

AL-Gard 3D Mask Guard

ALG-3DG

Helps ease breathing and prevents the mask from touching your mouth



3-Dimensional Design for Comfort + Function

3D ergonomic design **optimizes breathability** and **improves comfort** while separating the mask from contact with your mouth to **reduce condensation and allergic reactions**



Washable & Reusable

Reusable waterproof silicone material is hypoallergenic and can be conveniently washed with water and sanitized for long-lasting use



Secure Non-Slip Grip

Made specially for most disposable face masks with secure clips, the AL-Gard 3D Mask Guard will not slip out of place.



3D Mask Guard Fits All

Expert ergonomic design fits all users and is suitable for people of all ages and genders, including elderly and children



Quick & Easy Application Instructions



1

Align the 3D Mask Guard on the middle inside of the disposable face mask (facing you), with the clip at middle fold.



2

Securely slot and clip the ends of the 3D Mask Guard into the middle fold.



3

Expand the disposable mask and mould it to the concave shape of the 3D Mask Guard.



4

Wear the mask normally. Adjust the nose clip and attach ear loops.