

IBIS

- Unsurpassed adjustability
- Maximum user independence
- Accommodates changing needs
- Hand or foot propulsion



The IBIS Tilt-in-Space Manual Wheelchair takes the terms multi-functional and multi-adjustable to a whole new level. With its unique and versatile Sedeo Pro Seat System as well as its wide range of components and accessories, the IBIS is able to offer precise comfort and positioning, on-going adjustability for changing user needs, and precise chair settings for maximum user independence.

Un-surpassed, On-going Adjustability: The uniqueness of the IBIS begins with its Sedeo Pro Seat System. From the latin 'I Sit' the "Sedeo" Pro seat frame offers many built-in telescopic style adjustments that can be moved and set in place, making the IBIS a virtual "one size fits all" tilting wheelchair. These adjustments can be done immediately on the spot to ensure optimum user fit and allow for on-going adjustability if the user's needs change over time.

IBIS

IBIS

(with 12" rear wheels, attendant use)



Built-in Adjustability

- Seat width 16" - 22"
- Seat depth 16" - 23"
- Seat height 12.5" - 18.5"
- Overall chair width
- Back angle adjustment range of 39 degrees
- Armrest height & depth
- Static seat angle
- Rear wheel placement (centre of gravity)
- Many adjustments can be made with a "single tool"

Maximizing Independence

- The built-in adjustability of the IBIS goes to great lengths to ensure that the user has as much independence as possible. The chair's seat height can be set as low as 12.5" from the floor to allow for foot propelling.
- Hand propelling can also be optimized through choosing from multiple rear wheel sizes and rear wheel position settings. This ensures that the chair is properly balanced and the user can properly reach the rear wheels even in a tilted position.
- The unique Tilt-in-Space mechanism is designed to keep the user's knee height at a consistent level throughout the tilting range. This allows the user improved access to table tops and other work surfaces.
- By choosing the optional power tilt module, the user has the potential to change his or her position during the day with just the touch of a button. The hand control button can also be positioned for caregiver operation.

Considering the Caregiver

The IBIS also offers numerous features that benefit the caregiver:

- Adjustable height push handles ensure that the caregiver is comfortable when pushing the chair and operating the tilt-in-space mechanism.
- The gas spring assisted tilt-in-space mechanism allows for easy and smooth operation through a weight range of up to 350 lbs.



Many adjustments with one tool

Multi-adjustable Sedeo Pro Seat System



	IBIS
Seat width (adjustable)	*15" - 22"
Seat depth (adjustable)	**16" - 23"
Seat height (adjustable)	12.5" - 18.5"
Seat angle (adjustable)	0° - 6°
Armrest height (adjustable)	9" - 14"
Back angle (adjustable)	89° - 128°
Dynamic reclining back	Option
Folding back	Standard
Weight Capacity: Standard	275 lbs.
Weight Capacity: Upgrade	350 lbs.
Rear wheel sizes	12", 20", 22", 24"
Overall width (adjustable)	24.5" - 31"
Overall Length (excl. front rigging)	32"
Warranty	Base frame: Limited lifetime

* 15" seat width achieved with padded thigh supports/adductors option (see order form)

** possible tilt range limitations with seat depth over 20", when combined with low seat heights of 13.8" and less

Do You Need a Ramp?

