

ORION™ //



TILTING WHEELCHAIR
ORION II

 **FUTURE**
MOBILITY HEALTHCARE INC.



A



B



C



D



E



F

A. POSI -LOCK SYSTEM

GAS CYLINDER AND POSI-LOCK, LOCATED UNDERNEATH THE CHAIR, ALLOWS FOR A GENTLE, SLOW TILT FOR THE USER LOCKABLE AT INFINITE ANGLES.

B. TILT IN SPACE

ORION II OFFERS A TRUE 45 DEGREE TILT FOR OPTIMAL POSITIONING OF THE USER. IT HAS AN OPTIONAL 30 DEGREE OF RECLINE.

C. ATTENDANT FOOT TILT

ATTENDANT FOOT TILT ALLOWS THE WHEELCHAIR ATTENDANT TO EASILY TILT THE WHEELCHAIR BY PRESSING ON PEDAL TILTING IT UP OR DOWNWARDS.

D. ATTENDANT TRIGGER TILT

ALL BASE MODELS OF ORION II ARE EQUIPPED WITH THE TRIGGER TILT SYSTEM. THIS ALLOWS ATTENDANTS TO TILT THE WHEELCHAIR EASILY.

E. POWER TILT

POWER TILT IS AVAILABLE FOR EITHER THE ATTENDANT OR USER TO TILT THE WHEELCHAIR WITHOUT ANY ADDED MANUAL FORCE.

F. ADJUSTABLE AXLE AND SEAT HEIGHT

MULTIHOLE AXLE PLATE ALLOWS FOR VARIOUS REAR WHEEL POSITIONS GIVING THE USER CHOICE OF MORE OR LESS STABILITY. ALLOWS LOW SEAT HEIGHT OF 13 1/4".



TILTING

FEATURES ON THE ORION II WHEELCHAIR

- 30 DEGREE RECLINE OPTION
- LOW SEAT-TO-FLOOR HEIGHT
- STANDARD SEATING SYSTEM CONSISTING OF: TWO WAY CONTOURED PADDED FIXED HEADREST, ORION BACK WITH 1" DEPTH ADJUSTABILITY AND PRISM IDEAL CUSHION
- ALUMINUM SEAT PAN
- MULTI ADJUSTABLE AXLE PLATE
- REMOVABLE REAR WHEELS WITH CHROME OR PLASTIC COATED HANDRIMS
- ADJUSTABLE ARM HEIGHTS, GEL PADDING
- 70 DEGREE PINSTYLE FRONT RIGGING
- ALUMINUM OR COMPOSITE FOOTPLATES
- QUICK FLIP UP ANTI TIPPERS
- LOCKING COLLAR
- HEEL LOOPS
- ELEVATING LEGREST WITH CALF PAD, ALLOWS THE CLIENT POSITIONING OPTIONS FOR COMFORT AND PRESSURE RELIEF



CUSTOM ANYTHING!

CUSTOMER COMES FIRST. FUTURE MOBILITY HEALTHCARE INC. IS A CANADIAN MANUFACTURER OF WHEELCHAIRS, BACKS, CUSHIONS AND ACCESSORIES WITH OUR PARTS MANUFACTURED IN HOUSE. THIS GIVES US THE ABILITY TO CUSTOMIZE ANYTHING NEEDED FOR SPECIFIC CUSTOMER REQUESTS.

FUTURE MOBILITY IS PLEASE TO INTRODUCE THE ORION II, THE MOST INNOVATIVE AND AFFORDABLE TILT-IN-SPACE WHEELCHAIR ON THE MARKET TODAY. DEVELOPED FOR THOSE INDIVIDUALS REQUIRING ADJUSTABILITY TO REDISTRIBUTE THEIR BODY WEIGHT FOR OPTIMAL POSITIONING PURPOSES.














WHEELCHAIR



WWW.FUTUREMOBILITY.CA

ORION II

SPECIFICATIONS:

 SEAT WIDTH 15" - 24"	 SEAT TO FLOOR 13 1/4" - 20"	 REAR WHEEL 12", 20", 22" OR 24" COMPOSITE URETHANE
 SEAT DEPTH 15" - 24"	 OVERALL WIDTH SEAT WIDTH +9.5"	 CASTER SIZE 4", 5", 6", 7" AND 8" URETHANE
 FRONT RIGGING PINSTYLE, SWINGAWAY DUAL SWINGAWAY OR DUAL ELEVATING LEGRESTS	 ADJUSTABLE BACK ANGLE -7°, 0°, 7°, 14°, 21°	 BACK STYLE ORION BACK
 ARM STYLE ADJ. HEIGHT T-STYLE ADJ. HEIGHT PIN RELEASE T-STYLE ADJ. HEIGHT FLIP BACK	 PRODUCT WEIGHT 65 LBS (29 KG)	 WEIGHT CAPACITY 350 LBS (159 KG)
 ADP CODE Ontario WMA5FM005		

WARRANTY:

THE ORION II IS WARRANTED TO BE FREE FROM DEFECTS IN MATERIALS AND WORKMANSHIP FOR A PERIOD OF TWO (2) YEARS FROM THE DATE OF PURCHASE AND LIFETIME FOR THE FRAME. UPHOLSTERY FOR 90 DAYS, UNDER NORMAL USE BY ORIGINAL PURCHASER.

POWER TILT BATTERIES ARE WARRANTED FOR SIX (6) MONTHS WHILE THE CONTROL BOX, BATTERY CHARGER, ACTUATOR AND HAND CONTROL ARE WARRANTED FOR ONE (1) YEAR



FUTURE MOBILITY HEALTHCARE BELIEVES IN CONTINUED ENHANCEMENT, RESEARCH, INVENTION AND IMPROVEMENT OF MOBILITY PRODUCTS. OUR PRODUCTS CONTINUE TO IMPROVE TO MEET THE NEEDS OF OUR CUSTOMERS FOR A BETTER QUALITY OF LIFE.

FUTURE MOBILITY HEALTHCARE INC. 3223 ORLANDO DRIVE, MISSISSAUGA, ONTARIO, L4V 1C5
1-888-737-4011 FAX: 1-866-670-3377
WWW.FUTUREMOBILITY.CA
1 BUFFALO RIVER PLACE BUFFALO, NEW YORK, 14210, USA
1-855-9-FUTURE FAX: 716-783-9236,
WWW.FUTUREMOBILITY.COM
ALL RIGHTS RESERVED
PRINTED IN CANADA.
INNOVATIVE SEATING AND MOBILITY SOLUTIONS

