



DTK2

Deployment Trauma Kit



Product Overview

The DTK2 is designed to fit comfortably and securely behind ballistic plates or concealable armor. Undetectable to the wearer. The kit is made of durable RF Welded Vinyl and is completely waterproof.

Contents

Combat Gauze, Z-Fold Gauze x2, Occlusive Dressing x2, Decompression, NPA W/ 550, 3" Elastic Bandage, Nitrile Gloves, Casualty Card, 2" Safety Pins x2, Sharpie, Frog Tape, Eye Shield

Specs

Size: 9 x 11 x .75
Wt: 11.7 oz
Patents: 9980525, 9439819, and 8925115
Sealed, Plate Shaped, and 3/4 in thick

Features

- RF Welded Vinyl Container
- 100% Waterproof
- Reduced Bulk, Weight, and Stiffness
- Can be customized for unit buys
- Gov/LE Discounts Available

Additional Info

NSN: 6545-01-707-4139
GSA: GS-07F-5965P / Part # DTK-2
DAPA: Part # DTK-2
ECAT: Part # DTK-2
FedMall: Part # DTK-2

