



TRK

Team Response Kit



Product Overview

The Team Response Kit is a full-sized “CLS” waist pack with a padded hip belt and user-centric organization features. It has two modes of use: opening the top flap to access supplies while wearing the bag or placing the bag on the ground and opening it fully to access all supplies easily. The zipper-stop straps prevent accidental opening.

Contents

Combat Gauze x3, Tourniquet x3, Emergency Dressing 4in x3, Emergency Dressing Abdominal, Rolled Gauze x3, NPA x2, Lube x2, Chest Seal - Twin Pack x2, 14g Needle x2, Gloves x4 pairs, Eye Shield x2, Sharpie, TCCC Card x4, 2" Surgical Tape, Heat Sheet x2, SAM Splint, Triangular Cravat, Bandage x3, Strap Cutter, Scissor Leash, Shears - 7.5", Alcohol Prep Pad x6, Soft Litter

Specs

Size: 14" x 7" x 8"
Volume: 784 Cubic Inches
Pouch Weight: 1.28 lbs

Features

Carry handle on top, Loop panel for labels, Shear pocket with strap, Side compression straps, Spacer-mesh for ventilation, Front zippered pocket, Internal loops, Elastic loops inside lid, Internal slip pocket, Adjustable straps on bottom hold a pole-less litter or other gear, Suspender loops, Includes medical cross PVC patch, Grommets on the bottom for drainage, Made in USA, Berry Compliant

Additional Info

This product was built in partnership with High Speed Gear, a manufacturer, and distributor of 100 % American-made, Battle-Proven Tactical Gear.

