



## WCLB-2-PRO

### Wound CLUB™

### PRO Gen2



### Product Overview

The CLUB Pro Gen2 contains everything needed to run a 12 person wound packing class. One cube for every two students recommended.

### Contents

Each Pro Package comes with

- 1x Wound CLUB™
- 6x Wound Cubes™ Gen2
- 8x Frog Gauze™
- Advanced Blood Pumping System
- 1x Travel Case

### Specs

Material: Non-Toxic Silicone

Size: 4 x 4 x 4

Wt: 20 lbs

Patents Pending

### Features

- 2 in 1 Training Device (Junctional and Limb Trainer)
- Cube = Shoulder (Wound Packing Trainer)
- Tube = Upper Arm (Tourniquet, Wound Packing, and Pressure Dressing Trainer)
- Semi-transparent
- 4x Gunshot Wounds - Flat, Large, X-Large, and Bone In X-Large
- 2x Large Lacerations
- Realistic Blood Flow
- Realistic Tissue Loss and Tissue Density
- Anatomically Correct Humerus Bone
- Realistic Brachial Artery
- Compact Design (No wasted space)
- Durably built
- Affordable
- Blood simulant can be used with this product for added realism
- Gov/LE Discounts Available

