

## **Wound Simulators Product Information & Care Guide**

Due to the nature of soft silicone, debris will accumulate on Wound Simulation products.

## **Cleaning Instructions:**

Preferred method: Use warm water and a mild soap to clean after use.

Secondary method: Use 99% Isopropyl Alcohol.

Do not use abrasive cleaning products, as this may damage the product.

## **Durability Information:**

The high durometer of our Silicone Wound Cube allows us to stand behind a limited lifetime guarantee.

Some minor tearing or flaking around the wound pattern entry point will naturally occur and will not affect the use of the product.

## **Phokus Clear Blood:**

Phokus Clear Blood is a stainless lubricant designed for use with Phokus Wound Simulators.

Non-Toxic / Non-Hazardous











