



**MARLEN**  
**C A N A D A**

*Security by Design™*

[www.marlencanada.ca](http://www.marlencanada.ca)

2020 Edition



---

126-408 East Kent Avenue South  
Vancouver, BC V5X 2X7  
8am-4pm (Monday - Friday), PST.

**Orders can be placed via:**

Online: [nmswholesale.ca](http://nmswholesale.ca)

Email: [orders@nmswholesale.ca](mailto:orders@nmswholesale.ca)

Phone: 604.638.2761 | 1.844.379.9101

Fax: 604-879-3342



Product Overview	4
Barrier Ring	5
Protective Skin Barrier Sheets	5
<b>One-Piece Closed</b>	6
SuperFlat™ Closed No Tape	7
MaxPlus™ Closed No Tape	7
UltraMax™ DEEP Closed No Tape	8
UltraMax™ Closed No Tape	9
<b>One-Piece Drainable</b>	10
MegaMax™ Drainable No Tape with Kwick-Klose™	11
UltraMax™ DEEP Drainable No Tape with Kwick-Klose™	12
UltraMax™ Drainable No Tape with Kwick-Klose™	13
UltraLite™ Drainable with Tape & Kwick-Klose™	14
UltraLite™ Drainable with Tape & Snap-Lok™ Closure	15
MaxPlus™ Drainable No Tape Kwick-Klose™	16
MiniMax™ Drainable No Tape with Kwick-Klose™	16
SuperFlat™ MiniMax™ Drainable No Tape with Kwick-Klose™	17
SuperFlat™ Drainable No Tape with Kwick-Klose™	17
UltraExtra™ Drainable Tape with Snap-Lok™	18
UltraPlus™ Drainable with Tape Snap-Lok™	19
<b>One-Piece Urostomy</b>	20
UltraLite™ Urostomy with Tape	21
UltraMax™ Urostomy No Tape	22
UltraMax™ DEEP Urostomy No Tape	23
SuperFlat™ Urostomy No Tape	23
UltraExtra™ Urostomy with Tape	24
<b>Two-Piece Flanges and Pouches</b>	25
UltraMax™ Gemini™ No Tape Flange	26
UltraMax™ Gemini™ Pouch	27
<b>Two-Piece Reusable Ostomy Systems</b>	28
Flat All-Flexible	29
Convex All-Flexible	29
Deep Convex All-Flexible	30
Extra Large All-Flexible	30
Semi-Flexible	30
Two-Piece Drainable Pouch (Reusable)	31
Two-Piece Urostomy Pouch (Reusable)	32
<b>Ostomy Accessories</b>	33
Accessories	34
Urostomy Accessories	38
<b>Numerical Index</b>	39



# You're Always Covered Under the MARLEN Umbrella

You've always turned to MARLEN for the difficult to fit patient, the odd-sized stoma, the pouch size that nobody else seems to have. It's time to think of MARLEN CANADA as the universal solution for all your clients. We offer the largest combination of pouch types, convexities, sizes and body flanges in the industry.

## MegaMax™

- AquaTack™ Hydrocolloid Barrier
- No tape barrier
- Available in round or oval stoma openings in flat
- Transparent or opaque with Kwick-View™ window
- Pre-cut and cut-to-fit
- Belt tabs



## UltraMax™ DEEP

- Full barrier, no tape
- Flexible
- AquaTack™ Hydrocolloid Barrier
- Deep convex only
- Pre-cut and cut-to-fit
- Drainable pouch in transparent or opaque, no filter
- Closed pouch in transparent or opaque with filter
- Urostomy pouch, transparent only
- Drainable large flange option in transparent or opaque, no filter



## UltraMax™

- Flexible
- Full barrier, no tape
- Flat, or shallow cut-to-fit or pre-cut
- Transparent or opaque with Kwick-View™ window
- Closed, Drainable, Urostomy
- Closed with filter
- AquaTack™ Hydrocolloid Barrier



## UltraLite™

- Taped barrier
- Flexible
- AquaTack™ Hydrocolloid or Skin Shield™ Adhesive Skin Barrier
- Transparent or opaque with Kwick-View™ window
- Convex pre-cut only
- Flat, shallow, deep convex
- Filter only available in closed pouch
- Closed, Drainable, Urostomy



## MiniMax™

- Flat or shallow, cut-to-fit
- Filter
- Full barrier, no tape
- AquaTack™ Hydrocolloid Barrier
- Transparent or opaque with Kwick-View™ window
- 7.5" pouch length, 30% smaller than a conventional pouch



## SuperFlat™

- Extremely flat - 50% thinner than standard flange
- Very flexible
- AquaTack™ Hydrocolloid Barrier
- Standard length and 7.5" length
- With or without filter
- Closed, Drainable, Urostomy
- Cut-to-fit only



## MaxPlus™

- Full AquaTack™ Hydrocolloid Barrier, no tape
- Flat cut-to-fit to 80mm
- Transparent or opaque with Kwick-View™ window
- With or without filter
- Closed or Drainable
- Belt tabs
- Cutting surface has round and oval template



## UltraExtra™

- Skin Shield™ Adhesive Skin Barrier
- Flat or shallow convexity
- Cut-to-fit flat
- Round or oval pre-cut
- Transparent or opaque
- No filter
- Belt tabs
- 12" pouch length





## Barrier Ring

Picking the right barrier ring just got easier. MARLEN offers choices.

### MaxSeal™ Protective Barrier Ring

Easily moldable and offering maximum comfort, MaxSeal™ Protective Barrier Ring provides performance features similar to the leading barrier ring on the market at a value price. The ideal solution for improving the seal around the stoma. Resists fluids. Adheres well while also minimizing skin irritation.

DESCRIPTION	PRODUCT NO.
Large MaxSeal™ – 3.9" (100mm) outside diameter	9800
Regular MaxSeal™ – 2" (50mm) outside diameter	9700
5/32" (4mm) thickness	
Box of 10	



### UltraSeal™ Hydrocolloid Barrier

Bend, stretch, twist, pull. Flexible UltraSeal™ adjusts to accommodate most stomas, regardless of size and shape. Alcohol-free UltraSeal™ molds easily to fill in crevices and irregularities. The ring adheres well to both the body and the ostomy appliance, creating an ultra secure seal that prevents leakage and extends pouch wear time. UltraSeal™ removes easily, eliminating skin irritation.

DESCRIPTION	PRODUCT NO.
Large UltraSeal™ – 3.9" (100mm) outside diameter	9600 Box/5
Regular UltraSeal™ – 2" (50mm) outside diameter	9500 Box/10
3/32" (2mm) thickness	



## Protective Skin Barrier Sheets

MARLEN offers two types of barriers.

### AquaTack™ Hydrocolloid Barrier

- Absorbs skin moisture
- Adheres to moist skin
- Seals securely; removes easily from sensitive skin (controlled adhesion)



### Skin Shield™ Adhesive Skin Barrier

- Stays soft
- Elastic and conformable
- Will not absorb drainages
- Offers long-lasting protection





## One-Piece Closed

---

SuperFlat™ Closed No Tape	7
MaxPlus™ Closed No Tape	7
UltraMax™ DEEP Closed No Tape	8
UltraMax™ Closed No Tape	9





## SuperFlat™ Closed No Tape

*Marlen's Flattest, Most Discreet Pouch Ever*

- AquaTack™ Hydrocolloid Barrier
- Extremely flat - 50% thinner than standard flange
- Cut-to-fit 12-55mm
- Transparent or opaque with Kwick-View™ window
- Filter
- No belt tabs

Flat (No Convexity)	
DESCRIPTION	PRODUCT NO.
Opaque With Filter	83675
Transparent With Filter	82675
Box of 15	



*Cut-to-fit Flange*

## MaxPlus™ Closed No Tape

*One-Piece Disposable Ostomy System for Patients with LARGER Stomas*

- AquaTack™ Hydrocolloid Barrier
- Flat
- Cut-to-fit 12-80mm
- Transparent or opaque with Kwick-View™ window
- Filter
- Belt tabs
- Round or oval cutting guide on barrier surface

Flat (No Convexity)	
DESCRIPTION	PRODUCT NO.
Opaque With Filter	83680
Transparent With Filter	82680
Box of 15	




*Cut-to-fit Flange*

# UltraMax™ DEEP Closed No Tape

For MAXimum Choice and MAXimum Flexibility in a One-Piece System

- AquaTack™ Hydrocolloid Barrier
- Deep convexity
- Pre-cut 12-44mm
- Cut-to-fit 12-38mm
- Transparent or opaque with Kwick-View™ window
- Filter
- Belt tabs
- Also available in flat and shallow

 Deep Convexity		PRODUCT NUMBERS	
		OPAQUE	TRANSPARENT
STOMA OPENING			
Cut-to-fit	12-38 mm	83690	82690
1/2"	12 mm	83612	82612
5/8"	16 mm	83616	82616
3/4"	19 mm	83619	82619
7/8"	22 mm	83622	82622
1"	25 mm	83625	82625
1-1/8"	29 mm	83629	82629
1-1/4"	32 mm	83632	82632
1-3/8"	34 mm	83634	82634
1-1/2"	38 mm	83638	82638
1-5/8"	41 mm	83641	82641
1-3/4"	44 mm	83644	82644

Box of 15



*Cut-to-fit Flange*



*Pre-cut Flange*



# UltraMax™ Closed No Tape



For MAXimum Choice and MAXimum Flexibility in a One-Piece System

- AquaTack™ Hydrocolloid Barrier
- Flat or Shallow convexity
- Cut-to-fit flat 12-55mm
- Cut-to-fit convex 12-38mm
- Transparent or opaque with Kwick-View™ window
- Filter
- Belt tabs
- Also available in deep

Flat (No Convexity)		PRODUCT NUMBERS	
		OPAQUE	TRANSPARENT
STOMA OPENING			
Cut-to-fit	12-55 mm	83400	82400
1/2"	12 mm	83412	82412
5/8"	16 mm	83416	82416
3/4"	19 mm	83419	82419
7/8"	22 mm	83422	82422
1"	25 mm	83425	82425
1-1/8"	29 mm	83429	82429
1-1/4"	32 mm	83432	82432
1-3/8"	34 mm	83434	82434
1-1/2"	38 mm	83438	82438
1-5/8"	41 mm	83441	82441
1-3/4"	44 mm	83444	82444
1-7/8"	48 mm	83448	82448
2"	50 mm	83450	82450
2-1/8"	55 mm	83455	82455
2-1/4"	58 mm	83458	82458

Box of 15



*Cut-to-fit Flange*



*Pre-cut Flange*

Shallow Convexity		PRODUCT NUMBERS	
		OPAQUE	TRANSPARENT
STOMA OPENING			
Cut-to-fit	12-38 mm	83500	82500
1/2"	12 mm	83512	82512
5/8"	16 mm	83516	82516
3/4"	19 mm	83519	82519
7/8"	22 mm	83522	82522
1"	25 mm	83525	82525
1-1/8"	29 mm	83529	82529
1-1/4"	32 mm	83532	82532
1-3/8"	34 mm	83534	82534
1-1/2"	38 mm	83538	82538
1-5/8"	41 mm	83541	82541
1-3/4"	44 mm	83544	82544

Box of 15



*Cut-to-fit Flange*



*Pre-cut Flange*



## One-Piece Drainable

---

MegaMax™ Drainable No Tape with Kwick-Klose™	11
UltraMax™ DEEP Drainable No Tape with Kwick-Klose™	12
UltraMax™ Drainable No Tape with Kwick-Klose™	13
UltraLite™ Drainable with Tape & Kwick Klose™	14
UltraLite™ Drainable with Tape & Snap-Lok™ Closure	15
MaxPlus™ Drainable No Tape Kwick-Klose™	16
MiniMax™ Drainable No Tape with Kwick-Klose™	16
SuperFlat™ MiniMax™ Drainable No Tape with Kwick-Klose™	17
SuperFlat™ Drainable No Tape with Kwick-Klose™	17
Ultra™ Extra Drainable Tape with Snap-Lok™	18
UltraPlus™ Drainable with Tape Snap-Lok™	19





Marlen's largest pouch with a larger flange and larger cutting surface (up to 100mm)

- AquaTack™ Hydrocolloid Barrier
- Available in round or oval stoma openings in flat only
- Round Cut-to-fit (50-100mm)
- Oval Cut-to-fit (12-100mm)
- Transparent or opaque with Kwick-View™ window
- Belt Tabs

Flat (No Convexity)		PRODUCT NUMBERS	
		OPAQUE	TRANSPARENT
ROUND STOMA OPENING			
Cut-to-fit	12-100 mm	54000	54025
2"	50 mm	54001	54026
2-1/8"	54 mm	54002	54027
2-1/4"	57 mm	54003	54028
2-3/8"	60 mm	54004	54029
2-1/2"	63 mm	54005	54030
2-5/8"	67 mm	54006	54031
2-3/4"	70 mm	54007	54032
2-7/8"	73 mm	54008	54033
3"	76 mm	54009	54034
3-1/8"	79 mm	54010	54035
3-1/4"	82 mm	54011	54036
3-3/8"	85 mm	54012	54037
3-1/2"	88 mm	54013	54038
3-3/4"	95 mm	54014	54039
4"	100 mm	54015	54040

Box of 5



*Cut-to-fit Flange*

Flat (No Convexity)		PRODUCT NUMBERS	
		OPAQUE	TRANSPARENT
OVAL STOMA OPENING			
Cut-to-fit	12-100 mm	54000	54025
1" x 1-1/2"	25x38 mm	54016	54041
1" x 1-3/4"	25x44 mm	54017	54042
1-1/4" x 1-1/2"	32x38 mm	54018	54043
1-1/4" x 2"	32x50 mm	54019	54044
1-1/2" x 2"	38x50 mm	54020	54045
1-1/2" x 2-1/2"	38x62 mm	54021	54046
1-1/2" x 3"	38x76 mm	54022	54046
2" x 3"	50x76 mm	54023	54048
2" x 3-1/2"	50x87 mm	54024	54049

Box of 5




*Pre-cut Oval Flange*



# UltraMax™ DEEP Drainable No Tape with Kwick-Klose™

For MAXimum Choice and MAXimum Flexibility in a One-Piece System

- AquaTack™ Hydrocolloid Barrier
- Deep convexity
- Pre-cut convex 12-44mm
- Cut-to-fit convex 12-42mm & 12-64mm
- Transparent or opaque with Kwick-View™ window
- No filter
- Belt tabs
- Also available in flat and shallow

 Deep Convexity		PRODUCT NUMBERS	
		OPAQUE	TRANSPARENT
STOMA OPENING			
Cut-to-fit	12-42 mm	53690	52690
Cut-to-fit	12-64 mm	53699	52699
1/2"	12 mm	53612	52612
5/8"	16 mm	53616	52616
3/4"	19 mm	53619	52619
7/8"	22 mm	53622	52622
1"	25 mm	53625	52625
1-1/8"	29 mm	53629	52629
1-1/4"	32 mm	53632	52632
1-3/8"	34 mm	53634	52634
1-1/2"	38 mm	53638	52638
1-5/8"	41 mm	53641	52641
1-3/4"	44 mm	53644	52644
Box of 5			



*Cut-to-fit 12-42mm  
Flange*



*Cut-to-fit 12-64mm  
Flange*



*Pre-cut Flange*

# UltraMax™ Drainable No Tape with Kwick-Klose™



For MAXimum Choice and MAXimum Flexibility in a One-Piece System

- AquaTack™ Hydrocolloid Barrier
- Flat or shallow convexity
- Cut-to-fit flat 12-55mm
- Pre-cut flat 12-58mm
- Pre-cut convex 12-44mm
- Cut-to-fit convex 12-38mm
- Transparent or opaque with Kwick-View™ window
- No filter
- Belt tabs
- Also available in deep

Flat (No Convexity)		PRODUCT NUMBERS	
		OPAQUE	TRANSPARENT
STOMA OPENING			
Cut-to-fit	12-55 mm	53400	52400
1/2"	12 mm	53412	52412
5/8"	16 mm	53416	52416
3/4"	19 mm	53419	52419
7/8"	22 mm	53422	52422
1"	25 mm	53425	52425
1-1/8"	29 mm	53429	52429
1-1/4"	32 mm	53432	52432
1-3/8"	34 mm	53434	52434
1-1/2"	38 mm	53438	52438
1-5/8"	41 mm	53441	52441
1-3/4"	44 mm	53444	52444
1-7/8"	48 mm	53448	52448
2"	50 mm	53450	52450
2-1/8"	55 mm	53455	52455
2-1/4"	58 mm	53458	52458

Box of 5



Cut-to-fit Flange



Pre-cut Flange

Shallow Convexity		PRODUCT NUMBERS	
		OPAQUE	TRANSPARENT
STOMA OPENING			
Cut-to-fit	12-38 mm	53500	52500
1/2"	12 mm	53512	52512
5/8"	16 mm	53516	52516
3/4"	19 mm	53519	52519
7/8"	22 mm	53522	52522
1"	25 mm	53525	52525
1-1/8"	29 mm	53529	52529
1-1/4"	32 mm	53532	52532
1-3/8"	34 mm	53534	52534
1-1/2"	38 mm	53538	52538
1-5/8"	41 mm	53541	52541
1-3/4"	44 mm	53544	52544

Box of 5



Cut-to-fit Flange



Pre-cut Flange



# UltraLite™ Drainable with Tape & Kwick-Klose™

## One-Piece Disposable Ostomy System with a Slim, Contoured Design

- AquaTack™ Hydrocolloid Barrier or Skin Shield™ Adhesive Skin Barrier
- Flat, shallow or deep convexity
- Cut-to-fit flat 12-55mm
- Pre-cut flat 12-44mm
- Pre-cut convex 12mm-44mm
- Transparent or opaque with Kwick-View™ window
- No filter
- Belt tabs

Flat (No Convexity)		PRODUCT NUMBERS			
		SKIN SHIELD™ BARRIER		AQUATAK™ BARRIER	
STOMA OPENING		OPAQUE	TRANSPARENT	OPAQUE	TRANSPARENT
Cut-to-fit	12-55 mm	50400	51400	56400	57400
1/2"	12 mm	50412	51412	56412	57412
5/8"	16 mm	50416	51416	56416	57416
3/4"	19 mm	50419	51419	56419	57419
7/8"	22 mm	50422	51422	56422	57422
1"	25 mm	50425	51425	56425	57425
1-1/8"	29 mm	50429	51429	56429	57429
1-1/4"	32 mm	50432	51432	56432	57432
1-3/8"	34 mm	50434	51434	56434	57434
1-1/2"	38 mm	50438	51438	56438	57438
1-5/8"	41 mm	50441	51441	56441	57441
1-3/4"	44 mm	50444	51444	56444	57444

Box of 10



Cut-to-fit  
Flange

Pre-cut  
Flange

Shallow Convexity		PRODUCT NUMBERS			
		SKIN SHIELD™ BARRIER		AQUATAK™ BARRIER	
STOMA OPENING		OPAQUE	TRANSPARENT	OPAQUE	TRANSPARENT
1/2"	12 mm	50512	51512	56512	57512
5/8"	16 mm	50516	51516	56516	57516
3/4"	19 mm	50519	51519	56519	57519
7/8"	22 mm	50522	51522	56522	57522
1"	25 mm	50525	51525	56525	57525
1-1/8"	29 mm	50529	51529	56529	57529
1-1/4"	32 mm	50532	51532	56532	57532
1-3/8"	34 mm	50534	51534	56534	57534
1-1/2"	38 mm	50538	51538	56538	57538
1-5/8"	41 mm	50541	51541	56541	57541
1-3/4"	44 mm	50544	51544	56544	57544

Box of 10



Pre-cut Flange

Deep Convexity		PRODUCT NUMBERS			
		SKIN SHIELD™ BARRIER		AQUATAK™ BARRIER	
STOMA OPENING		OPAQUE	TRANSPARENT	OPAQUE	TRANSPARENT
1/2"	12 mm	50612	51612	56612	57612
5/8"	16 mm	50616	51616	56616	57616
3/4"	19 mm	50619	51619	56619	57619
7/8"	22 mm	50622	51622	56622	57622
1"	25 mm	50625	51625	56625	57625
1-1/8"	29 mm	50629	51629	56629	57629
1-1/4"	32 mm	50632	51632	56632	57632
1-3/8"	34 mm	50634	51634	56634	57634
1-1/2"	38 mm	50638	51638	56638	57638
1-5/8"	41 mm	50641	51641	56641	57641
1-3/4"	44 mm	50644	51644	56644	57644

Box of 10



Pre-cut Flange

# UltraLite™ Drainable with Tape & Snap-Lok™ Closure




## One-Piece Disposable Ostomy System with a Slim, Contoured Design

- AquaTack™ Hydrocolloid Barrier or Skin Shield™ Adhesive Skin Barrier
- Flat, shallow or deep convexity
- Cut-to-fit flat 12-55mm
- Pre-cut flat 12-44mm
- Pre-cut convex 12-44mm
- Transparent or opaque with Kwick-View™ window
- No filter
- Belt tabs

Flat (No Convexity)		PRODUCT NUMBERS			
		SKIN SHIELD™ BARRIER		AQUATAACK™ BARRIER	
STOMA OPENING		OPAQUE	TRANSPARENT	OPAQUE	TRANSPARENT
Cut-to-fit	12-55 mm	55000	55176	55088	55264
1/2"	12 mm	55001	55177	55089	55265
5/8"	16 mm	55002	55178	55090	55266
3/4"	19 mm	55003	55179	55091	55267
7/8"	22 mm	55004	55180	55092	55268
1"	25 mm	55005	55181	55093	55269
1-1/8"	29 mm	55006	55182	55094	55270
1-1/4"	32 mm	55007	55183	55095	55271
1-3/8"	34 mm	55008	55184	55096	55272
1-1/2"	38 mm	55009	55185	55097	55273
1-5/8"	41 mm	55010	55186	55098	55274
1-3/4"	44 mm	55011	55187	55099	55275

Box of 10



Shallow Convexity		PRODUCT NUMBERS			
		SKIN SHIELD™ BARRIER		AQUATAACK™ BARRIER	
STOMA OPENING		OPAQUE	TRANSPARENT	OPAQUE	TRANSPARENT
1/2"	12 mm	55030	55206	55118	55294
5/8"	16 mm	55031	55207	55119	55295
3/4"	19 mm	55032	55208	55120	55296
7/8"	22 mm	55033	55209	55121	55297
1"	25 mm	55034	55210	55122	55298
1-1/8"	29 mm	55035	55211	55123	55299
1-1/4"	32 mm	55036	55212	55124	55300
1-3/8"	34 mm	55037	55213	55125	55301
1-1/2"	38 mm	55038	55214	55126	55302
1-5/8"	41 mm	55039	55215	55127	55303
1-3/4"	44 mm	55040	55216	55128	55304

Box of 10



Deep Convexity		PRODUCT NUMBERS			
		SKIN SHIELD™ BARRIER		AQUATAACK™ BARRIER	
STOMA OPENING		OPAQUE	TRANSPARENT	OPAQUE	TRANSPARENT
1/2"	12 mm	55059	55235	55147	55323
5/8"	16 mm	55060	55236	55148	55324
3/4"	19 mm	55061	55237	55149	55325
7/8"	22 mm	55062	55238	55150	55326
1"	25 mm	55063	55239	55151	55327
1-1/8"	29 mm	55064	55240	55152	55328
1-1/4"	32 mm	55065	55241	55153	55329
1-3/8"	34 mm	55066	55242	55154	55330
1-1/2"	38 mm	55067	55243	55155	55331
1-5/8"	41 mm	55068	55244	55156	55332
1-3/4"	44 mm	55069	55245	55157	55333

Box of 10







## MaxPlus™ Drainable No Tape Kwick-Klose™

*One-Piece Disposable Ostomy System for Patients with LARGER Stomas*

- AquaTack™ Hydrocolloid Barrier
- Flat
- Cut-to-fit 12-80mm
- Transparent or opaque with Kwick-View™ window
- With or without filter
- Belt tabs

Flat (No Convexity)	
DESCRIPTION	PRODUCT NO.
Opaque With Filter	53680
Transparent With Filter	52680
Opaque Without Filter	53880
Transparent Without Filter	52880
Box of 5	



*Cut-to-fit Flange*

## MiniMax™ Drainable No Tape with Kwick-Klose™

*One-Piece Disposable Ostomy System For Patients Who Require a Smaller Pouch*

- AquaTack™ Hydrocolloid Barrier
- Flat or shallow convexity
- Cut-to-fit flat 12-58mm
- Cut-to-fit shallow 12-38mm
- Transparent or opaque with Kwick-View™ window
- With or without filter
- Belt tabs
- 7.5" pouch length

Flat (No Convexity)	PRODUCT NUMBERS	
	OPAQUE	TRANSPARENT
U-CUT STARTER HOLE		
With Filter	53600	52600
Without Filter	53800	52800
Box of 5		



*Cut-to-fit Flange*

Shallow Convexity	PRODUCT NUMBERS	
	OPAQUE	TRANSPARENT
U-CUT STARTER HOLE		
With Filter	53700	52700
Without Filter	53900	52900
Box of 5		



*Cut-to-fit Flange*



## SuperFlat™ MiniMax™ Drainable No Tape with Kwick-Klose™



*Marlen's Flattest, Most Discreet Pouch Ever*

- AquaTack™ Hydrocolloid Barrier
- Extremely flat - 50% thinner than standard flange
- Cut-to-fit 12-55mm
- Transparent or opaque with Kwick-View™ window
- With or without filter
- No belt tabs
- 7.5" pouch length

Flat (No Convexity)	
DESCRIPTION	PRODUCT NO.
Opaque With Filter	53670
Transparent With Filter	52670
Opaque Without Filter	53870
Transparent Without Filter	52870
Box of 5	



*Cut-to-fit Flange*

## SuperFlat™ Drainable No Tape with Kwick-Klose™

*Marlen's Flattest, Most Discreet Pouch Ever*

- AquaTack™ Hydrocolloid Barrier
- Extremely flat - 50% thinner than standard flange
- Cut-to-fit 12-55mm
- Transparent or opaque with Kwick-View™ window
- With or without filter
- No belt tabs
- 9" pouch length

Flat (No Convexity)	
DESCRIPTION	PRODUCT NO.
Opaque With Filter	53675
Transparent With Filter	52675
Opaque Without Filter	53875
Transparent Without Filter	52875
Box of 5	



*Cut-to-fit Flange*

# UltraExtra™ Drainable Tape with Snap-Lok™

## One-Piece Disposable System for Ileostomies, Colostomies and Urostomies

- Skin Shield™ Adhesive Skin Barrier
- Cut-to-fit flat 48-76mm
- Pre-cut flat 48-76mm
- Pre-cut shallow 48-76mm
- Transparent or opaque
- No filter
- Belt tabs
- 12" pouch length

Flat (No Convexity)		PRODUCT NUMBERS	
		STANDARD-SIZED POUCH	
STOMA OPENING		OPAQUE	TRANSPARENT
Cut-to-fit	48-76 mm	55506	55507
1-7/8"	48 mm	55100	55276
2"	50 mm	55101	55277
2-1/8"	54 mm	55102	55278
2-1/4"	57 mm	55103	55279
2-3/8"	60 mm	55104	55280
2-1/2"	63 mm	55105	55281
2-5/8"	67 mm	55106	55282
2-3/4"	70 mm	55107	55283
2-7/8"	73 mm	55108	55284
3"	76 mm	55109	55285
1" x 1-1/2"	25x38 mm Oval	55110	55286
1" x 1-3/4"	25x44 mm Oval	55111	55287
1-1/4" x 1-1/2"	32x38 mm Oval	55112	55288
1-1/4" x 2"	32x50 mm Oval	55113	55289
1-1/2" x 2"	38x50 mm Oval	55114	55290
1-1/2" x 2-1/2"	38x62 mm Oval	55115	55291
1-1/2" x 3"	38x76 mm Oval	55116	55292
2" x 3"	50x76 mm Oval	55117	55293
Box of 10			

Shallow Convexity		PRODUCT NUMBERS	
		STANDARD-SIZED POUCH	
STOMA OPENING		OPAQUE	TRANSPARENT
1-7/8"	48 mm	55129	55305
2"	50 mm	55130	55306
2-1/8"	54 mm	55131	55307
2-1/4"	57 mm	55132	55308
2-3/8"	60 mm	55133	55309
2-1/2"	63 mm	55134	55310
2-5/8"	67 mm	55135	55311
2-3/4"	70 mm	55136	55312
2-7/8"	73 mm	55137	55313
3"	76 mm	55138	55314
1" x 1-1/2"	25x38 mm Oval	55139	55315
1" x 1-3/4"	25x44 mm Oval	55140	55316
1-1/4" x 1-1/2"	32x38 mm Oval	55141	55317
1-1/4" x 2"	32x50 mm Oval	55142	55318
1-1/2" x 2"	38x50 mm Oval	55143	55319
1-1/2" x 2-1/2"	38x62 mm Oval	55144	55320
1-1/2" x 3"	38x76 mm Oval	55145	55321
2" x 3"	50x76 mm Oval	55146	55322
Box of 10			



Cut-to-fit Flange

Pre-cut Flange



Pre-cut Flange

# UltraPlus™ Drainable with Tape Snap-Lok™



## One-Piece Disposable System for Ileostomies, Colostomies and Urostomies

- Skin Shield™ Adhesive Skin Barrier with large tape-footprint
- Cut-to-fit flat 50-88mm
- Pre-cut flat 50-88mm
- Pre-cut oval flat 25-87mm
- Pre-cut convex 50-88mm
- Pre-cut oval convex 25-87mm
- Transparent or opaque
- No filter
- Belt tabs
- 13.75" pouch length

Flat (No Convexity)		PRODUCT NUMBERS	
		OPAQUE	TRANSPARENT
STOMA OPENING			
Cut-to-fit	50-88 mm	55352	55426
2"	50 mm	55353	55427
2-1/8"	54 mm	55354	55428
2-1/4"	57 mm	55355	55429
2-3/8"	60 mm	55356	55430
2-1/2"	63 mm	55357	55431
2-5/8"	67 mm	55358	55432
2-3/4"	70 mm	55359	55433
2-7/8"	73 mm	55360	55434
3"	76 mm	55361	55435
3-1/8"	79 mm	55362	55436
3-1/4"	82 mm	55508	55521
3-3/8"	85 mm	55509	55522
3-1/2"	88 mm	55510	55523
1" x 1-1/2"	25x38 mm Oval	55363	55437
1" x 1-3/4"	25x44 mm Oval	55364	55438
1-1/4" x 1-1/2"	32x38 mm Oval	55365	55439
1-1/4" x 2"	32x50 mm Oval	55366	55440
1-1/2" x 2"	38x50 mm Oval	55367	55441
1-1/2" x 2-1/2"	38x62 mm Oval	55533	55537
1-1/2" x 3"	38x76 mm Oval	55368	55442
2"x3"	50x76 mm Oval	55369	55443
2" x 3-1/2"	50x87 mm Oval	55370	55444

Box of 10



Cut-to-fit  
Flange



Pre-cut  
Flange

Shallow Convexity		PRODUCT NUMBERS	
		OPAQUE	TRANSPARENT
STOMA OPENING			
2"	50 mm	55371	55445
2-1/8"	54 mm	55372	55446
2-1/4"	57 mm	55373	55447
2-3/8"	60 mm	55374	55448
2-1/2"	63 mm	55375	55449
2-5/8"	67 mm	55376	55450
2-3/4"	70 mm	55377	55451
2-7/8"	73 mm	55378	55452
3"	76 mm	55379	55453
3-1/8"	79 mm	55380	55454
3-1/4"	82 mm	55512	55524
3-3/8"	85 mm	55513	55525
3-1/2"	88 mm	55514	55526
1" x 1-1/2"	25x38 mm Oval	55381	55455
1" x 1-3/4"	25x44 mm Oval	55382	55456
1-1/4" x 1-1/2"	32x38 mm Oval	55383	55457
1-1/4" x 2"	32x50 mm Oval	55384	55458
1-1/2" x 2"	38x50 mm Oval	55385	55459
1-1/2" x 2-1/2"	38x62 mm Oval	55534	55538
1-1/2" x 3"	38x76 mm Oval	55386	55460
2" x 3"	50x76 mm Oval	55387	55461
2" x 3-1/2"	50x87 mm Oval	55388	55462

Box of 10



Pre-cut  
Flange



## One-Piece Urostomy

---

UltraLite™ Urostomy with Tape	21
UltraMax™ Urostomy No Tape	22
UltraMax™ DEEP Urostomy No Tape	23
SuperFlat™ Urostomy No Tape	23
UltraExtra™ Urostomy with Tape	24



# UltraLite™ Urostomy with Tape



## One-Piece Disposable Ostomy System with a Slim, Contoured Design

- Push-Pull E-Z Drain Valve
- AquaTack™ Hydrocolloid Barrier or Skin Shield™ Adhesive Skin Barrier
- Flat, shallow, or deep convexity
- Cut-to-fit flat 12-44mm
- Pre-cut flat 12-44mm
- Pre-cut convex 12-44mm
- Transparent with cloth-like backing
- Belt tabs

Flat (No Convexity)		PRODUCT NUMBERS	
		SKIN SHIELD™ BARRIER	AQUATAACK™ BARRIER
STOMA OPENING		TRANSPARENT	TRANSPARENT
Cut-to-fit	12-44 mm	77000	77088
1/2"	12 mm	77001	77089
5/8"	16 mm	77002	77090
3/4"	19 mm	77003	77091
7/8"	22 mm	77004	77092
1"	25 mm	77005	77093
1-1/8"	29 mm	77006	77094
1-1/4"	32 mm	77007	77095
1-3/8"	34 mm	77008	77096
1-1/2"	38 mm	77009	77097
1-5/8"	41 mm	77010	77098
1-3/4"	44 mm	77011	77099

Box of 10

Shallow Convexity		PRODUCT NUMBERS	
		SKIN SHIELD™ BARRIER	AQUATAACK™ BARRIER
STOMA OPENING		TRANSPARENT	TRANSPARENT
1/2"	12 mm	77030	77118
5/8"	16 mm	77031	77119
3/4"	19 mm	77032	77120
7/8"	22 mm	77033	77121
1"	25 mm	77034	77122
1-1/8"	29 mm	77035	77123
1-1/4"	32 mm	77036	77124
1-3/8"	34 mm	77037	77125
1-1/2"	38 mm	77038	77126
1-5/8"	41 mm	77039	77127
1-3/4"	44 mm	77040	77128

Box of 10

Deep Convexity		PRODUCT NUMBERS	
		SKIN SHIELD™ BARRIER	AQUATAACK™ BARRIER
STOMA OPENING		TRANSPARENT	TRANSPARENT
1/2"	12 mm	77059	77147
5/8"	16 mm	77060	77148
3/4"	19 mm	77061	77149
7/8"	22 mm	77062	77150
1"	25 mm	77063	77151
1-1/8"	29 mm	77064	77152
1-1/4"	32 mm	77065	77153
1-3/8"	34 mm	77066	77154
1-1/2"	38 mm	77067	77155
1-5/8"	41 mm	77068	77156
1-3/4"	44 mm	77069	77157

Box of 10



# UltraMax™ Urostomy No Tape

For MAXimum Choice and MAXimum Flexibility in a One-Piece System

- Push-Pull E-Z Drain Valve
- AquaTack™ Hydrocolloid Barrier
- Flat or shallow convexity
- Pre-cut flat 12-55mm
- Cut-to-fit flat 12-55mm
- Pre-cut convex 12-44mm
- Cut-to-fit convex 12-38mm
- Transparent with cloth-like backing
- Belt tabs
- Also available in deep

Flat (No Convexity)		PRODUCT NUMBERS
STOMA OPENING		TRANSPARENT
Cut-to-fit	12-55 mm	72400
1/2"	12 mm	72412
5/8"	16 mm	72416
3/4"	19 mm	72419
7/8"	22 mm	72422
1"	25 mm	72425
1-1/8"	29 mm	72429
1-1/4"	32 mm	72432
1-3/8"	34 mm	72434
1-1/2"	38 mm	72438
1-5/8"	41 mm	72441
1-3/4"	44 mm	72444
1-7/8"	48 mm	72448
2"	50 mm	72450
2-1/8"	55 mm	72455
2-1/4"	58 mm	72458

Box of 5



Cut-to-fit Flange

Pre-cut Flange

Shallow Convexity		PRODUCT NUMBERS
STOMA OPENING		TRANSPARENT
Cut-to-fit	12-38 mm	72500
1/2"	12 mm	72512
5/8"	16 mm	72516
3/4"	19 mm	72519
7/8"	22 mm	72522
1"	25 mm	72525
1-1/8"	29 mm	72529
1-1/4"	32 mm	72532
1-3/8"	34 mm	72534
1-1/2"	38 mm	72538
1-5/8"	41 mm	72541
1-3/4"	44 mm	72544

Box of 5



Cut-to-fit Flange

Pre-cut Flange





## UltraMax™ DEEP Urostomy No Tape

*For MAXimum Choice and MAXimum Flexibility in a One-Piece System*

- Push-Pull E-Z Drain Valve
- AquaTack™ Hydrocolloid Barrier
- Deep convexity
- Cut-to-fit convex 12-38mm
- Pre-cut convex 12-44mm
- Transparent with cloth-like backing
- Belt tabs
- Also available in flat and shallow

Deep Convexity		PRODUCT NUMBERS
STOMA OPENING		TRANSPARENT
Cut-to-fit	12-38 mm	72690
1/2"	12 mm	72612
5/8"	16 mm	72616
3/4"	19 mm	72619
7/8"	22 mm	72622
1"	25 mm	72625
1-1/8"	29 mm	72629
1-1/4"	32 mm	72632
1-3/8"	34 mm	72634
1-1/2"	38 mm	72638
1-5/8"	41 mm	72641
1-3/4"	44 mm	72644

Box of 5



*Cut-to-fit Flange*



*Pre-cut Flange*

## SuperFlat™ Urostomy No Tape

*Marlen's Flattest, Most Discreet Pouch Ever*

- Push-Pull E-Z Drain Valve
- AquaTack™ Hydrocolloid Barrier
- Extremely flat - 50% thinner than standard flange
- Cut-to-fit 12-55mm
- Transparent or opaque with cloth-like backing
- No belt tabs

Flat (No Convexity)	
DESCRIPTION	PRODUCT NO.
Opaque	73875
Transparent	72875

Box of 5



*Cut-to-fit Flange*



# UltraExtra™ Urostomy with Tape

## One-Piece Disposable System for Urostomies

- Push-Pull E-Z Drain Valve
- Skin Shield™ Adhesive Skin Barrier
- Flat or shallow convexity
- Cut-to-fit flat 48-76mm
- Pre-cut flat 48-76mm
- Pre-cut flat oval 25-76mm
- Pre-cut shallow 48-76mm
- Pre-cut shallow oval 25-76mm
- Belt tabs
- 12" pouch length

Flat (No Convexity)		PRODUCT NUMBERS
STOMA OPENING		TRANSPARENT
Cut-to-fit	48-76 mm	77178
1-7/8"	48 mm	77012
2"	50 mm	77013
2-1/8"	55 mm	77014
2-1/4"	58 mm	77015
2-3/8"	60 mm	77016
2-1/2"	63 mm	77017
2-5/8"	67 mm	77018
2-3/4"	70 mm	77019
2-7/8"	73 mm	77020
3"	76 mm	77021
1" x 1-1/2"	25x37 mm Oval	77022
1" x 1-3/4"	25x44 mm Oval	77023
1-1/4" x 1-1/2"	32x38 mm Oval	77024
1-1/4" x 2"	32x50 mm Oval	77025
1-1/2" x 2"	38x55 mm Oval	77026
1-1/2" x 2-1/2"	38x62 mm Oval	77027
1-1/2" x 3"	38x76 mm Oval	77028
2" x 3"	50x76 mm Oval	77029
Box of 10		



*Cut-to-fit Flange*

Shallow Convexity		PRODUCT NUMBERS
STOMA OPENING		TRANSPARENT
1-7/8"	48 mm	77041
2"	50 mm	77042
2-1/8"	55 mm	77043
2-1/4"	58 mm	77044
2-3/8"	60 mm	77045
2-1/2"	63 mm	77046
2-5/8"	67 mm	77047
2-3/4"	70 mm	77048
2-7/8"	73 mm	77049
3"	76 mm	77050
1" x 1-1/2"	25x37 mm Oval	77051
1" x 1-3/4"	25x44 mm Oval	77052
1-1/4" x 1-1/2"	32x38 mm Oval	77053
1-1/4" x 2"	32x50 mm Oval	77054
1-1/2" x 2"	38x55 mm Oval	77055
1-1/2" x 2-1/2"	38x62 mm Oval	77056
1-1/2" x 3"	38x76 mm Oval	77057
2" x 3"	50x76 mm Oval	77058
Box of 10		



*Pre-cut Flange*





## Two-Piece Flanges and Pouches

---

UltraMax™ Gemini™ No Tape Flange	26
UltraMax™ Gemini™ Pouch	27



# UltraMax™ Gemini™ No Tape Flange

- AquaTack™ Hydrocolloid Barrier
- Flat or shallow convexity
- Cut-to-fit flat 12-55mm
- Pre-cut flat 12-38mm
- Cut-to-fit shallow 12-38mm
- Pre-cut shallow 12-38mm
- Belt tabs
- Raised coupling ring with “E-Z Guide™” Locating System
- Adhesive coupling



Flat (No Convexity)	STOMA OPENING
	Cut-to-Fit
STOMA OPENING WIDTH	PRODUCT NO.
w/ 1/2" (12mm-55mm) starter hole	20000
Box of 5	



Shallow Convexity	STOMA OPENING
	Cut-to-Fit
STOMA OPENING WIDTH	PRODUCT NO.
w/ 1/2" (12mm-38mm) starter hole	21000
Box of 5	



Flat (No Convexity)	STOMA OPENING	
	Pre-Cut	
W/PRE-CUT OPENINGS	PRODUCT NO.	
1/2"      12 mm	20012	
5/8"      16 mm	20016	
3/4"      19 mm	20019	
7/8"      22 mm	20022	
1"        25 mm	20025	
1-1/8"    29 mm	20029	
1-1/4"    32 mm	20032	
1-3/8"    34mm	20034	
1-1/2"    38 mm	20038	
Box of 5		



Shallow Convexity	STOMA OPENING	
	Pre-Cut	
W/PRE-CUT OPENINGS	PRODUCT NO.	
1/2"      12 mm	21012	
5/8"      16 mm	21016	
3/4"      19 mm	21019	
7/8"      22 mm	21022	
1"        25 mm	21025	
1-1/8"    29 mm	21029	
1-1/4"    32 mm	21032	
1-3/8"    34mm	21034	
1-1/2"    38 mm	21038	
Box of 5		

# UltraMax™ Gemini™ Pouch



With the E-Z Guide™ System, adhesive coupling

## UltraMax™ Gemini™ Closed-End

- Transparent or opaque with Kwick-View™ window
- Filter
- Belt tabs

DESCRIPTION	PRODUCT NO.
Opaque with Kwick-View™	25100
Transparent	25150
<b>Box of 30</b>	



## UltraMax™ Gemini™ Ileostomy (Drainable)

- Transparent or opaque with Kwick-View™ window
- Kwick-Klose™ closure
- With or without filter
- Belt tabs

DESCRIPTION	PRODUCT NO.
Opaque with Kwick-View™, with filter	25250
Transparent, with filter	25260
Opaque with Kwick-View™, without filter	25200
Transparent, without filter	25210
<b>Box of 10</b>	



## MARLEN UltraMax™ Gemini™ Urostomy

- Push-Pull E-Z Drain Valve
- Transparent or opaque with Kwick-View™ window
- Belt tabs

DESCRIPTION	PRODUCT NO.
Opaque with Kwick-View™	25300
Transparent	25350
<b>Box of 10</b>	





## Two-Piece Reusable Ostomy Systems

---

### Mounting Rings

Flat All-Flexible	29
Convex All-Flexible	29
Deep Convex All-Flexible	30
Extra Large All-Flexible	30
Semi-Flexible	30

### Two-Piece Drainable Pouch (Reusable)

Weight-Less™ Odour Ban™	31
Gas-Relief™ Odour Ban™	31
XTL* Odour Ban™	31
Solo™	31
White Neoprene Rubber	31

### Two-Piece Urostomy Pouches (Reusable)

All-Flexible™ Plastic	32
XTL* E-Z Drain™	32
E-Z Drain™	32
Latex Rubber	32
Solo Unit	32





## Flat All-Flexible

Accommodates stoma openings of 1/2" to 1-5/8" (1.27cm to 4.13cm) in 1/6" increments.

\*Please specify inside diameter (stoma size) when ordering.



### Basic All-Flexible


3-3/4" dia. (9.5 cm)

White Vinyl

WV-101

Green Neoprene Rubber

GN-101



### Oval Flat


3-3/4" x 2-3/4" (9.5 cm x 7.0 cm) -  
3/4" dia. (9.5 cm)

White Vinyl

WV-50

Green Neoprene Rubber

GN-50



### Compact Flat


3-3/4" x 3-1/8" (9.5 cm x 7.9 cm)

White Vinyl

WV-45

Green Neoprene Rubber

GN-45



### Petite Flat

3-3/4" x 2-7/8" (9.5 cm x 7.3 cm)

White Vinyl

WV-40


Green Neoprene Rubber

GN-40

## Convex All-Flexible

Accommodates stoma openings of 1/2" to 1-5/8" (1.27cm to 4.13cm) in 1/6" increments.

\*Please specify inside diameter (stoma size) when ordering.



### Oval Convex


3-3/4" x 2-3/4" (9.5 cm x 7.0 cm)

White Vinyl

WV-60

Green Neoprene Rubber

GN-60



### Compact Convex


3-3/4" x 3-1/8" (9.5 cm x 7.9 cm)

White Vinyl

WV-102

Green Neoprene Rubber

GN-102



### Petite Convex

3-3/4" x 2-7/8" (9.5 cm x 7.3 cm)

White Vinyl

WV-42

Green Neoprene Rubber


GN-42



## Deep Convex All-Flexible

Accommodates stoma openings of 1/2" to 1-5/8" (1.27cm to 4.13cm) in 1/6" increments.

\*Please specify inside diameter (stoma size) when ordering.



### Basic Deep Convex


3-3/4" dia. (9.5 cm)

White Vinyl

WV-80

Green Neoprene Rubber

GN-80



### Basic Deep Convex

3-3/4" dia. (9.5 cm)

White Vinyl

WV-85

Green Neoprene Rubber

GN-85

## Extra Large All-Flexible

Accommodates larger stoma openings and double-barreled stomas to 1-3/4" to 2-1/2" (4.45cm to 6.35cm) in 1/16" increments. \*Please specify inside diameter (stoma size) when ordering.



### "XTL" Oval Flat

4-1/2" x 3-1/2" (11.4 cm x 8.9 cm)

White Vinyl

XTL-55



### "XTL" Oval Convex

4-1/2" x 3-1/2" (11.4 cm x 8.9 cm)

White Vinyl

XTL-75

## Semi-Flexible

Accommodates stoma openings of 1/2" to 1-5/8" (1.27cm to 4.13cm) in 1/6" increments.

\*Please specify inside diameter (stoma size) when ordering.



### Basic Regular Convex

3-5/8" dia. (9.2 cm)

SF-10



### Compact Regular Convex

3-5/8" x 3-1/8" (9.2 cm x 7.9 cm)

SF-20



### Basic Medium Convex

4" dia. (10.2 cm)

SF-15



### Compact Medium Convex


3-5/8" x 3-1/8" (9.2 cm x 7.9 cm)

SF-65

## Two-Piece Drainable Pouch (Reusable)




\*Please specify Small, Medium, Large (size) pouch when ordering.



### Weight-Less™ Odour Ban™

Sizes: Infant, Pediatric, Small, Medium, Large


Clear	MDW-10-200
White	MDW-10-201
Beige	MDW-10-220



### Gas-Relief™ Odour Ban™

Sizes: Small, Medium, Large

Clear	GR-22-202
White	GR-22-203




### XTL\* Odour Ban™

Sizes: Medium, Large

White	XTL-MDW-20-210
-------	----------------

\*For Stoma Openings 1-3/4" to 2-1/2"



### Solo™

Sizes: Small, Medium, Large

<b>Note: When ordering specify pouch size and stoma opening</b>	SI-2001
---	---------




### White Neoprene Rubber

Small	MP-2-S
Large	MP-2-L

# Two-Piece Urostomy Pouch (Reusable)

## All-Flexible™ Plastic

Clear	MAF-12-114
White	MAF-12-115
Beige	MAF-12-222



Pouch Size	Length A	Width B	Capacity (approx.)
Infant	6-1/2" (16.5 cm)	3-1/2" (8.9 cm)	5 oz. (150 cc)
Child	17-5/8" (19.4 cm)	4-5/8" (12.4 cm)	10 oz. (300 cc)
Regular	8-5/8" (21.9 cm)	5" (12.7 cm)	15 oz. (450 cc)
Large	10" (25.4 cm)	5-3/8" (13.6 cm)	20 oz. (600 cc)

Box of 5

## XTL\* E-Z Drain™

\*For Stoma Openings 1-3/4" to 2-1/2"

Clear	XTL-EZD-24-212
-------	----------------




Pouch Size	Length A	Width B	Capacity (approx.)
Regular	8-7/8" (22.5 cm)	5" (12.7 cm)	14 oz. (420 cc)
Large	10-1/2" (26.7 cm)	5-3/8" (13.6 cm)	18 oz. (540 cc)

Box of 5

## E-Z Drain™

Clear	EZD-36-218
White	EZD-36-219
Beige	EZD-36-225




Pouch Size	Length A	Width B	Capacity (approx.)
Infant	6-1/2" (16.5 cm)	3-1/2" (8.9 cm)	5 oz. (150 cc)
Child	7-5/8" (19.4 cm)	4-5/8" (12.4 cm)	10 oz. (300 cc)
Regular	8-5/8" (21.9 cm)	5" (12.7 cm)	15 oz. (450 cc)
Large	10" (25.4 cm)	5-3/8" (13.6 cm)	20 oz. (600 cc)

Box of 5

## Latex Rubber


White / Child	MU-3R-C
White / Regular	MU-3R-R
White / Large	MU-3R-L



Pouch Size	Length A	Width B	Capacity (approx.)
Child	6-7/8" (17.5 cm)	3-1/8" (7.9 cm)	8 oz. (240 cc)
Regular	7-5/8" (19.4 cm)	5-3/8" (13.6 cm)	13 oz. (390 cc)
Large	10-3/8" (26.4 cm)	5-1/4" (13.3 cm)	24 oz. (720 cc)

## Solo Unit

Child	SU-3001-C
Regular	SU-3001-R
Large	SU-3001-L



Pouch Size	Length A	Width B	Capacity (approx.)
Child	7-5/8" (19.4 cm)	4-5/8" (12.4 cm)	10 oz. (300 cc)
Regular	8-5/8" (21.9 cm)	5" (12.7 cm)	15 oz. (450 cc)
Large	10" (25.4 cm)	5-3/8" (13.6 cm)	20 oz. (600 cc)





## Ostomy Accessories

---

Accessories	34
Urostomy Accessories	38




## Accessories



### Marlen Flange Xtender™

For people with active lifestyles or anyone concerned about added security - Box of 30

9200



### Night Drainage Kit

For use with MUPO-9, EZD-36 and Ultra Urostomy Pouches

MND-7

\*(Contains Screw Connector with Washer Seal)

\* MND-7A

Half-Gallon Plastic Container w/Special Cap

MND-7-4

Plastic Tubing (6 feet long)

MND-7-3

Rubber Tubing — (6 feet long)


MND-7-3R

\* Screw Connector with Washer Seal

310

Washer Seal Only (for Screw Connector)


315



### Drain Tube Extension

For use with EZD-36 and Ultra Urostomy Pouches

MDT-320



### Drain Tube Extension with Screw Type Fitting

For use with MAF-12, XTL-MAF-24 and MU-3R Appliances

MDT-320-A



### MaxSeal™ Protective Barrier Ring

Regular: outside diameter of 2" (50mm) - Box of 10

9700

Large: outside diameter of 3.9" (100mm) - Box of 10

9800




### UltraSeal™ Flexible Barrier Ring

Regular: outside diameter of 2" (50mm) - Box of 10

9500

Large: outside diameter of 3.9" (100mm) - Box of 5


9600



### AquaTack™ Protective Hydrocolloid Skin Barrier

4" x 4" (10 cm x 10 cm) with ½" (12 mm) starter hole - Box of 5

9100



### Skin Shield™ Protective Adhesive Skin Barrier

4" x 4" (10 cm x 10 cm) with ½" (12 mm) starter hole - Box of 5

9000



### Double-Faced Adhesive Tape Discs

3-7/8" outside diameter (9.9 cm) – Pre-cut stoma openings up to 1-5/8" (4.2 cm) - Box of 10  
*\*please specify inside diameter (stoma size) when ordering*

107



### Extra-Large Double-Faced Adhesive Tape Discs

Shaped to accommodate the Extra-Large XTL-55 and XTL-75 mounting rings - Box of 10

XTL-171



### Protex Powder (Gum Karaya)

4 oz. (120 cc) Bottle


P-116



### Small Protex Powder Pad (Karaya Seal)

Pre-cut with stoma opening sizes up to 1-1/2" (3.8 cm)  
*\*please specify inside diameter (stoma size) when ordering*


112



### Large Protex Powder Pad (Karaya Seal)

Outside diameter of 3-3/4" (9.5 cm) will accommodate openings up to 2-1/2" (6.4 cm) - Box of 5

PPC-175



### Protex Sheet

8" x 8" (20.3 cm x 20.3 cm) square


PS-176



### Plastic Foam Pads

4" diameter (10.2 cm) with pre-cut stoma openings

FP-601



### Rubber-O-Ring Seal

Pt. No. 103 is used with all mounting rings Except Extra-Large XTL Series

103

Pt. No. XTL-161 is used only with Extra-Large XTL-55 and XTL-75 Mounting Rings

XTL-161



### Lock Ring


605



### Elastic Appliance Belt (Adult)

Used with reusable mounting rings

104




### MARLEN Elastic Appliance Belt (Child)

Used with reusable mounting rings

105


## Accessories



### Ultra™ Elastic Waist Belt


Used with Ultra™ one-piece appliances for additional support.  
Individually packaged.

Regular	5004
Extra large	5004-XTL




### UltraMax™ Gemini™ Elastic Waist Belt

With the E-Z Guide™ System	29104
----------------------------	-------



### Urinary Drainage Fitting

108



### E-Z Drain Adapter

340-EZD



### Stoma Guide Strips

Pkg. of 100	603
-------------	-----



### Liquid Odour-Guard

4 oz. (120 cc) Bottle	503
-----------------------	-----



### Deodorizer

12 oz. (360 cc) Bottle	504
------------------------	-----



### MUCOSPERSE™ Mucous Dispersant

4 oz. (120 cc) Bottle	507
-----------------------	-----




### Eye Dropper Bottle

1 oz. (30 cc) Bottle	402
----------------------	-----




### Plastic Airing Hook

106



### 3M Micropore Tape

Box of 3 Rolls	612
----------------	-----



### Plastic Clamp Closure

Pouch closure for re-usable ileostomy (drainable) pouches.  
Individually packaged.


117



### Pouch Closure

One-piece pouch closure for re-usable ileostomy (drainable) pouches. After rolling up spout, sides fold and snap to secure.  
Individually packaged.


118



### Ultra™ Pouch Closure

One-piece pouch closure designed for the wider, Ultra™ drain spout. After rolling up spout, sides fold and snap to secure.  
Individually packaged.


5018



### Ultra™ Pouch Closure

Scissor-like action securely seals bottom of wider Ultra™ drain spout. Individually packaged.

5019



## Leg Bag Kit

With Regular size leg bag. For use with EZD-36 and Ultra Urostomy Pouches

MLB-300  
Complete Kit

\*(Contains Screw Connector with Washer Seal)

MLB-300A  
Complete Kit

### Kit Contents

- 1-Leg Bag
- 2-Adj. Velcro Straps
- 1-Two-Foot Tubing
- \*1-Screw Connector with Washer Seal


### Regular Size Leg Bag

Straps (2 required) 300-5  
300-10

Tubing (2 feet long) 300-25

\* Screw Connector with Washer Seal (For use with MAF-12, \* 310  
XTL-MAF-24 And MU-3R Appliances)

Washer Seal Only 315  
(for Screw Connector)



## Urinary Bedside Drain Bag Kit

For use with EZD-36 and Ultra Urostomy Pouches

MBD-8  
Complete Kit

\*(Contains Screw Connector with Washer Seal)

\* MBD-8A  
Complete Kit

### Kit Contents

- 1-Bedside Drain Bag
- 2-Plastic Hangers
- 1-Five-Foot Tubing
- \*1-Screw Connector with Washer Seal

### Bedside Drain Bag

Plastic Hangers (2 required) MBD-8-1  
MBD-8-2

Tubing (5 feet long) MBD-8-3

\* Screw Connector with Washer Seal \* 310  
(For use with MAF-12, XTL-MAF-24 And MU-3R Appliances)

Washer Seal Only (for Screw Connector) 315



# Numerical Index

---





## Numerical Index

ITEM #	PG #	ITEM #	PG #	ITEM #	PG #	ITEM #	PG #	ITEM #	PG #
103.....	<b>35</b>	21000 .....	<b>26</b>	50634 .....	<b>14</b>	52444 .....	<b>13</b>	53438 .....	<b>13</b>
104.....	<b>35</b>	21012 .....	<b>26</b>	50638 .....	<b>14</b>	52448 .....	<b>13</b>	53441 .....	<b>13</b>
105.....	<b>35</b>	21016 .....	<b>26</b>	50641 .....	<b>14</b>	52450 .....	<b>13</b>	53444 .....	<b>13</b>
106.....	<b>36</b>	21019 .....	<b>26</b>	50644 .....	<b>14</b>	52455 .....	<b>13</b>	53448 .....	<b>13</b>
107.....	<b>35</b>	21022 .....	<b>26</b>	51400 .....	<b>14</b>	52458 .....	<b>13</b>	53450 .....	<b>13</b>
108.....	<b>36</b>	21025 .....	<b>26</b>	51412 .....	<b>14</b>	52500 .....	<b>13</b>	53455 .....	<b>13</b>
112.....	<b>35</b>	21029 .....	<b>26</b>	51416 .....	<b>14</b>	52512 .....	<b>13</b>	53458 .....	<b>13</b>
117.....	<b>37</b>	21032 .....	<b>26</b>	51419 .....	<b>14</b>	52516 .....	<b>13</b>	53500 .....	<b>13</b>
118.....	<b>37</b>	21034 .....	<b>26</b>	51422 .....	<b>14</b>	52519 .....	<b>13</b>	53512 .....	<b>13</b>
300-5 .....	<b>38</b>	21038 .....	<b>26</b>	51425 .....	<b>14</b>	52522 .....	<b>13</b>	53516 .....	<b>13</b>
300-10 .....	<b>38</b>	25100 .....	<b>27</b>	51429 .....	<b>14</b>	52525 .....	<b>13</b>	53519 .....	<b>13</b>
300-25 .....	<b>38</b>	25150 .....	<b>27</b>	51432 .....	<b>14</b>	52529 .....	<b>13</b>	53522 .....	<b>13</b>
310.....	<b>34</b>	25200 .....	<b>27</b>	51434 .....	<b>14</b>	52532 .....	<b>13</b>	53525 .....	<b>13</b>
310.....	<b>38</b>	25210 .....	<b>27</b>	51438 .....	<b>14</b>	52534 .....	<b>13</b>	53529 .....	<b>13</b>
315.....	<b>34</b>	25250 .....	<b>27</b>	51441 .....	<b>14</b>	52538 .....	<b>13</b>	53532 .....	<b>13</b>
315.....	<b>38</b>	25260 .....	<b>27</b>	51444 .....	<b>14</b>	52541 .....	<b>13</b>	53534 .....	<b>13</b>
340-EZD ...	<b>36</b>	25300 .....	<b>27</b>	51512 .....	<b>14</b>	52544 .....	<b>13</b>	53538 .....	<b>13</b>
402.....	<b>36</b>	25350 .....	<b>27</b>	51516 .....	<b>14</b>	52600 .....	<b>16</b>	53541 .....	<b>13</b>
503.....	<b>36</b>	29104 .....	<b>36</b>	51519 .....	<b>14</b>	52612 .....	<b>12</b>	53544 .....	<b>13</b>
504.....	<b>36</b>	50400 .....	<b>14</b>	51522 .....	<b>14</b>	52616 .....	<b>12</b>	53600 .....	<b>16</b>
507.....	<b>36</b>	50412 .....	<b>14</b>	51525 .....	<b>14</b>	52619 .....	<b>12</b>	53612 .....	<b>12</b>
603.....	<b>36</b>	50416 .....	<b>14</b>	51529 .....	<b>14</b>	52622 .....	<b>12</b>	53616 .....	<b>12</b>
605.....	<b>35</b>	50419 .....	<b>14</b>	51532 .....	<b>14</b>	52625 .....	<b>12</b>	53619 .....	<b>12</b>
612.....	<b>36</b>	50422 .....	<b>14</b>	51534 .....	<b>14</b>	52629 .....	<b>12</b>	53622 .....	<b>12</b>
5004 .....	<b>36</b>	50425 .....	<b>14</b>	51538 .....	<b>14</b>	52632 .....	<b>12</b>	53625 .....	<b>12</b>
5004-XTL ..	<b>36</b>	50429 .....	<b>14</b>	51541 .....	<b>14</b>	52634 .....	<b>12</b>	53629 .....	<b>12</b>
5018 .....	<b>37</b>	50432 .....	<b>14</b>	51544 .....	<b>14</b>	52638 .....	<b>12</b>	53632 .....	<b>12</b>
5019 .....	<b>37</b>	50434 .....	<b>14</b>	51612 .....	<b>14</b>	52641 .....	<b>12</b>	53634 .....	<b>12</b>
9000 .....	<b>34</b>	50438 .....	<b>14</b>	51616 .....	<b>14</b>	52644 .....	<b>12</b>	53638 .....	<b>12</b>
9100 .....	<b>34</b>	50441 .....	<b>14</b>	51619 .....	<b>14</b>	52670 .....	<b>17</b>	53641 .....	<b>12</b>
9200 .....	<b>34</b>	50444 .....	<b>14</b>	51622 .....	<b>14</b>	52675 .....	<b>17</b>	53644 .....	<b>12</b>
9500 .....	<b>5</b>	50512 .....	<b>14</b>	51625 .....	<b>14</b>	52680 .....	<b>16</b>	53670 .....	<b>17</b>
9500 .....	<b>34</b>	50516 .....	<b>14</b>	51629 .....	<b>14</b>	52690 .....	<b>12</b>	53675 .....	<b>17</b>
9600 .....	<b>5</b>	50519 .....	<b>14</b>	51632 .....	<b>14</b>	52699 .....	<b>12</b>	53680 .....	<b>16</b>
9600 .....	<b>34</b>	50522 .....	<b>14</b>	51634 .....	<b>14</b>	52700 .....	<b>16</b>	53690 .....	<b>12</b>
9700 .....	<b>5</b>	50525 .....	<b>14</b>	51638 .....	<b>14</b>	52800 .....	<b>16</b>	53699 .....	<b>12</b>
9700 .....	<b>34</b>	50529 .....	<b>14</b>	51641 .....	<b>14</b>	52870 .....	<b>17</b>	53700 .....	<b>16</b>
9800 .....	<b>5</b>	50532 .....	<b>14</b>	51644 .....	<b>14</b>	52875 .....	<b>17</b>	53800 .....	<b>16</b>
9800 .....	<b>34</b>	50534 .....	<b>14</b>	52400 .....	<b>13</b>	52880 .....	<b>16</b>	53870 .....	<b>17</b>
20000 .....	<b>26</b>	50538 .....	<b>14</b>	52412 .....	<b>13</b>	52900 .....	<b>16</b>	53875 .....	<b>17</b>
20012 .....	<b>26</b>	50541 .....	<b>14</b>	52416 .....	<b>13</b>	53400 .....	<b>13</b>	53880 .....	<b>16</b>
20016 .....	<b>26</b>	50544 .....	<b>14</b>	52419 .....	<b>13</b>	53412 .....	<b>13</b>	53900 .....	<b>16</b>
20019 .....	<b>26</b>	50612 .....	<b>14</b>	52422 .....	<b>13</b>	53416 .....	<b>13</b>	54000 .....	<b>11</b>
20022 .....	<b>26</b>	50616 .....	<b>14</b>	52425 .....	<b>13</b>	53419 .....	<b>13</b>	54001 .....	<b>11</b>
20025 .....	<b>26</b>	50619 .....	<b>14</b>	52429 .....	<b>13</b>	53422 .....	<b>13</b>	54002 .....	<b>11</b>
20029 .....	<b>26</b>	50622 .....	<b>14</b>	52432 .....	<b>13</b>	53425 .....	<b>13</b>	54003 .....	<b>11</b>
20032 .....	<b>26</b>	50625 .....	<b>14</b>	52434 .....	<b>13</b>	53429 .....	<b>13</b>	54004 .....	<b>11</b>
20034 .....	<b>26</b>	50629 .....	<b>14</b>	52438 .....	<b>13</b>	53432 .....	<b>13</b>	54005 .....	<b>11</b>
20038 .....	<b>26</b>	50632 .....	<b>14</b>	52441 .....	<b>13</b>	53434 .....	<b>13</b>	54006 .....	<b>11</b>



## Numerical Index



ITEM #	PG #	ITEM #	PG #	ITEM #	PG #	ITEM #	PG #	ITEM #	PG #
54007	.....11	55006	.....15	55109	.....18	55176	.....15	55279	.....18
54008	.....11	55007	.....15	55110	.....18	55177	.....15	55280	.....18
54009	.....11	55008	.....15	55111	.....18	55178	.....15	55281	.....18
54010	.....11	55009	.....15	55112	.....18	55179	.....15	55282	.....18
54011	.....11	55010	.....15	55113	.....18	55180	.....15	55283	.....18
54012	.....11	55011	.....15	55114	.....18	55181	.....15	55284	.....18
54013	.....11	55030	.....15	55115	.....18	55182	.....15	55285	.....18
54014	.....11	55031	.....15	55116	.....18	55183	.....15	55286	.....18
54015	.....11	55032	.....15	55117	.....18	55184	.....15	55287	.....18
54016	.....11	55033	.....15	55118	.....15	55185	.....15	55288	.....18
54017	.....11	55034	.....15	55119	.....15	55186	.....15	55289	.....18
54018	.....11	55035	.....15	55120	.....15	55187	.....15	55290	.....18
54019	.....11	55036	.....15	55121	.....15	55206	.....15	55291	.....18
54020	.....11	55037	.....15	55122	.....15	55207	.....15	55292	.....18
54021	.....11	55038	.....15	55123	.....15	55208	.....15	55293	.....18
54022	.....11	55039	.....15	55124	.....15	55209	.....15	55294	.....15
54023	.....11	55040	.....15	55125	.....15	55210	.....15	55295	.....15
54024	.....11	55059	.....15	55126	.....15	55211	.....15	55296	.....15
54025	.....11	55060	.....15	55127	.....15	55212	.....15	55297	.....15
54026	.....11	55061	.....15	55128	.....15	55213	.....15	55298	.....15
54027	.....11	55062	.....15	55129	.....18	55214	.....15	55299	.....15
54028	.....11	55063	.....15	55130	.....18	55215	.....15	55300	.....15
54029	.....11	55064	.....15	55131	.....18	55216	.....15	55301	.....15
54030	.....11	55065	.....15	55132	.....18	55235	.....15	55302	.....15
54031	.....11	55066	.....15	55133	.....18	55236	.....15	55303	.....15
54032	.....11	55067	.....15	55134	.....18	55237	.....15	55304	.....15
54033	.....11	55068	.....15	55135	.....18	55238	.....15	55305	.....18
54034	.....11	55069	.....15	55136	.....18	55239	.....15	55306	.....18
54035	.....11	55088	.....15	55137	.....18	55240	.....15	55307	.....18
54036	.....11	55089	.....15	55138	.....18	55241	.....15	55308	.....18
54037	.....11	55090	.....15	55139	.....18	55242	.....15	55309	.....18
54038	.....11	55091	.....15	55140	.....18	55243	.....15	55310	.....18
54039	.....11	55092	.....15	55141	.....18	55244	.....15	55311	.....18
54040	.....11	55093	.....15	55142	.....18	55245	.....15	55312	.....18
54041	.....11	55094	.....15	55143	.....18	55264	.....15	55313	.....18
54042	.....11	55095	.....15	55144	.....18	55265	.....15	55314	.....18
54043	.....11	55096	.....15	55145	.....18	55266	.....15	55315	.....18
54044	.....11	55097	.....15	55146	.....18	55267	.....15	55316	.....18
54045	.....11	55098	.....15	55147	.....15	55268	.....15	55317	.....18
54046	.....11	55099	.....15	55148	.....15	55269	.....15	55318	.....18
54046	.....11	55100	.....18	55149	.....15	55270	.....15	55319	.....18
54048	.....11	55101	.....18	55150	.....15	55271	.....15	55320	.....18
54049	.....11	55102	.....18	55151	.....15	55272	.....15	55321	.....18
55000	.....15	55103	.....18	55152	.....15	55273	.....15	55322	.....18
55001	.....15	55104	.....18	55153	.....15	55274	.....15	55323	.....15
55002	.....15	55105	.....18	55154	.....15	55275	.....15	55324	.....15
55003	.....15	55106	.....18	55155	.....15	55276	.....18	55325	.....15
55004	.....15	55107	.....18	55156	.....15	55277	.....18	55326	.....15
55005	.....15	55108	.....18	55157	.....15	55278	.....18	55327	.....15



## Numerical Index

ITEM #	PG #	ITEM #	PG #	ITEM #	PG #	ITEM #	PG #	ITEM #	PG #
55328	15	55432	19	56400	14	57522	14	72619	23
55329	15	55433	19	56412	14	57525	14	72622	23
55330	15	55434	19	56416	14	57529	14	72625	23
55331	15	55435	19	56419	14	57532	14	72629	23
55332	15	55436	19	56422	14	57534	14	72632	23
55333	15	55437	19	56425	14	57538	14	72634	23
55352	19	55438	19	56429	14	57541	14	72638	23
55353	19	55439	19	56432	14	57544	14	72641	23
55354	19	55440	19	56434	14	57612	14	72644	23
55355	19	55441	19	56438	14	57616	14	72690	23
55356	19	55442	19	56441	14	57619	14	72875	23
55357	19	55443	19	56444	14	57622	14	73875	23
55358	19	55444	19	56512	14	57625	14	77000	21
55359	19	55445	19	56516	14	57629	14	77001	21
55360	19	55446	19	56519	14	57632	14	77002	21
55361	19	55447	19	56522	14	57634	14	77003	21
55362	19	55448	19	56525	14	57638	14	77004	21
55363	19	55449	19	56529	14	57641	14	77005	21
55364	19	55450	19	56532	14	57644	14	77006	21
55365	19	55451	19	56534	14	72400	22	77007	21
55366	19	55452	19	56538	14	72412	22	77008	21
55367	19	55453	19	56541	14	72416	22	77009	21
55368	19	55454	19	56544	14	72419	22	77010	21
55369	19	55455	19	56612	14	72422	22	77011	21
55370	19	55456	19	56616	14	72425	22	77012	24
55371	19	55457	19	56619	14	72429	22	77013	24
55372	19	55458	19	56622	14	72432	22	77014	24
55373	19	55459	19	56625	14	72434	22	77015	24
55374	19	55460	19	56629	14	72438	22	77016	24
55375	19	55461	19	56632	14	72441	22	77017	24
55376	19	55462	19	56634	14	72444	22	77018	24
55377	19	55506	18	56638	14	72448	22	77019	24
55378	19	55507	18	56641	14	72450	22	77020	24
55379	19	55508	19	56644	14	72455	22	77021	24
55380	19	55509	19	57400	14	72458	22	77022	24
55381	19	55510	19	57412	14	72500	22	77023	24
55382	19	55512	19	57416	14	72512	22	77024	24
55383	19	55513	19	57419	14	72516	22	77025	24
55384	19	55514	19	57422	14	72519	22	77026	24
55385	19	55521	19	57425	14	72522	22	77027	24
55386	19	55522	19	57429	14	72525	22	77028	24
55387	19	55523	19	57432	14	72529	22	77029	24
55388	19	55524	19	57434	14	72532	22	77030	21
55426	19	55525	19	57438	14	72534	22	77031	21
55427	19	55526	19	57441	14	72538	22	77032	21
55428	19	55533	19	57444	14	72541	22	77033	21
55429	19	55534	19	57512	14	72544	22	77034	21
55430	19	55537	19	57516	14	72612	23	77035	21
55431	19	55538	19	57519	14	72616	23	77036	21

## Numerical Index

ITEM #	PG #	ITEM #	PG #	ITEM #	PG #	ITEM #	PG #	ITEM #	PG #
77037	21	77122	21	82619	8	83641	8	SF-15	30
77038	21	77123	21	82622	8	83644	8	SF-20	30
77039	21	77124	21	82625	8	83675	7	SF-65	30
77040	21	77125	21	82629	8	83680	7	SI-2001	31
77041	24	77126	21	82632	8	83690	8	SU-3001-C	32
77042	24	77127	21	82634	8	EZD-36-218	32	SU-3001-L	32
77043	24	77128	21	82638	8	EZD-36-219	32	SU-3001-R	32
77044	24	77147	21	82641	8	EZD-36-225	32	WV-40	29
77045	24	77148	21	82644	8	FP-601	35	WV-42	29
77046	24	77149	21	82675	7	GN-40	29	WV-45	29
77047	24	77150	21	82680	7	GN-42	29	WV-50	29
77048	24	77151	21	82690	8	GN-45	29	WV-60	29
77049	24	77152	21	83400	9	GN-50	29	WV-80	30
77050	24	77153	21	83412	9	GN-60	29	WV-85	30
77051	24	77154	21	83416	9	GN-80	30	WV-101	29
77052	24	77155	21	83419	9	GN-85	30	WV-102	29
77053	24	77156	21	83422	9	GN-101	29	XTL-55	30
77054	24	77157	21	83425	9	GN-102	29	XTL-75	30
77055	24	77178	24	83429	9	GR-22-202	31	XTL-161	35
77056	24	82400	9	83432	9	GR-22-203	31	XTL-171	35
77057	24	82412	9	83434	9	MAF-12-114	32	XTL-EZD-24-212	
77058	24	82416	9	83438	9	MAF-12-115	32	.....	32
77059	21	82419	9	83441	9	MAF-12-222	32	XTL-MDW-20-210	
77060	21	82422	9	83444	9	MBD-8	38	.....	31
77061	21	82425	9	83448	9	MBD-8-1	38		
77062	21	82429	9	83450	9	MBD-8-2	38		
77063	21	82432	9	83455	9	MBD-8-3	38		
77064	21	82434	9	83458	9	MBD-8A	38		
77065	21	82438	9	83500	9	MDT-320	34		
77066	21	82441	9	83512	9	MDT-320-A	34		
77067	21	82444	9	83516	9	MDW-10-200	31		
77068	21	82448	9	83519	9	MDW-10-201	31		
77069	21	82450	9	83522	9	MDW-10-220	31		
77088	21	82455	9	83525	9	MLB-300	38		
77089	21	82458	9	83529	9	MLB-300A	38		
77090	21	82500	9	83532	9	MND-7	34		
77091	21	82512	9	83534	9	MND-7-3	34		
77092	21	82516	9	83538	9	MND-7-3R	34		
77093	21	82519	9	83541	9	MND-7-4	34		
77094	21	82522	9	83544	9	MND-7A	34		
77095	21	82525	9	83612	8	MP-2-L	31		
77096	21	82529	9	83616	8	MP-2-S	31		
77097	21	82532	9	83619	8	MU-3R-C	32		
77098	21	82534	9	83622	8	MU-3R-L	32		
77099	21	82538	9	83625	8	MU-3R-R	32		
77118	21	82541	9	83629	8	P-116	35		
77119	21	82544	9	83632	8	PPC-175	35		
77120	21	82612	8	83634	8	PS-176	35		
77121	21	82616	8	83638	8	SF-10	30		





*Security by Design*™

[www.nmswholesale.ca](http://www.nmswholesale.ca) | [info@nmswholesale.ca](mailto:info@nmswholesale.ca) | [orders@nmswholesale.ca](mailto:orders@nmswholesale.ca)

126 - 408 East Kent Avenue South, Vancouver, BC V5X 2X7

T: 604.638.2761 | 1.844.379.9101 | F: 604.879.3342

