


 **LIPOSUCTION** **FAT GRAFTING** **CONTOURING** **AUGMENTATION** **ADIPOSE TISSUE** **RECONSTRUCTION****M.D. RESOURCE**
medical device resource corporation**INNOVIA MEDICAL COMPANIES**

1.888.229.2875 | 651.789.3939
www.innoviamedical.com/shippertmedical



Innovia Medical is a **family of specialty surgical companies** that are passionately committed to partnering with medical professionals to help elevate the delivery of patient care and improve clinical outcomes.

ELEVATING CARE. IMPROVING OUTCOMES.

Innovia Medical companies develop and manufacture a portfolio of product solutions designed with both the medical professional and their patients in mind, to deliver care that:

- Saves time for both the medical professional and patient
- Delivers a consistent and predictable care experience
- Improves patient confidence with reduced recovery times and less discomfort

Innovia Medical does this through a commitment to innovation leadership, unsurpassed product quality and a dedication to provide responsive, personalized customer service.

INNOVATION LEADERSHIP.

Innovation leadership is fueled through our dedication to collaboration with medical professionals. We regularly seek feedback and input from physicians, nurses, and other medical professionals through product research, testing and development.

This allows us to truly understand the solutions medical professionals need, and guides our innovative culture focused on inventing new, industry-leading solutions and reinventing current solutions that make ordinary products extraordinary, across the specialty surgery spectrum.

UNSURPASSED PRODUCT QUALITY.

Our products deliver proven solutions and dependability which is the result of rigorous quality control and a relentless commitment to ongoing manufacturing improvements. We use only the highest quality materials and subject all products to rigorous inspection protocols. This is why we're proud to have so many customers who have continued to rely on our solutions for more than forty years.

RESPONSIVE, PERSONALIZED CUSTOMER SERVICE.

From our sales representatives to our customer service team, every associate at Innovia Medical strives to deliver an unrivaled customer experience. This begins by ensuring that each customer receives the product that is right for their specific needs. We provide assistance through technical support and live customer service representatives with extensive product knowledge, available throughout the day. We deliver quick and thorough responses, every time—from order and shipping confirmations to delivery date communications and more ensuring continuous satisfaction.

We know you're committed to providing the best patient care experience possible and that's why we commit to being your trusted partner helping you every step of the way. At Innovia Medical, we are not just products, we are part of the patient care solution.

Learn more at www.innoviamedical.com.



TABLE OF CONTENTS

ANESTHETIC INFILTRATION	18	FILTRON® ACCESSORIES	5
AQUAVAGE	6	INTERNAL NASAL SPLINTS	26
AQUAVAGE ACCESSORIES	7	LIPOSUCTION CANNULAS AND HANDLES	10
ASPIRATORS	12	NASAL PACKS	25
CANNULA CLEANING	9	NASAL SPLINTING AND PACKING GUIDE	30
CANNULA HOLE PATTERNS	11	SYRINGES AND LOCKS	17
DENVER SPLINT®	28	SEROMACATH® WOUND DRAINAGE	21
DISPOSABLE CANNULAS	20	THERMOPLAST	27
FAT GRAFTING ACCESSORIES	16	TUBING	14
FAT GRAFTING CANNULA KITS	8	VERONIQUE® COMPRESSION GARMENTS	22
FAT GRAFTING GUIDE	15	WASTE CANISTERS	12
FILTRON®	4		

To simplify procurement, Shippert Medical products are available through distributors and Group Purchasing Organization contracts.

Medline, Henry Schein, McKesson, Cardinal Health and Owens & Minor are some of the participating companies. Please call customer service with any questions.

HOW TO ORDER



Shippert Medical and MD Resource have united under Innovia Medical to offer the most comprehensive fat grafting product options in the industry.

PHONE: 1.888.229.2875 | 651.789.3939
FAX: 1.888.229.1941 | 651.789.3979
ONLINE: www.innoviamedical.com/shippertmedical
EMAIL: orders@innoviamedical.com

FILTRON® with Filter Technology

The FILTRON® does not change the normal flow of the physician's technique. Standard tumescent solution is infused. The FILTRON® unit is placed in-line between the harvest cannula and the waste canister/suction source. Suction pressure is reduced to 15-18"Hg and fat is drawn into the mesh filter within the canister. Waste (blood, epinephrine, Xylocaine) and non-viable cells are pulled through the filter and are suctioned out to the waste canister. When the desired volume is collected, the transfer tube is attached and the FILTRON® unit is inverted. The quality fat is then drawn directly into a luer lock syringe and is ready to re-inject.

FEATURES


- Closed, single-use system
- Sterile, disposable
- Separates tissue and fluids by filtration
- Filters and removes waste, non-viable cells and excess fluid during harvest
- No decant time required - fat is immediately ready for reinjection
- Compatible with most aspirators and instruments
- Eliminates need for centrifuge
- Five volumes: 100cc, 250cc, 500cc, 1000cc, 2000cc

FILTER TECHNOLOGY

FAT HARVESTED

QUALITY FAT
RETAINED IN FILTER

WASTE FILTERED AND
SUCTIONED OUT



TOP CHOICE FOR VIABLE FAT GRAFTING!

FILTRON®

Unit volume is total quality fat that is ready for reinjection.
The waste fraction is already removed.

Each Unit Includes:

- 6 ft Harvest Tubing (I.D. 5/16" and O.D. 1/2")
- Transfer Tube with luer lock tip for removal of fat

PART NUMBER	VOLUME	QTY
3-TT-FILTRON 100	100cc	1
3-TT-FILTRON 250	250cc	1
3-TT-FILTRON 500	500cc	1
3-TT-FILTRON 1000	1000cc	1
3-TT-FILTRON 2000	2000cc	1



FILTRON® Stand

- Operating room canister stand with four universal brackets
- Devices can be hung upright for harvest and upside down for fat removal
- Five legs for stability
- To maintain sterility, a sterile drape or mayo stand cover is placed over the top of the stand prior to hanging a FILTRON® device
- **3-T-STAND-TALL** - Extended height model raises the sterile field. Specially designed for use with the FILTRON®

PART NUMBER	DESCRIPTION	HEIGHT	QTY
3-T-STAND-DLX	Standard Height	34 (in)	1
3-T-STAND-TALL	Extended Height	52 (in)	1



FILTRON® Stand Add-on Bracket

- Optional bracket to attach FILTRON® units to sterile drape and O.R. stand
- Creates one more sterile barrier between stand, drape and fat collection unit
- Devices can be hung upright for harvest or upside down for fat removal
- Reusable
- Autoclavable

Standard Instructions:

1. Place sterile drape over O.R. floor stand
2. Attach sterile add-on bracket
3. Hang sterile FILTRON® unit on bracket

PART NUMBER	QTY
3-TT-BRACK-SMT	1



FILTRON® IV Pole Bracket

- Universal bracket with clamp to fit any IV pole
- Devices can be hung upright for harvest or upside down for fat removal
- Reusable
- Autoclavable

PART NUMBER	QTY
3-TT-BRACK-HK	1



AquaVage with Decantation Technology

The AquaVage does not change the normal flow of the physician's technique. Standard tumescent solution is infused. The AquaVage is placed in-line between the harvest cannula and waste canister/suction source. Suction pressure is reduced to 15-18" Hg. Lactated Ringers solution is placed in the bottom of the AquaVage canister to enable separation of the fat as it flows in. When fat hits the Lactated Ringers, it separates from the waste portion and floats to the top. Once gravity decantation process is complete, fat is removed via the extraction tube with a Luer Lock or Toomey syringe.

DECANTATION TECHNOLOGY




DURABLE
CANISTER

FEATURES

- Closed, single-use system
- Sterile, disposable
- Separates tissue and fluids by decantation
- Compatible with most aspirators and instruments
- Eliminates need for centrifuge
- Two volumes: 1200cc and 2000cc

ACENTRIC FUNNEL DESIGN



NEW
BARBED SPOUT

HOLES FOR
LIQUID DRAINAGE

aquaVage®



1200cc

2000cc

AquaVage

Each Unit Includes:

- 6' harvest tubing
- Transfer Tube with luer lock tip for removal of fat

PART NUMBER	VOLUME	QTY
AV1200	1200cc	1
AV2000	2000cc	1



AquaVage Stand

The AquaVage Stand is made completely of stainless steel allowing for reuse. It is composed of 2 detachable legs and a main plate that holds the AquaVage.

- Autoclavable

PART NUMBER	DIMENSIONS	QTY
AVSTAND	7 ¾ W x 8 ⅝ L x ½ H (in)	1



SHIPPERT MEDICAL AND MD RESOURCE HAVE UNITED UNDER INNOVIA MEDICAL TO OFFER THE MOST COMPREHENSIVE FAT GRAFTING PRODUCT OPTIONS IN THE INDUSTRY.

Fat Grafting Cannula Kits

MOST POPULAR



Starter Kit

3-CK-FGST 8 piece complete cannula kit for those starting fat grafting. Large and small volume. Includes our most popular infiltration, harvest and reinjection cannulas.

- Autoclavable

Kit Includes:

- BAN-16gaugex15cm - Infiltration LL cannula
- BAN-14gaugex30cm - Infiltration LL cannula
- 24H-3mmx30cm - Harvest cannula on Blue Aluminum handle
- 4H-3mmx30cm - Harvest cannula on Blue Aluminum handle
- 1H2-3mmx15cm - ReInjection LL cannula (body)
- 1H2-14gaugex8cm - ReInjection LL cannula (body)
- 1H2-16gaugex8cm - ReInjection LL cannula (face)
- 1H2-18gaugex8cm - ReInjection LL cannula (lips, eyelids)



Small Volume Kit

3-CK-FGSV 4 piece complete kit with infiltration, harvest, and reinjection cannulas. Includes our most common cannulas for small volume fat grafting; face, hands, and corrections.

- Autoclavable

Kit Includes:

- BAN-16gaugex15cm - Infiltration LL cannula
- 24H-3mmx15cm - Harvest LL cannula
- 1H2-16gaugex8cm - ReInjection LL cannula (face)
- 1H2-18gaugex8cm - ReInjection LL cannula (lips, eyelids)



Large Volume Kit

3-CK-FGLV 3 piece complete kit with infiltration, harvest, and reinjection cannulas. Includes our most common cannulas for large volume fat grafting; buttocks, breast, and body contouring.

- Autoclavable

Kit Includes:

- BAN-14gaugex30cm - Infiltration LL cannula
- 4H-3mmx30cm - Harvest cannula on Blue Aluminum handle
- 1H2-3mmx15cm - ReInjection LL cannula



Breast Kit

3-CK-FGB 4 piece kit for fat grafting to the breast. Includes our most common infiltration, harvest, and reinjection cannulas.

- Autoclavable

Kit Includes:

- BAN-14gaugex30cm - Infiltration LL cannula
- 4H-3mmx30cm - Harvest cannula on Blue Aluminum handle
- 1H2-3mmx15cm - ReInjection LL cannula
- 1H2-14gaugex8cm - ReInjection LL cannula

Cannula Cleaning and Sterilization

Manual Cleaning Instructions

- 1. Prepare cleaning solution per the manufacturer’s instructions
- 2. Immerse instrument into cleaning solution and individually clean using a soft bristled brush to remove all visual soil and contaminants for a minimum of two (2) minutes
- 3. Individually rinse the instruments with running cold tap water for a minimum of one (1) minute
- 4. Dry with soft, lint-free cloth and clean lumens with clean compressed air

Sterilization

Sterilize in a pre-vacuum steam sterilization cycle at 132°C (270°F) for four (4) minutes



Cannula Cleaning Brushes

Stainless steel rod with medical grade nylon bristles. Non-sterile.

PART NUMBER	SIZE	QTY
3-C-BR2	2mmx42cm (0.07 x 16.5in)	1
3-C-BR3	3mmx48cm (0.11 x 19in)	1
3-C-BR4	4mmx48cm (0.15 x 19in)	1
3-C-BR5	5mmx48cm (0.2 x 19in)	1
3-C-BR6	6mmx48cm (0.23 x 19in)	1
3-C-BR8	8mmx48cm (0.31 x 19in)	1
3-C-BR10	10mmx48cm (0.39 x 19in)	1

HOW TO ORDER CANNULAS

A multitude of hole patterns, diameters and length combinations are available. Use the example below to help build your cannula. Go to www.shippertmedical.com to order today.

Luer Lock: **3** - **C** - **BAN** - **3** x **15**
Part Number Hole Pattern Diameter Length

Handle: **3** - **H** - **BA** - **MC** - **4** x **30**
Part Number Hole Pattern Diameter Length

Diameter

18 gauge

16 gauge

14 gauge

2mm

3mm

4mm

5mm

6mm

Length

8cm

10cm

12cm

15cm

20cm

25cm

30cm

35cm

40cm

45cm

Blue Aluminum Handle with Cannula

3-H-BA-Hole Pattern-DiameterxLength

1" Blue Aluminum Handle with Cannula of choice

- Autoclavable
- Limited Availability



Luer Lock Handle

3-H-BAL

1" Blue Aluminum Handle with Luer Lock Hub

- Autoclavable



Luer Lock Cannulas

3-C-Hole Pattern-DiameterxLength

- Stainless Steel
- Autoclavable



Toomey 60cc Cannulas

3-C-T60-Hole Pattern-DiameterxLength

60cc hub fits 60cc Toomey syringes and handles

- Stainless Steel
- Autoclavable



Biplane® Handle with Cannula

3-H-SB-Hole Pattern-DiameterxLength

- Unique ergonomic design
- Reduces fatigue and repetitive stress
- Lightweight

Specially designed for FAT GRAFTING

- Cannula runs entire length of handle
- Smooth channel allows thorough cleaning with brush
- Autoclavable



Hole Patterns

INFILTRATION

BAN 6 Hole Body Spray

GAN 6 Hole Garden Spray

MH 6 Multi Hole Spray

HARVEST

4H Large Volume 4 Hole Harvester

24H Small Volume 24 Hole Harvester

REINJECTION

1H2 One Hole Type II (Coleman)

1H3 One Hole Type III (Spoon Tip)

1HV One Hole V Type

TRANSFER

OP Open-Ended

PLEASE ALLOW 5-7 BUSINESS DAYS FOR CANNULA ORDERS

**No returns on cannulas orders.*

STANDARD

MC Mercedes

TP Triport

SP Spatula

2SP 2 Port Spatula

LV Las Vegas

MAN Bucket

MS Modern Spiral

TV Toledo V

BL Blue

DS Disrupter

BB Becker Basket

BTD Becker Tear Drop



Waste Canisters

Implosion resistant waste collection canisters. No inserts, direct and easy view of measurement. Disposable.

- Both volumes are compatible with canister mount
- Liner-free
- Seal-able ports
- High flow pathogen filter
- Easy visualization of aspirate

PART NUMBER	VOLUME	QTY
LS7020	2000cc	10/box
LS7030	3000cc	10/box

MOST POPULAR!



Overflow Fat Trap

The Overflow Fat Trap prevents moisture from entering the LS2 Aspirator. Under normal operating conditions, this filter only needs to be changed during regular maintenance.

PART NUMBER	QTY
165-023.1	1



Portable Aspirator

Inexpensive suction source for those starting a new practice or doing in-office procedures.

- 1/8 horsepower pump delivers 0-26 in Hg
- 120 VAC, 60Hz, 5 Amps
- Extra strength premium suction feet attach the unit safely and securely to most surfaces
- Dual protection system with float shut-off in canister lid and in filter
- Thermal limit auto reset
- Velcro strap holds cord securely to make portability easier
- Weight 12.6 lbs. (5.7 kg)
- Length 13 1/4" (33.7 cm), Width 8" (20.3 cm), Height 9 3/4" (24.8 cm)

PART NUMBER	DESCRIPTION	QTY
3-S-PM60	Easy Vac aspirator with hospital grade power cord and one 800cc canister	1
3-S-1690	Hydrophobic In-line Filter	10/box
3-S-502519-10	800cc Disposable Canister with lid	10/box
3-S-1906-12	1200cc Disposable Canister with lid	12/box

FRONT VIEW



BACK VIEW



Touch Up Portable Aspirator

Small, quiet, and portable. Class 2, type B, power aspirator with a detachable cordset, foot-pedal control and vacuum adjustment capability. Canisters are secured to the back with the supplied canister ring.

- 1/8 HP Vacuum Pump
- Vacuum Adjustable from 0 to 27 inches of Hg
- Flow Rate 2.5 CFM free air
- Power 110 or 220 VAC
- Sound Level 38dB +/- 3dB measured 3ft. from unit front
- Enclosure Scratch free and dent resistant Kydex® case
- Dimensions: 9 W x 14 1/2 D x 10 1/2 H (in)
- Weight: 20 lbs

PART NUMBER	DESCRIPTION
TOUCHUP	Easy Vac aspirator with hospital grade power cord and one 800cc disposable canister



LS2 Aspirator

Up to 29 in Hg. 110 pounds. Graphite veins self lubricate. Sealed oil-less vacuum system, virtually maintenance-free. A taller case is provided for convenience and ease for canister monitoring. Space on "table-top" designed for equipment such as a PAL System console and/or an infiltration pump. A rear-mounted IV pole holds two canisters along with tumescent solution bags. Base with grips for easy moving. Lockable canisters. Bright digital display. Universal machine.

- Two ½ Hp Vacuum pump with 2 maintenance free, non-polluting, oil-less rotary vane pumps.
- Vacuum adjustable from 0-29" Hg
- Flow rate 4.2 SCFM
- Cooling 2 high-speed low noise fans
- Power 115 or 230 VAC 50/60 Hz
- Sound level <50 dBA
- Powder coated sheet aluminum with double sound insulation
- Chassis Steel frame with vibration - free parallel motor mounts
- Thermal protector

Includes:

- Two universal power cords
- Hospital grade power cord
- Pneumatic foot-pedal
- IV pole with canister mounts
- Storage rack
- 3 feet of tubing
- IV hooks
- Biofilter
- 2000cc disposable waste canisters - 10/box

PART NUMBER	DIMENSIONS	WEIGHT	QTY
LS2	15 ½ L x 15 ½ W x 24 H (in) 394 L x 394 W x 610 H (mm)	110 lbs (20 kg)	1



Replacement Absolute Biofilter

Replacement Absolute Biofilter. Each 99.999% DOP rated for 0.02 micron particles. Prevents potential pathogens aerosolized by your aspirator from being exhausted into the O.R. air. 15 hour replacement interval.

PART NUMBER	QTY
FL3000	1
FL3000-05	5/box

NX3 Aspirator

Super quiet piston pump aspirator. At only 45 pounds, it is lightweight yet still offers the capabilities of less mobile, heavier units. Competes with HK, WJ, Mentor standard pumps. Smaller scale in O.R. suite.

- ½ Hp Vacuum pump with maintenance free, non-polluting, dual-rocking piston
- Vacuum adjustable from 0-27" Hg
- Flow rate 4.2 SCFM
- Cooling 1 high-speed low noise fan
- Power 115 or 230 VAC 50/60 Hz
- Sound level <50 dBA
- Kydex® with double sound insulation
- Steel frame with vibration free motor mounts
- Thermal protector

Includes:

- Two universal power cords
- Hospital grade power cord
- Pneumatic foot-pedal
- IV pole with canister mounts
- Storage rack
- 3 feet of tubing
- IV hooks
- Biofilter
- 2000cc disposable waste canisters - 10/box

PART NUMBER	DIMENSIONS	WEIGHT	QTY
NX3	14 L x 16 W x 14 ½ H (in) 355 L x 1041 W x 368 H (mm)	45 lbs (20 kg)	1





Infiltration Tubing with Pump Segment

Single Spike

Auto-Fuse Infiltration tubing with a luer connector on one end and a single spike on the other end. Compatible with all peristaltic head infiltration pumps.

- Sterile
- Disposable

PART NUMBER	OUTER DIAMETER	LENGTH	QTY
AFTD	3/8" (.375")	13 feet	10/box



Double Spike

Auto-Fuse Infiltration tubing with a luer connector on one end and a double spike on the other end. Compatible with all peristaltic head infiltration pumps.

- Sterile
- Disposable

PART NUMBER	OUTER DIAMETER	LENGTH	QTY
AFTDY	3/8" (.375")	13 feet	10/box



High Pressure Aspiration Tubing - Liposuction

Compatible with any aspirator on the market. Extremely subtle with excellent drapability.

- Thick walls to withstand high pressure suction
- PVC
- Latex Free
- Sterile
- Disposable

PART NUMBER	INNER DIAMETER	OUTER DIAMETER	LENGTH	QTY
LS7010	3/8" (.375")	5/8" (.625")	8 feet	10/box
LS7012	3/8" (.375")	5/8" (.625")	10 feet	10/box



Low Pressure PVC Aspiration Tubing - Fat Grafting

- For use with Low Pressure Suction 1/2 atmosphere, 10-15" Hg, or 50% power
- Thinner walls than standard liposuction tubing
- PVC
- Latex Free
- Sterile
- Disposable

PART NUMBER	INNER DIAMETER	OUTER DIAMETER	LENGTH	QTY
SFRT10	1/4" (.25")	3/8" (.375")	10 feet	10/box
SFRT12	1/4" (.25")	3/8" (.375")	12 feet	10/box



Silicone Tubing

A 50 foot very soft formulated silicone. This tubing is designed for the K Pump or any peristaltic pump. Non-sterile. Intended to be cut into smaller sections for repetitive injection.

- Autoclavable

PART NUMBER	INNER DIAMETER	OUTER DIAMETER	LENGTH	QTY
AFTR-BULK	.188"	3/8" (.375")	50 feet	1



PAL Microaire® Tapered Aspiration Tubing

- Individual
- Sterile
- 12' tapered tubing
- Disposable

PART NUMBER	DESCRIPTION	LENGTH	QTY
PAL-1200	Tubing	12 feet	10/box
3-PAL-CPLR	Luer Adapter - Allows use of LL cannula with a PAL hand piece	N/A	1

FAT GRAFTING GUIDE

			<100cc	Eyelids	Lips	Face	Body (Small)	Body (Large)	Breast	Buttocks
Infusion Cannulas	3-C-BAN-16x15		Luer Lock 6 Hole Spray 16 gaugex15cm	•	•	•	•			
	3-C-GAN-16x15									
	3-C-BAN-14x30						•	•	•	•
	3-C-GAN-14x30									
Harvest Cannulas	3-C-24H-3x15		24 hole 3mmx15cm	•	•	•	•			
	3-H-BA-24H-3x15 (handle)									
	3-H-BA-4H-3x30 (handle)		24 hole 3mmx30cm				•	•	•	•
Collection Units	3-S-V60		Manual VacLok® Syringe	•	•	•				
	3-TT-FILTRON 100		100cc		•	•	•			
	3-TT-FILTRON 250		250cc				•			
	3-TT-FILTRON 500		500cc				•	•	•	
	3-TT-FILTRON 1000		1000cc					•	•	
	AV1200		1200cc					•	•	
	AV2000		2000cc					•	•	•
	3-TT-FILTRON 2000		2000cc					•	•	•
Injection Cannulas	3-C-1H2-18x8		Luer Lock 1 Hole, Type II 18 gaugex8cm	•	•	•				
	3-C-1H2-16x8		Luer Lock 1 Hole, Type II 16 gaugex8cm	•			•			
	3-C-1H2-14x8		Luer Lock 1 Hole, Type II 14 gaugex8cm				•		•	
	3-C-1H2-3x15		Luer Lock 1 Hole, Type II 3mmx15cm					•	•	•

This is a guide only. Treatment plan is determined by healthcare professional.



VacLok® Syringes For Manual Harvest

Allow clinicians to lock the plunger in multiple positions for a negative pressure vacuum. This feature minimizes hand fatigue while increasing procedural efficiencies

- For manual fat harvest without outside suction source
- Luer Lock tip
- Built-in locking mechanism to maintain negative pressure at small increments
- Polycarbonate barrel for clear visibility
- ABS/PC/Silicone
- Sterile
- Disposable

PART NUMBER	VOLUME	QTY
3-S-V10	10cc	25/box
3-S-V20	20cc	20/box
3-S-V30	30cc	20/box
3-S-V60	60cc	25/box



Disposable Luer to Luer Transfer

- Female Luer Lock to Luer Lock transfer system
- Large finger grips protect against contamination
- Polypropylene
- Sterile
- Disposable
- Used to transfer tissue or emulsify between syringes

PART NUMBER	QTY
3-S-13901	50/box



Reusable Luer to Luer Transfer

- Female Luer Lock to Luer Lock transfer system
- Used to transfer tissue or emulsify between syringes
- Stainless Steel
- Autoclavable


PART NUMBER	SIZE	QTY
3-T-TLL3	3mm	1
3-T-TLL4	4mm	1



Reusable Toomey to Luer Lock Syringe Adapter

- 60cc Toomey Monoject syringe hub to Luer Lock Syringe
- Stainless Steel
- Used to transfer tissue between syringes
- Autoclavable

PART NUMBER	QTY
3-S-T60-LL-SS	1



BD Luer-Lok™ Syringes

- BD Syringes
- Luer Lock
- Sterile

PART NUMBER	VOLUME	QTY
3-S-LL1-20	1cc	20/box
3-S-LL3-20	3cc	20/box
3-S-LL5-20	5cc	20/box
3-S-LL10-20	10cc	20/box
3-S-LL20-20	20cc	20/box
3-S-LL30-20	30cc	20/box
3-S-LL60-20	60cc	20/box

Monoject™ Syringes

PART NUMBER	DESCRIPTION	VOLUME	QTY
3-S-T10-20	Slip Tip, non-sterile	10cc	20/box
3-S-T35-20	Irrigation Tip, sterile	35cc	20/box
3-S-T60-20	Toomey Tip, sterile	60cc	20/box
3-S-MLL60-20	Luer Lock Tip, sterile	60cc	20/box

Syringe Racks

Plastic Syringe Racks are made of a special polypropylene, steam autoclavable up to 121°C (250°F), and will be able to stay submerged and maintain stability in water. Easy to assemble. Stackable.

PART NUMBER	SIZE	HOLE	HOLES DIAMETER	FOR SYRINGE	QTY
3-S-R5	5cc	90	13mm	5cc, 3cc & 1 cc BD	1
3-S-R10	10cc	60	16mm	10cc BD	1
3-S-R30	30cc	24	25mm	30cc BD, 20cc Monoject	1
3-S-R60	60cc	21	30mm	60cc BD, 35cc & 60cc Monoject	1

Dimension of assembled units for all specifications:

24.6 L x 10.4 W x 6.4 H (mm)

9.7 L x 4.09 W x 2.52 H (in)

Luer Lock Syringe Tip Caps

Red polypropylene luer lock cap for standard luer lock syringes.

- Packaged individually
- Sterile
- Disposable

PART NUMBER	QTY
3-S-65001	100/box

Horseshoe and Clicker Locks

Locks are used to maintain negative pressure in syringes during manual tissue aspiration.

Horseshoe

- Aluminum
- Autoclavable

PART NUMBER	DESCRIPTION	FOR SYRINGE	QTY
3-S-HL10M	10cc Horseshoe/Twist Lock	10cc Monoject Syringe	1
3-S-HL10	10cc Horseshoe/Twist Lock	10cc BD Syringe	1
3-S-HL60	60cc Horseshoe/Twist Lock	60cc BD or Monoject Syringe	1

Clicker

- Aluminum with stainless steel spring
- Autoclavable

PART NUMBER	DESCRIPTION	FOR SYRINGE	QTY
3-S-CL10M	10cc Clicker/Snapper Lock	10cc Monoject Syringe	1
3-S-CL10	10cc Clicker/Snapper Lock	10cc BD Syringe	1
3-S-CL60	60cc Clicker/Snapper Lock	60cc BD and Monoject Syringe	1



Horseshoe



Clicker



Infiltration K-Pump

Quiet, efficient, vibration free peristaltic infiltration pump. The control knob allows setting from 0 to 1000ml per minute flow rate. Infiltration tubing (AFTD/AFTDY) can be easily mounted with the push of a lever that lets you insert silicone tubing section quickly and accurately, and locking in place. Dual pneumatic foot-pedal capability, a three roller peristaltic head; a safety lever and additional power connector for a Tumescant Measuring Device (TMD). Pump is available in either latching or non-latching operation.

PART NUMBER	DESCRIPTION	DIMENSIONS	WEIGHT	QTY
KPUMP	Latching	13 L x 10 W x 9 H (in) 330 L x 254 W x 229 H (mm)	9 lbs	1
KPUMP-NL	Non-Latching	13 L x 10 W x 9 H (in) 330 L x 254 W x 229 H (mm)	9 lbs	1



Tumescant Measuring Device (TMD)

Its convenient pole top mount delivers effortless measurement of total tumescent delivery. Weight lost from the hanging tumescent fluid bag is digitally displayed with +/- 1 ml accuracy. The very low current circuitry enables the TMD to be "piggy backed" to our K Pump or it can be used as a "stand-alone" device.

PART NUMBER	QTY
TMD	1



Replacement Pancake Foot Pedal

Standard pneumatic, pancake foot-pedal with 10 ft hose.

PART NUMBER	QTY
AS4500-LOW	1



Replacement K-Pump Power Supply

This Power Supply is used for both the TMD and the K-Pump. Only one power supply is needed if you have both. Instead you can purchase an umbilical cord that matches your set up to power both TMD and the K-Pump with one power supply.

PART NUMBER	QTY
115-003	1



3-H-IHSL



3-H-IHSS



3-H-BALI

Infiltration Handles

3-H-IHSL Large 1" Blue Aluminum Infiltration Handle
Stainless Steel Luer Lock fittings on both ends.
• Autoclavable

3-H-IHSS Small ¾" Blue Aluminum Infiltration Handle
Stainless Steel Luer Lock Fittings on both ends.
• Autoclavable

3-H-BALI Shippert Black Aluminum Infiltration Handle, ¾" diameter
On/Off Push Button and Stainless Steel Luer Lock fittings on both ends.
• Autoclavable



3-C-GAN



3-C-BAN



3-C-MH

Tumescent Infiltration Cannulas


Luer Lock Cannulas for infusion/infiltration. Select diameter and length.
• Stainless Steel
• Autoclavable

PART NUMBER	DESCRIPTION
3-C-GAN	6 Hole Garden Spray
3-C-BAN	6 Hole Body Spiral
3-C-MH	6 Multi Hole Spray

Most Common Diameters/Lengths:

3-C-holepattern-14gaugex30cm (body/large areas)

3-C-holepattern-16gaugex15cm (face/small areas)



HK Monty Infiltration Cannulas Sets

Luer Lock Cannulas Sets for infusion/infiltration. Stainless Steel.
• Autoclavable

Each Set Includes:

- 2 Full Monty - holes along 90% of cannula length, from the distal end
- 2 Half Monty - holes along 50% of cannula length, from the distal end
- 1 Tip Monty - holes along 20% of cannula length, from the distal end

PART NUMBER	SIZE	QTY
HK-IC14MSet	14 gaugex20cm	5 cannulas
HK-IC16MSet	16 gaugex15cm	5 cannulas
HK-IC18MSet	18 gaugex15cm	5 cannulas

Sterile Disposable Luer Lock Cannulas

- Sterile
- Disposable
- One patient use
- Reduce risk of spreading infectious disease

- No misplacing instruments at multiple O.R.s
- No large instrument investment
- Lipocyte friendly hole patterns

- Reduce cross-contamination
- No cleaning time or labor
- No autoclave



Sterile Disposable Luer Lock Cannulas

Hub with standard Luer Lock connection. Sterile. 10/box.

4 Hole
3-CD-4H-3x15

Description
4 Hole Harvester, 3mmx15cm

24 Hole
3-CD-24H-3x15

Description
24 Hole Harvester, 3mmx15cm

One Hole Type II
3-CD-1H2-3x15
3-CD-1H2-16x8
3-CD-1H2-18x9

Description
1 Hole II Injector, 3mmx15cm
1 Hole II Injector, 16gaugex8cm
1 Hole II Injector, 18gaugex9cm

Open-Ended
3-CD-OP-3x20

OP, 3mmx20cm
To remove fat from collection devices

Sterile Disposable - Handle with Cannula


Lightweight white handle with thumb groove. Sterile. 6/box.

4 Hole
3-H-D-4H-3x13 4 Hole, 3mmx13cm
3-H-D-4H-3x28 4 Hole, 3mmx28cm
3-H-D-4H-4x13 4 Hole, 4mmx13cm
3-H-D-4H-4x28 4 Hole, 4mmx28cm

24 Hole
3-H-D-24H-3x13 24 Hole, 3mmx13cm
3-H-D-24H-3x28 24 Hole, 3mmx28cm
3-H-D-24H-4x13 24 Hole, 4mmx13cm
3-H-D-24H-4x28 24 Hole, 4mmx28cm

Triport
3-H-D-TP-4x13 Triport, 4mmx13cm
3-H-D-TP-4x28 Triport, 4mmx28cm

Mercedes
3-H-D-MC-3x23 Mercedes, 3mmx23cm
3-H-D-MC-4x28 Mercedes, 4mmx28cm



Sterile Disposable Micro Re-injection Cannula

- Sterile
- Blunt atraumatic tip
- One side port hole pattern

- Luer Lock hub
- 5 year shelf life

- For fat or filler injection
- 20/box

PART NUMBER	DESCRIPTION
3-CD-1811	18 gaugex11cm
3-CD-187	18 gaugex7cm
3-CD-217	21 gaugex7cm**
3-CD-229	22 gaugex9cm
3-CD-227	22 gaugex7cm

**most popular

PART NUMBER	DESCRIPTION
3-CD-225	22 gaugex5cm
3-CD-237	23 gaugex7cm
3-CD-235	23 gaugex5cm
3-CD-245	24 gaugex5cm
3-CD-2438	24 gaugex3.8cm

PART NUMBER	DESCRIPTION
3-CD-255	25 gaugex5cm**
3-CD-265	26 gaugex5cm
3-CD-2638	26 gaugex3.8cm
3-CD-2625	26 gaugex2.5cm
3-CD-274	27 gaugex4cm**



SeromaCath® Wound Drainage Kits

Sterile wound drainage systems for post-operative seromas following mastectomy, axillary dissection, groin dissection, latissimus dorsi flap, TRAM flap, liposuction, abdominoplasty, breast implants, and other procedures.

- Connecting drainage tube for a bulb suction reservoir
- Bulb sold separately
- Sterile

PART NUMBER

SIZE

QTY

GR-1275

2.75" Blunt for Seromas surrounding Breast implants
2" (5cm) 14G Teflon catheter/stylet needle/plug with two side holes
2.75" (7cm) 17G blunt metal catheter with two side holes
Connecting drainage tube for a bulb suction reservoir
GR-2001 Bulb sold separately

5 kits/
box

GR-1035

3.5" for Deep Seromas common with abdominoplasty
3.5" (8.9cm) 14G Teflon catheter/stylet needle/plug with six side holes
Connecting drainage tube for a bulb suction reservoir
GR-2001 Bulb sold separately

5 kits/
box

GR-1002

2" for Common Seromas
2" (5cm) 14G Teflon catheter/stylet needle/plug with six side holes
Connecting drainage tube for a bulb suction reservoir
GR-2001 Bulb sold separately

5 kits/
box

**SAPPHIRE WOUND DRAINAGE BULB
MUST BE PURCHASED SEPARATELY.**



GR-2001

SeromaCath® Sapphire Wound Drainage Bulb

- High silicone content gives it the highest suction of any commercially available 100cc silicone reservoir
- Offers 55% more suction than the popular "JP Bulb"
- For use in general, orthopedic, urologic and neurosurgery
- Markers at 25, 50, 75 and 100cc for fast and easy fluid measurement
- Anti-reflux valve helps prevent back flow
- A strap with stainless steel clip makes attaching the Sapphire to a gown or clothing convenient and discreet
- Sterile

PART NUMBER

SIZE

QTY

GR-2001

100ml

1 bulb


**NECESSARY FOR USE WITH SEROMACATH®
WOUND DRAINAGE KITS ABOVE.**



3" Breast Wrap

V-510-3


- Adjustable downward compression
- Adjustable front hook & loop closures
- Nylon/Spandex
- White
- One size



6" Breast Wrap

V-510-6


- Adjustable downward compression
- Adjustable front hook & loop closures
- Nylon/Spandex
- White
- One size



9" Breast Wrap

V-510-9

- Adjustable downward compression
- Adjustable front hook & loop closures
- Nylon/Spandex
- White
- One size



Vianna Lightweight Support Bra

V-458


- Seamless molded contouring cups
- Front hook & eye closure
- Smooth, seamless back
- Nylon/Spandex
- White, champagne or black
- XS-3XL



6" Single Panel Binder

V-350-6

- Adjustable front hook & loop closure
- Reinforced back stays
- Strong, durable stretch
- White
- XS-4XL



9" Double Panel Binder

V-350-9


- Two overlapping 6" wide panels form 9" wide binder
- Adjustable front hook & loop closure
- Reinforced back stays
- Strong, durable stretch
- White
- XS-4XL



3" Breast Wrap with Lace

V-510J

- Adjustable downward compression
- Adjustable front hook & loop closures
- Elegant front lace trim
- Nylon/Spandex
- White or black
- One size



Versatile Stabilizing Breast Wrap

V-510-X

- Customized band positioning
- 2" wide customizing bands
- 3" wide anchoring strap
- Adjustable front hook & loop closures
- Nylon/Spandex
- White
- One size



Universal Facial Band

V-210

- Adjustable/detachable anchor straps
- Soft, durable stretch
- Hook & loop closures
- White or champagne
- One size



Universal Facial Band with Cold Hot Compresses

V-210K

- Detachable cotton knit pouches with cold/hot compresses
- Soft, durable stretch
- Adjustable/detachable anchor straps
- Hook & loop closures
- White
- One size



Universal Facial Band

V-210N

- Adjustable contouring face and neck bands
- Soft, durable stretch
- Hook & loop closures
- White or champagne
- Small-XL



Occipital Universal Facial Band

V-210OB

- Adjustable/detachable anchor straps
- Soft, durable stretch
- Hook & loop closures
- White or champagne
- One size



Universal Facial Band with Cotton Lining

V-210U

- Adjustable/detachable anchor straps
- Soft cotton lining
- Soft, durable stretch
- Hook & loop closures
- White or champagne
- One size



Universal Split Ear Facial Band

V-210S

- Adjustable/detachable anchor straps
- Split ear design
- Soft, durable stretch
- Hook & loop closures
- White or champagne
- One size

To maintain sanitary requirements:

NO RETURNS ON GARMENT ORDERS

Garments are not guaranteed in stock. Please order one week prior to procedure.

All Veronique Garments are **LATEX-FREE**.

VERONIQUE FEMALE SIZE CHART

SIZE	WAIST	HIP	BUST	BREAST	UNDER BICEP CIRCUMFERENCE
X-Small	21-31"	31-33"	27-29"	25-26"	9.5-10.25"
Small	24-26"	34-36"	30-32"	27-30"	10.5-11.25"
Medium	27-29"	37-39"	33-36"	31-34"	11.5-12.25"
Large	30-32"	40-42"	37-40"	35-38"	12.5-13.25"
XL	33-35"	43-45"	41-44"	39-42"	13.5-14.25"
2XL	36-39"	45-48"	45-48"	43-46"	14.5-15.25"
3XL	40-43"	49-52"	49-52"	47-49"	15.75-16.75"
4XL	44-47"	53-56"	53-56"	50-52"	17.00-18.00"

VERONIQUE MALE SIZE CHART

SIZE	WAIST	HIP	CHEST
Small	31-33"	36-38"	35-37"
Medium	34-36"	39-41"	38-41"
Large	37-39"	42-44"	42-45"
XL	40-43"	45-47"	46-49"
2XL	44-48"	48-51"	50-53"
3XL	49-53"	52-55"	54-58"
4XL	54-56"	56-59"	59-62"



Zippered High Back Torso Brief with Bra

V-B85

- Front zipper with cotton guard tape
- Soft molded cotton knit bra without underwire
- Targeted compression panels
- Adjustable hook & eye shoulder straps and crotch closure
- Breathable powernet fabric retains compressive memory
- White, champagne or black
- XS-4XL



Zippered High Back Girdle with Bra-Above Knee

V-B853

- Front zipper with cotton guard tape
- Soft molded cotton knit bra without underwire
- Targeted compression panels
- Adjustable hook & eye shoulder straps
- Reinforced open crotch
- Breathable powernet fabric retains compressive memory
- White, champagne or black
- XS-4XL



Zippered High Back Girdle with Bra-Below Knee

V-B855

- Front zipper with cotton guard tape
- Soft molded cotton knit bra without underwire
- Targeted compression panels
- Adjustable hook & eye shoulder straps
- Reinforced open crotch
- Gripper lace on leg
- Breathable powernet fabric retains compressive memory
- White, champagne or black
- XS-4XL



Zippered Above Knee Girdle

V-753

- Side separating zippers with cotton guard tape
- Abdominal compression panels
- 3"-wide waistband
- Gripper lace on leg
- Reinforced open-crotch
- Breathable powernet fabric retains compressive memory
- White, champagne or black
- XS-4XL



Zippered Abdominal Girdle

V-850

- Side zippers with cotton guard tape
- Abdominal compression panels
- Adjustable detachable suspenders
- Gripper lace on leg
- Reinforced open-crotch
- Breathable powernet fabric retains compressive memory
- White, champagne or black
- XS-4XL



Zippered Body Girdle

V-853

- Side zippers with cotton guard tape
- Abdominal compression panels
- Adjustable detachable suspenders
- Gripper lace on leg
- Reinforced open-crotch
- Breathable powernet fabric retains compressive memory
- White, champagne or black
- XS-4XL



Adjustable Compression Vest

V-640


- Adjustable front hook & eye closure
- Wide, adjustable hook & loop shoulder straps
- 2" wide waist band
- Breathable powernet fabric retains compressive memory
- White or black
- S-4XL




Zippered Abdominal Garment

V-840


- Side zippers with cotton guard tape
- Targeted compression panels in key surgical sites
- Wide, adjustable hook & loop shoulder straps
- Reinforced open-crotch
- Breathable powernet fabric retains compressive memory
- White or black
- S-4XL



Compressed



Hydrated



Anatomical with Airway

AbsorbENT™ Nasal Packs

Nasal and sinus packs utilize three-dimensional interconnected open cells that are soft and pliable for epistaxis management. Sterile Polyvinyl Alcohol (PVA) foam. Single use.

- String attached for retention and removal
- Soft, lint and fiber free
- Gentle, but firm compression to nasal mucosa
- Rapid wicking properties. Even and fast expansion

Note: Dimensions listed are for hydrated state.

PART NUMBER	SIZE	DIMENSIONS	QTY
RH-2315-20	Small	5 ½ L x 2 W x 1 ½ H (cm)	20
RH-7401-10	Small	4 ½ L x 2 W x 1 ½ H (cm)	10
RH-7408-10	Small with Airway	4 ½ L x 2 W x 1 ½ H (cm)	10
RH-7404-10	Large	8 L x 2 W x 1 ½ H (cm)	10
RH-7405-10	Large with Airway	8 L x 2 W x 1 ½ H (cm)	10
RH-7409-10	Large Anatomical	8 L x 3 W x 1 ½ H (cm)	10
RH-7410-10	Large Anatomical with Airway	8 L x 3 W x 1 ½ H (cm)	10
RH-7411-10	Large Anatomical - Thin	8 L x 3 W x 1 H (cm)	10
RH-7406-10	Large - Thin	8 L x 2 W x 1 H (cm)	10
RH-2310-10	Extra Large	10 L x 2 ½ W x 1 ½ H (cm)	10
RH-2312-10	Extra Large without String	10 L x 2 ½ W x 1 ½ H (cm)	10

11S-NS300-08AS	Small	3 L x 1 W x 2 H (cm)	8
11S-NS500-08AS	Medium	5 L x 1 W x 2 H (cm)	8
11S-NS800-08AS	Large	8 L x 1 W x 2 H (cm)	8
11S-NS1000-08AS	Extra Large	10 L x 1 W x 2 H (cm)	8

Rhino Rocket® Nasal Packs with Applicator

- Sterile PVA Expandacell® foam
- Epistaxis management after trauma, nasal or sinus surgery
- Gentle, but firm mucosal compression
- Tamponade of anterior or middle segment epistaxis
- Approximately 30 seconds to position and insert
- When in contact with moisture, foam swells to six times its compressed diameter
- String attached for retention and easy removal
- Dimensions listed are for hydrated state

Slimline

PART NUMBER	SIZE	DIMENSIONS	QTY
11S-S0300-08AS	Small	3 L x 2 W x 1 H (cm)	8
11S-S0500-08AS	Medium	5 L x 2 W x 1 H (cm)	8
11S-S0800-08AS	Large	8 L x 2 W x 1 H (cm)	8
11S-S0800-4PAS	Large Pairs	8 L x 2 W x 1 H (cm)	4 Pairs
11S-S1000-08AS	Extra Large	10 L x 2 W x 1 H (cm)	8

Wideline

11W-W0300-08AS	Small	3 L x 2 W x 1 ½ H (cm)	8
11W-W0500-08AS	Medium	5 L x 2 W x 1 ½ H (cm)	8
11W-W0800-08AS	Large	8 L x 2 W x 1 ½ H (cm)	8
11W-W1000-08AS	Extra Large	10 L x 2 W x 1 ½ H (cm)	8

Starter Pack

11S-STARTERRHINO	Variety pack for new users. One of each size.	8/box
-------------------------	--	-------



PATENTED APPLICATOR

- Minimizes blood containment concerns
- Protects foam from moisture
- Aids in insertion




Breeze™ Nasal Airway Splint

Unique airway geometry creates a pressure gradient force that increases airflow. Suture indicators are located at anterior tip to aid suture placement. Lubricious coating aids insertion and removal. Sterile, single-use.

- Increase airflow
- Soft medical-grade silicone
- Commonly referred to as "Doyle Splint"

PART NUMBER	DESCRIPTION	QTY
SP-79300	With Airway, Sterile	1 Pair
SP-79300-05	With Airway, Sterile	5 Pairs
SP-79300-50	With Airway, Sterile	50 Pairs
SP-79200-B	With Airway, Non-Sterile	50 Pairs
SP-79400	Without Airway, Sterile	1 Pair



Nasal Septal Button

One-piece design with 2 flange leaves for maximum security and cleanliness. Designed to reduce or eliminate many symptoms of nasal septal perforations. Easily trimmed to ensure proper fit. Medical-grade silicone. 6mm posts on all sizes. Sterile, single-use.

PART NUMBER	DESCRIPTION	DIMENSIONS	QTY
SP-78100	Small	3cm x 4mm	1
SP-78105	Medium	5cm x 4mm	1
SP-78107	Large	7cm x 4mm	1

***Pre-marked in 1cm increments for easy trimming. Only applies to large size.*



Bi-Valve Nasal Splint

Non-stick surface reduces adhesion and/or hematoma. Medial slit facilitates insertion and removal. Holes for easier suturing. Sterile, single use. Clear fluoroplastic, medical-grade material.

PART NUMBER	DESCRIPTION	QTY
SP-27000	Standard Thin (.25mm)	6 Pairs/box
SP-27005	Standard Thick (.5mm)	6 Pairs/box
SP-27010	Oversized Thin (.25mm)	6 Pairs/box
SP-27015	Oversized Thick (.5mm)	6 Pairs/box

Thermoplast Splints

- Rigid, opaque, lightweight material
- Soften with hot water and mold over taped nose

Options Available:


- Pre-cut Denver Splint® sizes
- 3x3" Flats
- 9x12" Sheets to cut your own shape



Pre-Cut Denver Splint® Series 4000

PART NUMBER	DESCRIPTION	QTY
10-4021-05K	Ivory 1/16"	Kit 5/box
10-4021-10NK	Ivory 1/16"	Non-Kit 10/box
10-4031-05K	Ivory 3/32"	Kit 5/box
10-4031-10NK	Ivory 3/32"	Non-Kit 10/box
10-4022-05K	Blush 1/16"	Kit 5/box
10-4022-10NK	Blush 1/16"	Non-Kit 10/box
10-4032-05K	Blush 3/32"	Kit 5/box
10-4032-10NK	Blush 3/32"	Non-Kit 10/box

INDICATE SIZE ON SERIES 4000



Kit Includes:

- One Denver Splint®
- Two Pre-cut Protecto Strips, sterile
- One alcohol wipe, sterile
- One skin prep, sterile
- One dorsal pad



3x3" Squares

PART NUMBER	DESCRIPTION	QTY
AQ-F-05-21	Ivory 1/16"	5/box
AQ-F-10-21	Ivory 1/16"	10/box
AQ-F-05-31	Ivory 3/32"	5/box
AQ-F-10-31	Ivory 3/32"	10/box
AQ-F-05-22	Blush 1/16"	5/box
AQ-F-10-22	Blush 1/16"	10/box
AQ-F-05-32	Blush 3/32"	5/box
AQ-F-10-32	Blush 3/32"	10/box

9x12" Sheets

PART NUMBER	DESCRIPTION	QTY
AQ-9x12-21	Ivory 1/16"	1
AQ-9x12-31	Ivory 3/32"	1
AQ-9x12-22	Blush 1/16"	1
AQ-9x12-32	Blush 3/32"	1

Denver Nasal Splint®

- Stabilization of the nose after surgery
- Reduction of edema after nasal surgery or trauma
- Moldable aluminum retains nasal bones in position
- Layers give cushioned pressure on subcutaneous tissue
- Dorsal pad reduces dead space and possibility of polybeak
- Four sizes: Mini, Petite, Small/Medium, Large
- Kit or Non-Kit



FEATURES

- Adhesive
- Moldable aluminum
- No water mess
- Approximately 60 seconds to apply
- Inconspicuous
- Water resistant

Denver Splint® - One Piece

Denver Splint® Series 2000

- Soft moldable aluminum and 1/8" ventilating foam
- Length of wear 5-7 days


PART NUMBER	DESCRIPTION	DIMENSIONS	QTY
10-2000-05KM	Mini	50mm x 28mm	Kit 5/box
10-2000-05KP	Petite	50mm x 30mm	Kit 5/box
10-2000-05KS	Small/Medium	65mm x 35mm	Kit 5/box
10-2000-05KL	Large	69mm x 43mm	Kit 5/box
10-2000-10NKM	Mini	50mm x 28mm	Non-kit 10/box
10-2000-10NKP	Petite	50mm x 30mm	Non-kit 10/box
10-2000-10NKS	Small/Medium	65mm x 35mm	Non-kit 10/box
10-2000-10NKL	Large	69mm x 43mm	Non-kit 10/box

Denver Splint® Series 3000 - No Adhesive

- Soft moldable aluminum with a thin foam pad
- Similar to "Brown Nasal Splint"

PART NUMBER	DESCRIPTION	DIMENSIONS	QTY
10-3000-05KM	Mini	50mm x 28mm	Kit 5/box
10-3000-05KP	Petite	50mm x 30mm	Kit 5/box
10-3000-05KS	Small/Medium	65mm x 35mm	Kit 5/box
10-3000-05KL	Large	69mm x 43mm	Kit 5/box
10-3000-10NKM	Mini	50mm x 28mm	Non-kit 10/box
10-3000-10NKP	Petite	50mm x 30mm	Non-kit 10/box
10-3000-10NKS	Small/Medium	65mm x 35mm	Non-kit 10/box
10-3000-10NKL	Large	69mm x 43mm	Non-kit 10/box

Denver Splint® - Two Piece



Kit Includes:

- One Denver Splint®
- Two Pre-cut Protecto Strips, sterile
- One alcohol wipe, sterile
- One skin prep, sterile
- One dorsal pad

Splinting Steps - 2 Peice



1. Alcohol & Skin Prep



2. Protecto Tape



3. Dorsal Pad



4. Adhesive Velcro



5. Aluminum Splint

Denver Splint® Series 1500

- Soft moldable aluminum and Velcro™
- Length of wear 7-10 days

PART NUMBER	DESCRIPTION	DIMENSIONS	QTY
10-1500-05KM	Mini	50mm x 28mm	Kit 5/box
10-1500-05KP	Petite	50mm x 30mm	Kit 5/box
10-1500-05KS	Small/Medium	65mm x 35mm	Kit 5/box
10-1500-05KL	Large	69mm x 43mm	Kit 5/box
10-1500-10NKM	Mini	50mm x 28mm	Non-kit 10/box
10-1500-10NKP	Petite	50mm x 30mm	Non-kit 10/box
10-1500-10KNS	Small/Medium	65mm x 35mm	Non-kit 10/box
10-1500-10NKL	Large	69mm x 43mm	Non-kit 10/box

Denver Splint® Series 5000

- Soft moldable aluminum and Velcro™ with attached dorsal pad
- Length of wear 7-10 days

PART NUMBER	DESCRIPTION	DIMENSIONS	QTY
10-5000-05KM	Mini	50mm x 28mm	Kit 5/box
10-5000-05KP	Petite	50mm x 30mm	Kit 5/box
10-5000-05KS	Small/Medium	65mm x 35mm	Kit 5/box
10-5000-05KL	Large	69mm x 43mm	Kit 5/box
10-5000-10NKM	Mini	50mm x 28mm	Non-kit 10/box
10-5000-10NKP	Petite	50mm x 30mm	Non-kit 10/box
10-5000-10NKS	Small/Medium	65mm x 35mm	Non-kit 10/box
10-5000-10NKL	Large	69mm x 43mm	Non-kit 10/box

Denver Splint® Series 1800

- Soft moldable aluminum and Velfoam
- Length of wear 7-10 days

PART NUMBER	DESCRIPTION	DIMENSIONS	QTY
10-1800-05KM	Mini	50mm x 28mm	Kit 5/box
10-1800-05KP	Petite	50mm x 30mm	Kit 5/box
10-1800-05KS	Small/Medium	65mm x 35mm	Kit 5/box
10-1800-05KL	Large	69mm x 43mm	Kit 5/box
10-1800-10NKM	Mini	50mm x 28mm	Non-kit 10/box
10-1800-10NKP	Petite	50mm x 30mm	Non-kit 10/box
10-1800-10NKS	Small/Medium	65mm x 35mm	Non-kit 10/box
10-1800-10NKL	Large	69mm x 43mm	Non-kit 10/box

BEST SELLER

Sinus Packing And Splinting Cross Reference

Shippert Medical offers equivalents to competitor nasal packs - plus many more sizes! This chart offers easy cross referencing and ordering. All nasal and sinus packs are made from sterile, polyvinyl alcohol (PVA) foam.

SHIPPERT PART NUMBER	COMPETITORS PART NUMBER	DESCRIPTION
Nasal M-Type		
RH-7401-10	440400/400400	M-Type Small
RH-7408-10	440409/400409	M-Type Small with Airway
RH-7404-10	440402/400402	M-Type Large
RH-7405-10	440410/400410	M-Type Large with Airway
Nasal P-Type		
RH-2310-10	440406 (Pope)/400406	P-Type XL
RH-2315-20	400421 (Weimert)	P-Type Small
Nasal D-Type		
RH-7409-10	440424 (Doyle)/400424	D-Type Wideline Anatomical
RH-7411-10	440413 (Doyle)/400413	D-Type Slimline Anatomical
Sinus Nasal K & T Types		
RH-7530-10	400530 (Kennedy)	K-Type, Ultra-Slim
RH-7532-10	400428 (Toffel)	T-Type, Slim
<i>*All Shippert packs have a attached string</i>		
Ear Packs		
OT-5201-50	400141	Ear Speedi-Wick 15mm
OT-5202-50	400146	Ear Speedi-Wick 25mm
25-1500FNE-50AS	425008	Ear Fenestrated 15mm
OT-5204-10	400355	Ear Post-Op Small
OT-5203-10	400162	Ear Post-Op Large
Denver® Nasal Splints		
10-2000-05KL	1528136	Denver Splint, 2000 Kit Large
10-2000-05KS	1528126	Denver Splint, 2000 Kit Sm/Med
10-2000-05KP	1528116	Denver Splint, 2000 Kit Petite
10-2000-05KM	1528116	Denver Splint, 2000 Kit Mini
10-2000-10NKL	1528137/1528138	Denver Splint, 2000 Non-Kit Large
10-2000-10NKS	1528120/1528121	Denver Splint, 2000 Non-Kit Sm/Med
10-2000-05NKP	1528110/1528111	Denver Splint, 2000 Non-Kit Petite
10-2000-05NKM	1528110/1528111	Denver Splint, 2000 Non-Kit Mini
Thermoplast Nasal Splints		
10-4021-05KL	1529120	Denver Splint, 4000 Kit Large
10-4021-05KS	1529110	Denver Splint, 4000 Kit Sm/Med
10-4021-05KP	1529100	Denver Splint, 4000 Kit Petite
10-4021-10NKS	1529010	Denver Splint, 4000 Non-Kit Sm/Med
10-4021-10NKP	1529000	Denver Splint, 4000 Non-Kit Petite
10-4022-10NKL	1529020	Denver Splint, 4000 Non-Kit Large



SHOP ONLINE or CALL
www.innoviamedical.com/shippertmedical
 1.888.229.2875 | 651.789.3939



Shippert Medical and MD Resource have united under Innovia Medical to offer the most comprehensive fat grafting product options in the industry.

Tissu-Trans®
 Plastic and Cosmetic Surgery


Summit Medical, an Innovia Medical Company
 815 Vikings Parkway
 Suite 100
 St. Paul, MN 55121 | USA
 1.888.229.2875 | 651.789.3939
www.innoviamedical.com/shippertmedical

CAUTION: Federal (USA) law restricts these devices for sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in product labelling supplied with each device. Specifications and designs are subject to changes without notice or obligation on the part of the manufacturer.

All brands are trademark or registered trademarks of their respective owners.
 Denver Splint, Tissu-Trans and Rhino Rocket are registered trademarks of Shippert Medical. AbsorbENT and Breeze are all trademarks of Summit Medical.

©2019 Innovia Medical. All rights reserved. M-PLAS-5000 | RevA | SEPT2019

6248 South Troy Circle
Suite A
Centennial, CO 80111



1.888.229.2875 | 651.789.3939
www.innoviamedical.com/shippertmedical

SHIPPERT MEDICAL AND MD RESOURCE HAVE UNITED UNDER INNOVIA MEDICAL TO OFFER THE MOST COMPREHENSIVE FAT GRAFTING PRODUCT OPTIONS IN THE INDUSTRY.