



# MY DRUGTEST

4 EASY STEPS  
TO A  
PEACE OF MIND





# MY DRUGTEST

**1. DETACH CAP**

**AND DIP THE**

**STRIP INTO THE**

**CUP OF URINE**

**2. TAKE THE**

**STRIP OUT AFTER**

**THE 'C' LINE**

**APPEARS**

**3. REPLACE CAP**

**AND LAY ON**

**A FLAT SURFACE**


**4. READ RESULTS**

**AT 5 MINUTES**

**HC**




**NEGATIVE**




**NEGATIVE**



**POSITIVE**



**INVALID**



**INVALID**