



tilibit nanosystems GmbH  
Lichtenbergstraße 8  
DE 85748 Garching b. München  
Germany  
info@tilibit.com

## Prefab. nanostructure PF-2, cuboid with large aperture, staple mixture

---

**Conc.: 475 nM**  
**Vol.: 1053 µl**

---

**Store at -20°C**

**Product No.:** PF-2-WS-0

**Lot No.:** L-PF-2-WS-1

**Description:** This staple mixture contains all 208 oligonucleotide sequences necessary to form the prefab. nanostructure PF-2, cuboid with small aperture design. The oligonucleotides are dissolved in ultrapure water at an individual sequence concentration of 475 nM.

**Detailed usage recipes:**

Castro CE, et al: "A primer to scaffolded DNA origami" — Nature Methods. 2011 Mar; 8(3): 221-9

