



tilibit nanosystems GmbH
Lichtenbergstraße 8
DE 85748 Garching b. München
Germany
info@tilibit.com

Prefab. nanostructure PF4 Gear, staple mixture

Conc.: 425nM
Vol.: 1177 µl

Store at -20°C

Product No.: PF4-2-0

Lot No.: L-PF4-2-1

Description: This staple mixture contains all 234 oligonucleotide sequences necessary to form the prefab. nanostructure PF4 gear. The oligonucleotides are dissolved in ultrapure water at an individual sequence concentration of 425 nM.

Detailed usage recipes:

Castro CE, et al: "A primer to scaffolded DNA origami" — Nature Methods. 2011 Mar; 8(3): 221-9

