



tilibit nanosystems GmbH  
Lichtenbergstraße 8  
DE 85748 Garching b. München  
Germany  
info@tilibit.com

## Prefab. nanostructure PF5, high asymmetry design, staple mixture

---

**Conc.: 500nM**  
**Vol.: 1000 µl**

---

**Store at -20°C**

**Product No.:** PF5-2-0

**Lot No.:** L-PF5-2-1

**Description:** This staple mixture contains all 161 oligonucleotide sequences necessary to form the prefab. nanostructure PF5, high asymmetry design. The oligonucleotides are dissolved in ultrapure water at an individual sequence concentration of 500 nM.

**Detailed usage recipes:**

Castro CE, et al: "A primer to scaffolded DNA origami" — Nature Methods. 2011 Mar; 8(3): 221-9

