

PRODUCT BOOKLET génia





NOTES	
	W. States and Control

CONTENTS

E-collars	04
Protective Equipment	12
Mano Médical	20
Dressings	24
Bandages	37
Urinary Catheter Range	42
Disposable Syringes & Needles Range	46
Vétopop Apparel	52

Génia USA Inc. 1411 Broadway, 16th Floor New York, NY 10018 – USA

contact@geniausa.com






NEOCOL®

Transparent E-Collar with neoprene border

- Transparent polypropylene for better visibility
- Quick-fastening system with secured double-tab closure
- Separate fastenings for collar
- Outer anti-scratch engraving
- Neoprene border acts as "shock absorber"
- Easy storage system and color-coded sizes





Code	Size
306420 •	7.5 cm
306422 •	10 cm
306424	12.5 cm
306426 •	15 cm
306428 •	20 cm
306430 •	25 cm
306432 •	30 cm
306438	Set 2 of each size from 7.5 to 30 cm























FUNCOL®

Transparent E-collar with paw-print design

- Quick-fastening system with secured double-tab closure
- Outer anti-scratch engraving
- Scalloped border acts as "shock absorber"
- Reduced echo and static
- Easy storage system and color-coded sizes

The range is expanding! Sizes 35 cm and 40 cm available





Code	Size
306200 •	7.5 cm
306202 •	10 cm
306204 •	12.5 cm
306206 •	15 cm
306218* •	17.5 cm
306208 •	20 cm
306210 •	25 cm
306212 •	30 cm
306213* •	35 cm
306214* ●	40 cm
306216	Set 2 of each size

*Not included in the set.

from 7.5 to 30 cm



www.geniausa.com

7



GREENCOL

Transparent ecofriendly E-collar

- Transparent polypropylene for better visibility
- Quick-fastening system with secured double-tab closure
- Outer anti-scratch engraving
- Scalloped border acts as "shock absorber"
- Reduced echo and static
- Easy storage system







Code	Size
306300	7.5 cm
306302	10 cm
306304	12.5 cm
306306	15 cm
306308	20 cm
306310	25 cm
306312	30 cm
306316	Set 2 of each size from 7.5 to 30 cm





VELFAST®

Transparent E-collar with quick Velcro attachment system

- Transparent and silent peripheral neoprene edging absorbs impacts and noise
- Atraumatic neoprene edging around the E-collar for more comfort and quicker adjustments
- Simple and quick Velcro attachment with excellent tear resistance
- Velcro strip adjustment for different neck circumference
- Outer anti-scratch engraving and reinforced with neoprene lining
- Color-coded according to size





Code	Size
306450 •	7.5 cm
306452 •	10 cm
306454 •	12.5 cm
306456 •	15 cm
306458 •	20 cm
306460 •	25 cm
306462 •	30 cm
306468	Set 2 of each size from 7.5 to 30 cm





COSYPAD® NEW GENERATION


Versatile medical grade anti-microbial mat

- Keeps patients clean, comfortable, and dry
- Can be used before, during and after surgery
- Non-absorbent liquids pass through the anti-microbial mesh
- Reduces hypothermia unlike towels and bare stainless steel
- Reduces pressure on muscles/joints to prevent bedsores
- Woven with one, continuous thread of polypropylene
- Rinse down or machine wash, drip dry, and repeat













Continuous sewing	Code	Dimensions	Weight
0 · · · · · · · · · · · · · · · ·	327200B	13 x 20 in	200 g
	327201B	15 x 28 in	382 g
	327203B	22 x 28 in	505 g

327205B

24 x 45 in



www.geniausa.com

884 q



VACUUMAT®

Vacuum immobilization mattress

- Made of sturdy polyurethane material
- 15% more microbeads to perfectly conform to the animal's shape
- Vacuum system allows complete immobilization in the desired position ideal for radiology, ultrasound, MRI, CT scan and surgery
- Thermal insulation, flexible, comfortable, and sturdy material
- High quality safety valve, which can maintain air for more than 72 hours
- Radiolucent and easy to clean









DRYSOFT®

Highly absorbent pads made of super absorbent gel

- Five layers of fast absorbing padding = no leaks!
- Absorbs up to twice as much as similar product / you will use less
- Anti-microbial polymer gel increases absorption and reduces odor








DRAPVET

Génia disposable surgical drapes, non-woven, sterile

- One smooth waterproof layer
- One non-woven, pre-sterilized, and absorbent layer that is soft like fabric
- 53g/m²





Code	Dimensions
317178	18 x 30 in
317180	30 x 35 in
317182	49 x 49 in
317184	79 x 118 in





CARBON TECH

The Carbon Tech V2 Heating Mat is designed to be particularly robust while also providing an optimal heat distribution evenly over its entire surface.

- Maintains patient temperature during anesthesia, so recovery is simple and safe
- Gently warms the patient while minimizing the risk of burns
- High density carbon fabric/flexible structure
- Liquid-resistant coating that is soft and easy to clean
- Integrated probe with temperature sensor
- Integrated temperature limiter (107.6°F 42°C)
- 3 power levels: 50%, 80% and 100%
- 5 years guaranteed from scratches and bites!







LASER FLEX LT16

Cold laser therapy, typical applications in veterinary medicine:

- Pain treatment
- Facilitation post-surgical healing, aiding wound healing
- Inflammatory symptom reduction
- Decrease arthritis pain
- Edema reduction
- Pre-established protocols for the different pathologies of cats, dogs, horses and NACs (reptiles, small mammals, and birds)
- Protocols adapted according to the color of the animal's coat
- Ability to modify or create new, personalized protocols
- Integration of unique protocols for different phases of continuous, low, medium, or high frequency treatment






(J Z -S S Ш \square

ANIPLAST® ALGI

Super absorbent calcium and sodium alginate dressing with silver ions

- Silver ions with antibacterial properties
- Forms a protective gel on wound surfaces
- Hemostatic





Code	Dimensions
327020	4 x 4 in
327022	4 x 8 in




C S ഗ Ш

ANIPLAST® SURGI

Non-adhesive absorbent dressing

- Non-woven absorbent pad that prevents maceration and infection
- Micro-perforated polyethylene sheets on each surface which limits maceration





Code	Dimensions
327030	4 x 4 in
327032	4 x 8 in



ANIPLAST® HYDRO

Absorbent healing dressing

- Gelatin free hydrocolloid layer applied to a polyurethane film to form a gel on the wound
- Maintains a moist and warm environment on the surface of the wound which promotes granulation and epithelialization along the wound's self-debridement
- Bacteria resistant





Code	Dimensions
327014	4 x 4 in
327016	8 x 8 in

29




S C (N) S Ш Ω

ANIPLAST® TULLO

100% cotton absorbent gauze infused with white paraffin

- Made from reticular cotton fabric that allows exudate to drain from the wound ensuring no maceration
- Pro-inflammatory activity that promotes granulation






Code	Dimensions
327006	4 x 4 in
327008	4 x 8 in



ANIPLAST® SPONGI

Super absorbent dressing

- Non-adhesive, semi-occlusive foam dressing on a thick water-resistant polyurethane film
- Permeable to gas and water vapors but resistant to bacteria
- Limits maceration and infection
- Keeps the wound surface moist and warm which encourages wound's self-debridement, granulation, and epithelialization





Code Dimensions	
327000	4 x 4 in
327002	8 x 8 in



THE CHOICE IS BASED ON THE VETERINARIAN'S CLINICAL ASSESSMENT

TYPE OF WOUND		DRESSINGS	DRESSING CHANGE
	Non-infected wound. Surgical wound. Healing by first intention. Low exudate wound.	SURGI	When saturated. From every day to every few days.
	Wound in debridement stage. High or moderate exudate and/or cavitary wound.	ALGI Routh side against the wound	When saturated. From every day to every few days.
	Wound in debridement stage or in granulation stage. Moderate exudate wound and/or without cavitation.	HYDRO	Every three days on average. Change if discoloured or swollen.
	Wound in granulation stage. Low exudate wound and/or burns.	TULLO	From every day to every few days.
	Wound from debridement to epithelialisation. Medium/high exudate wound Aural hematoma.	SPONGI Foam side against the wound	When saturated or every five days.

ANIPLAST® SET

Set of ANIPLAST® dressings for small and large animals



Code	Content	
327040	Small animals set 4 of each dressing 4 x 4 in	
327042	Large animals set 2 of each dressing 4 x 8 in (SURGI, ALGI, TULLO) 8 x 8 in (HYDRO, SPONGI)	





ADERPLAST®

White adhesive elastic bandage

- Aerated strip and 100% woven cotton
- Flexible elastic band for effective tightening
- Body made from highly adhesive surgical rubber
- Protection from siliconized paper, low noise with edging for easy removal
- Good skin tolerance
- Excellent resistance to repeated movement



. Aerated strip 100 % woven cotton

		- 1
Alerman		Abelle
	1	
Sinia	Sil -	
	- general	

Code	Dimensions
307250	1.2 x 98 in
307252	2.4 x 98 in
307254	3.1 x 98 in
307256	4 x 98 in





SOFT SPLINT ROLL

Aluminum roll splints with foam on both sides

- Flexibility to fit any shape
- Soft and resilient for efficient immobilization
- Lightweight, compact, and comfortable
- Washable closed-cell foam
- Easy to cut, reducing sharp edges





Code	Dimensions
307830 •	1.8 x 10 in
307832 •	2.2 x 14 in
307834 •	3.5 x 16 in





FUN STRAP®

Durable medical-grade bandages

- Self-adhesive with no sticking to hands or fur
- Fun colors and designs for pet owners
- Odor protection to prevent biting and licking
- Easy to adjust and tear, no scissors necessary










Starting strip



Code	Dimensions
307240	2 x 177 in
307241	2.8 x 177 in
307242	3.9 x 177 in





UROLENE® CAT

Premium quality sterile feline urinary catheter

- 360-degree rotating collar with 4 stitching holes
- Rounded end with 2 side holes
- Atraumatic for easy insertion into the urethra
- Luer-Lock hub making it easy to flush or collect with a Luer-Lock syringe
- Metal stylet provided for easier insertion
- DEHP-free & radiotransparent
- Sterile and individually wrapped











Code	Model	Dimensions
306900 •	With stylet	3FG - 1.0 x 110 mm
306902 •	With stylet	4FG - 1.3 x 140 mm
306904 •	Without stylet	3FG - 1.0 x 110 mm
306908 •	Without stylet	4FG - 1.3 x 140 mm



\square Ш \vdash ш 工 U R I N A R C A C A RANG



UROLENE® DOG

Premium quality sterile canine urinary catheter

- Rounded end with 2 side holes
- Atraumatic for easy insertion into the urethra
- Luer-Lock hub making it easy to flush or collect with a Luer-Lock syringe
- DEHP-free & radiotransparent
- Sterile and individually wrapped







Code

Dimensions

	306920 •	4FG - 1.3 x 300 mm
	306922 •	5FG - 1.7 x 350 mm
187. Annual Control of the Control o	306924 •	6FG - 2.0 x 450 mm
Treas secon than they write a cateful of Remissibilities for hand / Sonda vesciciate cateful (Sonda of Sonda of	306926 •	8FG - 2.6 x 500 mm
urnaria perro	306928 •	10FG - 3.3 x 600 mm





DISPOJECT® SAFETY INSULIN SYRINGES

Three parts 40 IU insulin syringe of 1 ml with a 29G extra fine needle

- Sterile, disposable, and latex free
- Innovative protective system to lock the needle after use for safety and ease of mind





Code Capacity Graduation
107136 1 mL 1 IU



MEIIII Landsid. Lind III

DISPOJECT® DISPOSABLE SYRINGES

Disposable three-piece polypropylene syringes

- Latex free, silicone-free, DEHP-free, styrene free and BPA free
- Smooth glide and coned tip to prevent loss of product
- Tipped plunger to ensure 1 ml size high accuracy
- Sterilized with ethylene oxide

1 mL Tuberculin model







Code	Capacity	Hub	Graduation
107140* •	1 ml	Concentric Luer	1/10 mL
107142	2 ml	Concentric Luer	1/10 mL
107144 •	5 ml	Concentric Luer	1/5 mL
107146 •	10 ml	Eccentric Luer	1/5 mL
107148 •	20 ml	Eccentric Luer	1 mL
107150 •	30 ml	Eccentric Luer	1 mL
107152	50 ml	Eccentric Luer	1 mL
107154 •	50 ml	Feeding	1mL

*Tuberculin model.





DISPOSABLE NEEDLES

Stainless steel, sterile, disposable, hypodermic needles

- Long and tri-beveled made for canine and feline skin
- Compatible with Luer or Luer-Lock hubs
- Resin hub for extra strength to prevent bending
- Easy to open paper packaging
- Sterilized with ethylene oxide







Code	Dimensions	Dimensions
147752 •	16 x 0.5 mm	25G x 5/8
147754 •	16 x 0.6 mm	23G x 5/8
147762 •	16 x 0.8 mm	21G x 5/8
147764 •	25 x 0.8 mm	21G x 1
147766 •	38 x 0.8 mm	21G x 1 1/2"
147770 •	25 x 0.9 mm	20G x 1
147772 •	38 x 0.9 mm	20G x 1 1/2
147776	38 x 1.1 mm	19G x 1 1/2
147780 •	38 x 1.2 mm	18G x 1 1/2"






COOLING JACKET

Cooling jacket for dogs

- Improves your pet's well-being and reduces stress
- With its modern micro-perforated cooling fibers, this jacket keeps your pet cool for several hours
- Fastens with Velcro to put on easily
- Reusable
- Hand wash
- Available in 4 sizes









Moisten with cool water and wring out before use

Code	Dimensions	Length
700700	S	14 in
700701	M	16 in
700702	L	20 in
700703	XL	24 in



NOTES



NOTES



