




# PICCPerfect<sup>®</sup>

## Antimicrobial PICC Line Cover


### FEATURES



Two access points for lower or higher insertion site



Bonded anti-slip edges



Antimicrobial, moisture-wicking and ultra-soft fabric

- Patented double-layer design holds tubing and lumen(s)
- Bonded anti-slip edges help maintain cover in place
- Multiple colors & sizes available for adults and children
- Non-adhesive, latex-free

### BENEFITS

- Supports PICC and midline care at home or on-the-go
- Holds tubing and lumen(s) between double-layer design
- Helps prevent possible line entanglement during activities of daily living
- Helps prevent discomfort caused by the contact of lines with the skin
- Increase satisfaction and comfort with a PICC (Based on 1,000+ reviews)



# PICCPe<sup>®</sup>fect

## Antimicrobial PICC Line Cover

### AVAILABLE COLORS AND SIZES

Below find the available sizes and colors of the PICCPe<sup>®</sup>fect.

ITEM NAME	SIZE	SKU
PICCPe <sup>®</sup> fect, Midnight Black	XS	PLC4-1MID-XS
PICCPe <sup>®</sup> fect, Midnight Black	S	PLC4-1MID-S
PICCPe <sup>®</sup> fect, Midnight Black	M	PLC4-1MID-M
PICCPe <sup>®</sup> fect, Midnight Black	L	PLC4-1MID-L
PICCPe <sup>®</sup> fect, Midnight Black	XL	PLC4-1MID-XL
PICCPe <sup>®</sup> fect, Slate	XS	PLC4-1SLT-XS
PICCPe <sup>®</sup> fect, Slate	S	PLC4-1SLT-S
PICCPe <sup>®</sup> fect, Slate	M	PLC4-1SLT-M
PICCPe <sup>®</sup> fect, Slate	L	PLC4-1SLT-L
PICCPe <sup>®</sup> fect, Slate	XL	PLC4-1SLT-XL

### SIZE GUIDE

Measure the circumference of the arm at the PICC insertion site using a flexible tape measure, and select the recommended size in the chart.

*Ex: if circumference is 8", then use size XS. If in between two sizes, use the smaller size.*

RANGE	SIZE
7.5" - 9"	XS
9 - 10.5"	S
10.5" - 12.5"	M
12.5" - 14.5"	L
14.5" - 18"	XL

### COMPOSITION

85% Polyester  
10% Spandex  
5% Recycled S Cafe

### CARE

The suggested use of the PICCPe<sup>®</sup>fect is 1 week before washing. Please refer to the care instructions on the product tag on the inner layer.

### ★★★★★ CLINICIAN REVIEW

"Patients are less likely to pick at dressing edges and it keep the dressing more secure."  
- Infusion Clinic, Clearwater, FL 01/27/2023



### FOR MORE INFORMATION

Scan the QR code for a demo video of the PICCPe<sup>®</sup>fect.