



BRAVO Classic

BRAVO Twsit

BRAVO Series™

Powered by APOLLO&ATLAS

Tremor Assistance & Quantitative Evaluation System



Tremor is severely affecting the patients and caregivers in many aspects:



Loss of independence
on Eating



Avoid Social Contact



Caregiving Stress



Psychological Burden

Innovative Therapy and Symptom Evaluation System in The Intelligent Digital Health Era

Solving functional problem and realize symptom evaluation by self-adaptive micro robot


BRAVO Series™

Powered by APOLLO&ATLAS

Tremor Assistance & Quantitative Evaluation System



Eating aids for patients with hand tremors



Tremor Reduciton



Eat independently with a BRAVO spoon



Unable to lift food up with a normal spoon

A patient with Parkinson's Disease in Finland



BRAVO Twist™

Powered by APOLLO&ATLAS

Tremor Assistance & Quantitative Evaluation System

Innovative Therapy and Symptom Evaluation System in The Intelligent Digital Health Era

Solving functional problem and realize symptom evaluation by self-adaptive micro robot

BRAVO Twist™

Powered by APOLLO&ATLAS



For patients with Parkinson's Disease or Essential Tremor



Micro-robotic technology realizes automatic eating assistance function



Real-time tremor monitoring and quantitative evaluation



Check symptom data and reports through Mobile APP



User-friendly Accessibility



Health Care Funding Support in Europe



Innovative Therapy and Symptom Evaluation System in The Intelligent Digital Health Era

Solving functional problem and realize symptom evaluation by self-adaptive micro robot



BRAVO Twist™

Powered by APOLLO&ATLAS


Tremor Assistance & Quantitative Evaluation System



Eating aids for patients with hand tremors



Tremor Reduciton



Checking hand tremor data through Mobile APP(GYENNO SPOON)



Disease Condition Assessment



Disease Progression Tracking




Treatment Effectiveness Assessment

Innovative Therapy and Symptom Evaluation System in The Intelligent Digital Health Era

More user-friendly and auxiliary functions to improve the patient experience

The Twist

Enjoy noodles with the help of the rotating fork attachment




Left/right hand Auto-recognition

Always keep the attachment facing upwards by auto-recognition of holding hands




Innovative Therapy and Symptom Evaluation System in The Intelligent Digital Health Era

More user-friendly and auxiliary functions to improve the patient experience




Magnetic Attachments

Washable & Easy to install and detach




Built-in Rechargeable Battery

One week of use with a single charge



Protective Case

Cover and go



Protective Strap

Keep safe from accidental dropping

GYENNO SCIENCE



BRAVO Classic™

Powered by APOLLO&ATLAS

Tremor Assistance

Innovative Therapy and Symptom Evaluation System in The Intelligent Digital Health Era

Solving functional problem and realize symptom evaluation by self-adaptive micro robot

BRAVO Classic™

Powered by APOLLO&ATLAS



For patients with Parkinson's Disease or Essential Tremor



Micro-robotic technology realizes automatic eating assistance function



User-friendly Accessibility




Health Care Funding Support in Europe




Innovative Therapy and Symptom Evaluation System in The Intelligent Digital Health Era

More user-friendly and auxiliary functions to improve the patient experience




Magnetic Attachments

Washable & Easy to install and detach




Built-in Rechargeable Battery

One week of use with a single charge



Protective Case

Cover and go



Protective Strap


Keep safe from accidental dropping

Comparison



BRAVO Twsit





Eating Aids




Tremor Data Collection





Twist




Magnetic Attachments
(Spoon & Fork)



Protective Cover



Protective Strap



BRAVO Classic

