

HBV LONGITUDINAL PANEL

(Catalog No. 6517)

For Research Use Only

Application Notes: This 13 member panel was collected from a 18 year old Caucasian Male.

PANEL NUMBER	DATE OF DRAW	ABBOTT HBSAG S/CO	ABBOTT HBC TOTAL S/CO	ABBOTT HBC IGM S/CO	DIASORIN HBEAG S/CO	DIASORIN ANTI-HBE S/CO	ABBOTT ANTI-HBS S/CO	ORTHO ANTI-HCV ECI	ORTHO ANTI-HIV ½ ECI
6517-01	6/21/01	>40.8	0.06	>6.4	0.4	0.5	0.3	25.1	0.08
6517-02	7/4/01	0.02	0.05	>6.4	0.1	0.2	0.8	24.9	0.06
6517-03	7/19/01	0.31	0.02	>6.0	0.3	0.06	1.25	24.5	0.07
6517-04	8/1/01	0.05	0.02	>6.0	0.3	0.05	1.77	23.5	0.07
6517-05	8/15/01	0.05	0.07	5.0	0.4	0.01	3.1	23.4	0.14
6517-06	8/30/01	0.04	0.05	3.2	0.2	0.02	5.0	22.0	0.07
6517-07	9/13/01	0.02	0.04	2.0	0.3	0.01	9.9	23.9	0.08
6517-08	9/27/01	0.02	0.06	1.6	0.2	0.02	7.5	19.9	0.08
6517-09	10/11/01	0.17	0.06	1.8	0.3	0.02	16.0	20.8	0.08
6517-10	10/25/01	0.04	0.06	1.4	0.2	0.02	15.7	28.8	0.06
6517-11	11/8/01	0.08	0.06	1.3	0.3	0.02	15.7	29.0	0.07
6517-12	11/21/01	0.08	0.05	1.3	0.3	0.03	20.9	29.8	0.06
6517-13	12/4/01	0.05	0.04	1.2	0.3	0.02	>21.0	30.0	0.13

*Donor is co-infected with Hepatitis C.