

THANK YOU for your interest in the Dignity™ AllDay 400 Wheelchair. A specialty chair designed around 3 pillars of technology: **Independence, comfort, and mobility.**

The result is a safer chair that improves user quality of life, and significantly reduces the need for toilet transfers with its patented drop-seat technology.



"I DON'T REALLY KNOW HOW MY WIFE AND I WOULD GET ALONG WITHOUT THE DIGNITY™ ALLDAY CHAIR BECAUSE I AM NO LONGER ABLE TO EASILY TRANSFER HER."

-DAVID J. HUSBAND OF USER

The chair is not only safer for users, but reduces the need for caregivers during toilet procedures; leaving users with a greater sense of dignity.

Clinical trials prove that the use of a Dignity™ AllDay 400 Wheelchair reduces dangerous lifting by more than 50% making it ideal for both age-in-place and institutional settings such as hospitals, assisted living facilities and skilled nursing facilities.



5 Year Frame Warranty



Our standard Dignity™ AllDay 400 Wheelchair frames and wheels have passed the most stringent FDA required strength testing. Users of The Dignity™ AllDay 400 Wheelchair who weigh more than three hundred pounds will indicate upon the placement of an order their need for a specially reinforced bariatric toileting chair. The Dignity™ AllDay 400 is not intended for shower use.

Dignity and TruGlides are trademarks of Majestic Medical LLC. © Majestic Medical LLC 2020. All rights reserved. Printed in the U.S.A. Please visit dignityallday.com for warranty and liability information.

TECHNICAL SPECIFICATIONS

Frame Type & Material	Solid fixed J-frame configuration* with 875" O.D. 1020 CDS DOM Steel
Seat Width	Standard – 16", 18" and 20" (Bariatric – 22", 24" and 25")
Seat Depth	18"
Seat Height	Adjustable – 16" - 20" in one-inch increments
Backrest Height	16"
Overall Chair Width	Add 9.5" to seat width, e.g. 20" seat has 29.5" overall chair width
Overall Chair Height (to top of handles bars)	Adjustable – 34" - 38" in one-inch increments
Chair Weight Limits	Standard max weight 300 lbs. Reinforced (bariatric) max weight 600 lbs.
Overall Chair Weight	65 lbs.
Front Wheel Assemblies	TruGlides™ shock absorbing/self-righting/anti-flutter safety wheels, height adjustable in one-inch increments
Front Wheel Casters	Standard – 6" x 1" (Bariatric – 6" x 2")
Rear Wheels	24" x 1" (Six-spoke design w/ urethane tire)
Rear Axle	Height adjustable in one-inch increments
Armrest	Solid, fixed w/padded rest
Leg Rests	Standard: two-pin, adjustable length only
Commode Bucket	A detachable and portable commode bucket for users with urgency issues
Frame Color	 Majestic Blue

CONTACT US TODAY TO LEARN MORE

Phone | 1-800-940-4326

Email | info@MajesticMedicalLLC.com

Web | dignityallday.com

Majestic Medical, LLC
Sioux Falls, South Dakota 57105

Independence for the wheelchair user and caregiver alike!



MAINTAIN YOUR DIGNITY



INTRODUCING
THE DIGNITY™ ALLDAY
400 WHEELCHAIR


INDEPENDENCE



A: Drop seat in the up position for all-day use.



B: Drop seat lowered for transfer-free toileting.

 **See how it works:**
DignityAllDay.com

The Dignity™ AllDay 400 Wheelchair was designed as a manual fixed-frame wheelchair for all-day use with an amazing feature — it has a lever-controlled drop-seat system that allows the user to back over and use a toilet without ever needing to do a transfer.

Our chair easily fits over standard toilets. The result is a safer toileting process that provides greater user independence.

Additionally, it can greatly reduce, and in some cases eliminate, the need for supervision and handling by caregivers. If the average person uses the restroom 6-7 times per day, using the Dignity™ AllDay 400 would **eliminate over 2000 transfers annually!**

MOBILITY

The Dignity™ AllDay 400 Wheelchair is the complete package with the innovation of the TruGlides™ Caster System featuring anti-shock and anti-flutter technology. These patented front suspension casters are vital in enhancing user mobility, comfort, and safety.



Vibrations from rough surfaces can contribute to user fatigue, and be the cause of discomfort and even injury. The TruGlides™ Caster front suspension greatly reduces wheelchair vibration leading a smooth and comfortable ride on all surfaces. TruGlides™ helps casters align automatically, reducing wheel flutter and drag, for better maneuverability, ultimately giving the user a safer and more comfortable ride.



Beyond excellent mobility, the Dignity™ AllDay 400 Wheelchair is built to last. Testing exceeded the ISO 7176-8 requirements for static, impact, and fatigue strength. The results are directly related to structural stability in chair design and wheel performance.

COMFORT



Our memory foam cushions are designed for maximum comfort during long periods of use. These 2-1/2" thick cushions are covered with medical quality vinyl, these cushions are durable and easy to clean with mild soap and water or a disinfectant solution.

Note: It is recommended that wheelchair users work with a seating specialist prior to using a new wheelchair. Also, memory foam cushioning can also help minimize the development of pressure sores.

