



# ARYSE

## MALLET FINGER SPLINT

### Features


- Bilateral
- Lightweight
- Simple to apply

### Directions for Use

1. Place affected finger in splint
2. With medical tape (not included), tape splint around finger until finger is secure and comfortable

### Scope of Application

- Finger sprains
- Mallet finger
- Stable finger fractures
- Volar plate injuries
- Among others



<u>Finger Splint</u>	<u>Size</u>	<u>Circumference of Injured Joint</u>	<u>Length of Tip of Finger to Injured Joint</u>
100	X Small	up to 1.75"	up to 2"
101	Small	1.75 - 2"	up to 2"
102	Medium	2 - 2.25"	up to 2.25"
103	Large	2.25 - 2.375"	up to 2.5"
104	X Large	2.375 - 2.5"	up to 2.25"
105	2X Large	2.5 - 2.75"	up to 2.375"
106	3X Large	2.75 - 3"	up to 2.5"
107	4X Large	3" and up	up to 2.75"

\*Measure circumference of injured joint first and then length of tip of finger to injured joint

1801 N. 1st Street Suite 2 Lincoln, NE 68528

www.aryse.com



# ARYSE

## MALLET FINGER SPLINT

### Features


- Bilateral
- Lightweight
- Simple to apply

### Directions for Use

1. Place affected finger in splint
2. With medical tape (not included), tape splint around finger until finger is secure and comfortable

### Scope of Application

- Finger sprains
- Mallet finger
- Stable finger fractures
- Volar plate injuries
- Among others



<u>Finger Splint</u>	<u>Size</u>	<u>Circumference of Injured Joint</u>	<u>Length of Tip of Finger to Injured Joint</u>
100	X Small	up to 1.75"	up to 2"
101	Small	1.75 - 2"	up to 2"
102	Medium	2 - 2.25"	up to 2.25"
103	Large	2.25 - 2.375"	up to 2.5"
104	X Large	2.375 - 2.5"	up to 2.25"
105	2X Large	2.5 - 2.75"	up to 2.375"
106	3X Large	2.75 - 3"	up to 2.5"
107	4X Large	3" and up	up to 2.75"

\*Measure circumference of injured joint first and then length of tip of finger to injured joint

1801 N. 1st Street Suite 2 Lincoln, NE 68528

www.aryse.com



**ARYSE**

---

**MALLET FINGER SPLINT**

**Warranty**

All products, with the exception of exo-performance products, sold by ARYSE are warranted against defects in material and workmanship for a period of one (1) year from the date of shipment. All exo-performance products sold by ARYSE are warranted against defects in material and workmanship for a period of six (6) months. If you believe any ARYSE product you have purchased has a defect in material or workmanship or has failed during normal use within the warranty period, please contact ARYSE at 855-877-2886 for assistance. Please call your supplier with any questions unrelated to the warranty. If product repair or replacement is necessary, the customer will be solely responsible for all shipping charges, insurance, and proper packaging to prevent breakage in transit. This warranty does not apply to defects resulting from any customer actions, such as mishandling, operation outside of design limits, misapplication, improper repair, or unauthorized modification. No other warranties are expressed or implied.

**Cleaning and General Care**

Hand wash with lukewarm water and mild soap. Rinse well making sure to remove all soap residue as it can cause skin irritation and deterioration of the support. Air dry. Do not expose to direct heat or sunlight while drying to prevent shrinkage and damage to your support.

**Notice**

While every effort has been made to obtain maximum strength, function, and durability, there is no guarantee that further injury will be prevented through the use of this device. Intended for single patient use only. Intended to be dispensed by the prescription or order of a licensed practitioner. If you experience pain while using this product, discontinue use and consult your doctor.

© 2021 Cruze Distribution LLC. All rights reserved.

---

1801 N. 1st Street Suite 2 Lincoln, NE 68528

[www.aryse.com](http://www.aryse.com)



**ARYSE**

---

**MALLET FINGER SPLINT**

**Warranty**

All products, with the exception of exo-performance products, sold by ARYSE are warranted against defects in material and workmanship for a period of one (1) year from the date of shipment. All exo-performance products sold by ARYSE are warranted against defects in material and workmanship for a period of six (6) months. If you believe any ARYSE product you have purchased has a defect in material or workmanship or has failed during normal use within the warranty period, please contact ARYSE at 855-877-2886 for assistance. Please call your supplier with any questions unrelated to the warranty. If product repair or replacement is necessary, the customer will be solely responsible for all shipping charges, insurance, and proper packaging to prevent breakage in transit. This warranty does not apply to defects resulting from any customer actions, such as mishandling, operation outside of design limits, misapplication, improper repair, or unauthorized modification. No other warranties are expressed or implied.

**Cleaning and General Care**

Hand wash with lukewarm water and mild soap. Rinse well making sure to remove all soap residue as it can cause skin irritation and deterioration of the support. Air dry. Do not expose to direct heat or sunlight while drying to prevent shrinkage and damage to your support.

**Notice**

While every effort has been made to obtain maximum strength, function, and durability, there is no guarantee that further injury will be prevented through the use of this device. Intended for single patient use only. Intended to be dispensed by the prescription or order of a licensed practitioner. If you experience pain while using this product, discontinue use and consult your doctor.

© 2021 Cruze Distribution LLC. All rights reserved.

---

1801 N. 1st Street Suite 2 Lincoln, NE 68528

[www.aryse.com](http://www.aryse.com)