

### **Introducing:** CLINITEST Rapid COVID-19 Antigen Self-Test



CLINITEST®
Rapid COVID-19 Antigen Self-Test
Corean tela-ske:
3 intriviors Tables with halfer and Tips
Intractions for tips

SIEMENS
Healthineers

C T

Simulation Tables with halfer and Tips
Intractions for tips

SIEMENS

Realthineers

SIEMENS

T

#### **Self-Test kits includes:**

- Test devices (cassettes)
- Sterile swabs
- Instructions for use
- Prefilled extraction tubes with buffer solution and tips

For the moments that matter.

300 million+ CLINITEST COVID-19 tests already acquired outside U.S. and counting...

### **Product Details**<sup>1</sup>



- OTC (non-prescription); self-guided; use on asymptomatic individuals
- Qualitative; visually interpreted, no instrument required
- Omicron SARS-CoV-2 variant unlikely to impact performance of CLINITEST Rapid COVID-19
   Antigen (Self-)Test<sup>2</sup>
- Nasal swab samples
- Results in 15 minutes
- 5 tests / box option available for U.S. market initially; 1 test / box option added over time
- Room temperature or refrigerated storage
- 11-month shelf life (from date of original manufacturing)
- 86.5% Sensitivity; 99.3% Specificity
- Dedicated Customer Support Line
- U.S. webpage containing product details and How-To tutorial video; website accessible via QR code printed on outside of test box
- This test is also authorized for non-prescription home use with self-collected anterior nasal (nares) swab samples from individuals aged 14 years or older, or adult-collected anterior nasal (nares) swab samples from individuals aged 2 years or older, with or without symptoms or other epidemiological reasons to suspect COVID-19 when tested twice over three days with at least 24 hours (and no more than 48 hours) between tests.

### **CLINITEST Rapid COVID-19 Antigen Self-Test<sup>1</sup>**



### **Click** to See How It Works:



<sup>1.</sup> Based on FDA Emergency Use Authorization

<sup>2.</sup> Based on initial assessment conducted by Siemens Healthineers

## 3 Reasons Why CLINITEST Rapid COVID-19 Antigen Self-Test <sup>1</sup>



1

### **Brand**

trust, from a brand already relied on globally to support the world's toughest health challenges



2

### Quality

Accurate results, test after test:

11-Month Shelf Life<sup>1</sup>

86.5% Sensitivity 99.3% Specificity

Omicron SARS-CoV-2 variant unlikely to impact performance<sup>1</sup>

**300 million+** CLINITEST COVID-19 tests already acquired outside U.S.



3

### Supply

# Major, Sustained Supply

Tens of millions of tests made available to the U.S. market within weeks of EUA

Major sustained supply available after that



