LM-101



The Marco Lensmeter 101, like all Marco lensmeters, is well suited for measuring soft and hard contact lenses as well as conventional spectacle lenses. The Marco Lensmeter 101 is a compact deluxe model with a full 90-degree inclination, which facilitates the measurement of contact lenses when used in the vertical position. Standard equipment on the 101 includes such deluxe features as a prism compensator, external power and axis readings, LED illumination, AC or battery-powered operation, and an American-style cross-line target. The oversized lens table permits the measurement of conventional lenses from 30 to 90 millimeters in diameter. The instrument is well suited for laboratory or office use.



Play it safe.

Protect your warranty. Buy your Marco product only from an authorized Marco distributor. Call Marco for the name of the distributor nearest you.

LENSMETER 101

SPECIFICATIONS Stock #4004

Compact Deluxe Model with Full 90 Degree Inclination

• FEATURES

LED Illumination system

AC Adapter or Battery powered operation

Full 90-degree inclination

Prism compensator

External power and axis readings American-style cross-line target Acceptable for laboratory use

Incomparable construction
Discriminating optics

Will read hard or soft contact lenses and conventional lenses from 24mm through 90mm in diameter

TECHNNICAL DAT

Type Ocular type; external reading

Target American Cross-line

Vertex Power Range: +25 to -25 diopters

tep: 0.125 diopters up to ±3 diopters

0.25 diopters beyond ±3 diopters

Cylindrical Axis Range: 0 to 180 degrees

Step: 1 degree

Prismatic Power Range: 5 prism diopters in field of view

+5 to -5 diopters

(20 prism diopters by prism compensator)

Step: 1 prism diopter

Eyepiece Focusing

Diameter of Lenses

Accommodated 24 to 90mm

Lamp AC 100V (LED illumination)

Target Illumination Switchable Dual color, Green/Orange LED (2V 25mA)

DC 5V 1.2A (AC adapter 100~240V 50~60Hz 0.5A)

2 Each 1.5V D size battery (included)

Tilting Angle 0 to 90 degrees

• GENERAL

Power source

Dimensions At 0 degrees: 170mm (W) x 450mm (D) x 224mm (H)

At 90 degrees: 170mm (W) x 342mm (D) x 470mm (H)

Weight Approximately 14 lbs.

Color/Finish Blac

Standard Accessories Vinyl dust cover, 7mm nosepiece, prism compensator

Optional Accessories Manual marking system (for Lab use)

Storage Temperature 10C –40C, Humidity 30%–85%

Dust-free environment



