

## DELUXE 2 STANDS



*A new, innovative design complements the proven, 30-year mechanical operating system.*



*An electronic release button for slit lamp arm adjustments and a chair elevation switch are both conveniently located on the lower slit lamp arm.*



*An efficient counterbalanced refractor arm provides one-handed, immediate three-way locking of the refractor.*



*An attractive new overhead lamp features multiple adjustments, rotation stops and an on/off switch.*



*A full-featured flat overlay switch panel controls the slit lamp, chart projector, accessory instrument, overhead room lights, lamp brightness, chair elevation, and voltage for a corded instrument.*



*Three rechargeable wells and three 115V electrical receptacles.*



*Scuff-resistant, high-impact base cover is easily cleaned.*



*Standard Marco ten-step, hand-finished painting process.*

**DELUXE**

The industry's standard in efficient, instrument-delivery systems just got better. The new modern **DELUXE 2 STAND** clearly offers more pure quality and cosmetic appeal than any other instrument stand on the market. Innovative features such as an electronic slit lamp arm release button, overhead room light switch, and a laborious ten-step hand-finished painting process make the Deluxe 2 Stand a unique and differentiated product. While maintaining our usual high standards of function, quality and durability, the new Deluxe 2 Stand also remains competitively priced. For a compelling side-by-side demonstration against any other instrument stand, contact your local authorized Marco distributor.



*Deluxe 2 Lower Slit Lamp Arm*



*Deluxe 2 Stand full-featured instrument control panel, Stock #1206*



*Non-console Deluxe 2 Stand instrument control panel, Stock #1213*




Ask your Marco Distributor about The Marco Package Plus Program.



**PLAY IT SAFE.**

Protect your warranty. Buy your Marco product only from an authorized Marco distributor. Call Marco for the name of the distributor nearest you.

**MEASUREMENTS**



PO55000 BR-DLX2STD

**T H E L E A D E R I N V I S I O N D I A G N O S T I C S™**

