

VeraPlug™

Punctal Occlusion System

The **VeraPlug™** is designed to provide a simple and effective treatment for chronic dry eye.

- Conforms perfectly to the puncta for **unparalleled retention**
- Low-profile dome for **enhanced patient comfort**
- Proprietary shaft design for **easy insertion** and **proper anatomic fit**
- Available in four sizes



Choose the **sterile preloaded** configuration for the ultimate in convenience or **nonsterile bulk** for increased value.

SIZE CHART

| PUNCTAL OPENING | VERAPLUG™ SIZE |
|-----------------|----------------|
| 0.4mm to 0.6mm | Small |
| 0.6mm to 0.7mm | Medium |
| 0.7mm to 0.8mm | Large |
| 0.8mm to 1.0mm | X-Large |



AVAILABLE CONFIGURATIONS

Sterile Preloaded

Individually packaged and preloaded on a sterile, disposable inserter. 1 pair per box.

| SIZE | 1 PAIR |
|---------|----------|
| Small | VSP-1001 |
| Medium | VSP-1002 |
| Large | VSP-1003 |
| X-Large | VSP-1004 |

Nonsterile Bulk

Packaged nonsterile, 1 pair per poly bag, 10 pairs per box. Inserters included.

| SIZE | 10 PAIRS |
|---------|----------|
| Small | BNS-1001 |
| Medium | BNS-1002 |
| Large | BNS-1003 |
| X-Large | BNS-1004 |