

OS-100



OS-200



OS-300



OS-400



KEY FEATURES:

- » Four different stool styles to meet every need
 - » OS-100 Round stool with no back
 - » OS-200 Round stool with back
 - » OS-300 Contoured stool with back
 - » OS-400 Patient stool with auto locking wheels
- » Extended pneumatic elevation (8 inches)
- » Premium casters roll well on carpet and / or tile
- » Choice of any Topcon upholstery color

Performance & Specifications

Dimensions	OS-100	OS-200	OS-300	OS-400
Seat Size	15in (381mm)	15in (381mm)	16in x 15in (406mm x 381mm)	18in x 17in (457mm x 432mm)
Seat Foam Thickness	3in (76mm)	3in (76mm)	4in (102mm)	3in (76mm)
Retracted Seat Height	19in (483mm)	19in (483mm)	20in (507mm)	19in (483mm)
Extended Seat Height	27in (686mm)	27in (686mm)	28in (712mm)	27in (686mm)
Back Size	NA	10in x 15in (254mm x 381mm)	10in x 15in (254mm x 381mm)	12in x 17in (305mm x 432mm)
Base Footprint	24in x 24in (610mm x 610mm)	24in x 24in (610mm x 610mm)	24in x 24in (610mm x 610mm)	24in x 24in (610mm x 610mm)
Elevation Control	Side lever	Side lever	Side lever	Side lever
Caster Type	2in dual	2in dual	2.5in dual	2in dual/locking
Weight Capacity	250lbs (113.4kg)	250lbs (113.4kg)	250lbs (113.4kg)	250lbs (113.4kg)
Shipping Weight	21lbs (9.5kg)	28lbs (12.7kg)	31lbs (14kg)	31lbs (14kg)

Movements

Elevation Stroke	8in (203mm)	8in (203mm)	8in (203mm)	8in (203mm)
Back Adjustment Range	NA	3in (76mm)	3in (76mm)	3in (76mm)


Appearance

Upholstery Color	Your choice	Your choice	Your choice	Your choice
Cylinder Color	Black	Black	Chrome	Chrome
Base Color	Black	Black	Chrome	Chrome
Caster Color	Black	Black	Chrome	Chrome




Upholstery Selections

Standard Colors



Designer Series Colors



Topcon Medical Systems, Inc.
 111 Bauer Drive, Oakland, NJ 07436, U.S.A.
 t 800.223.1130 f 201.599.5248 • www.topconmedical.com
 © 2012 Topcon Medical Systems, Inc. MCA# 711