

JP90



Specifications

Yields approx. 30-120 cases per day

Dimensions: Height 6" | Width 2 1-4" | Depth 11"

Weight: 2 lbs | .90 kilos

Use with JP Bullets and JP Flask

Kits available for Unbreakable & iFlex

Patent pending

JP90 HANDHELD INJECTOR

saves time and money

- Injects: **Clasps, 1-3 tooth unilaterals, up to 5 consecutive tooth partial/saddle, saddles for combination cases, and repairs**
- Lightweight, compact, minimal maintenance
- Starter kits include 2 waxed cases for practice
- JP Bullets (cartridges) available in all shades of Unbreakable, iFlex, & Karadent
- **External heater: TCS Digital Furnace**
- 1 year warranty
- Made in the USA

KIT



Thermoplastic clasps



Saddles for combination cases



1-3 Tooth unilaterals



Up to 5 consecutive teeth partial saddles

STARTER KIT INCLUDES



1 JP90 Handheld Injector

TCS Digital Furnace included in item #s 7020-JPI & 7020-JPU



10 large JP Bullets



10 small JP Bullets



JP flask



TCS putty 400 gms



2 waxed cases for practice



5 metal sprues



1 matrix pusher



1 timer



866 • 426 • 2970 | 562-426-2970
www.tcsdentalinc.com

Item Codes	Description
7001-JPU	JP90 starter kit - Unbreakable™ bullets
7020-JPU	JP90 starter kit w/ digital furnace (110V) Unbreakable™ bullets
7001-JPI	JP90 starter kit - iFlex™ bullets
7020-JPI	JP90 starter kit w/ digital furnace (110V) iFlex™ bullets
7000-JP	JP90 handheld injector
7005-JP	JP flask