

Karadent™

Semi-rigid microcrystalline polymer for full dentures, semi flexible partial dentures, clear clasps, frameworks & nightguards




*Karadent™
clear framework partial*



MONOMER FREE &
HYPOALLERGENIC

NO MORE BROKEN APPLIANCES



Prescribe Karadent™ for a comfortable, aesthetic, break resistant, and durable appliance.

TIPS

1. For best results use lightweight alginate impression materials.
2. Adjusting Karadent™ is quick and easy. You can adjust using the same techniques as acrylic dentures.
3. Use a carbide bur for making small adjustments and Green Mounted Stones to smoothen area.
4. For a touch-up polish, Blue Hi-Shine polishing paste (item #4215-01) is recommended.
5. Recommend TCS Dental Appliance Cleaner to your patients for optimal maintenance.

Specify Karadent™ on your Lab Slip

Karadent™ by tcs®

Karadent™ by tcs™ is a monomer free, semi-rigid microcrystalline resin formulated for complete dentures, semi-flexible partials, clear clasps, clear frameworks and nightguards. Karadent™ combines the durable strength of nylons with the simplicity of acrylics to create a comfortable, aesthetic and reliable appliance.

Indications

- › Full dentures, partial dentures, clear clasps, frameworks & nightguards
- › Economical, aesthetic solution for missing teeth
- › Patients with concerns over metal restorations
- › Patients with allergies to monomer and/or nickel
- › Transitional restoration during the healing period for implants
- › Flexible partials when more stability is needed due to poor tissue support or few remaining teeth.
- › Flexible partials when more stability is needed for free-end saddles with teeth missing from the 1st bicuspid.

Benefits

- › Highly aesthetic
- › Maximum strength and durability
- › No volume shrinkage- exceptional fit
- › Minimal water absorption
- › Quick and easy chairside adjustments
- › Easy to repair, rebase, and relin with acrylic (laboratory procedure)

Insurance Codes:

D5110 Maxillary Complete Denture

D5120 Mandibular Complete Denture

Shades:



Standard Pink | Clear (clasps & framework)

