





Total Intravenous Anaesthesia Sets





High flow distal Codan valve (up to 280 mL/min) and anti-siphon valves (to prevent free-flow and back-flow) for the syringe pumps included



Microbore tubing used in all Tivatek™ kits ensures a low priming volume and reduces drug wastage



Anti-kink bridge as additional safety feature



High specification non-kink tubing offers a reduced risk of occlusion – even with several tight knots in the infusion line, drugs continue to run at flow rates set



2-part TIVA Infusion Sets with integral needle-free access ports and back check valve now available

Ordering information

Total Intravenous Anaesthesia Sets - clinical dead space 0.68 mL

Product code	Product description	Priming volume mL	Units per box
2-way			
TT202E	2-way 2.0m TIVA Infusion Set with high flow valve	2.2	50
TT252FE	2-way 2.5m 2-part TIVA Infusion Set with 2 needle-free access ports and valve	3.0	50
TT302E	2-way 3.0m TIVA Infusion Set with high flow valve	3.0	50
3-way			
TT223E	3-way 2.2m TIVA Infusion Set with high flow valve	4.3	50
TT223ER	3-way 2.2m TIVA Infusion Set with high flow valve and red connectors	4.3	50
TT253E	3-way 2.5m TIVA Infusion Set with high flow valve	4.9	50
TT253FE	3-way 2.5m 2-part TIVA Infusion Set with 2 needle-free access ports and valve	5.7	50
TT303E	3-way 3.0m TIVA Infusion Set with high flow valve	5.7	50
4-way			
TT224ER	4-way 2.2m TIVA Infusion Set with high flow valve and red connectors	6.4	25
TT304	4-way 3.0m TIVA Infusion Set with high flow valve and Tivalok™	8.3	25

Manufactured in the UK



Fannin (UK) Ltd.
Westminster Industrial Estate
Repton Road, Measham,
Swadlincote, Derbyshire,
DE12 7DT
T: +44 (0) 1530 514566
E: serviceuk@fannin.eu

Fannin Ltd.
Fannin House
South County Business Park
Leopardstown
Dublin 18
T: +353 (0) 1290 7000
E: info@fannin.eu