

Aston Shower Commode

A versatile solution for showering and toileting needs.



- Seat is moulded for comfort
- Armrests flip up for side transfers
- Backrest adjusts tool-free
- Simple e-clip height adjustment



Product Features

- Can be used solely as a commode where access to toilet is difficult
- Also assists in transporting users with restricted mobility from toilet to shower, or vice versa, to limit moving and handling transfers
- Self-propelling rear wheel kit available for easy attachment to enable independent use
- Polyurethane (PU) seat, seat-infill, back and armrests are water-impermeable, easy to clean and comfortable
- Bright blue upholstery provides strong colour contrast for visually impaired users
- Seat is moulded for comfort and warmth and has a front cut-away and rear recess to facilitate personal hygiene
- Lightweight, rustproof aluminium frame for ease of manoeuvring/handling and longevity
- Simple e-clip height adjustment to suit different toilet heights
- Tool-free height and angle adjustable footplates accommodate user support requirements
- Footplates swing in/out to facilitate transfers and are easily removable for ease of storage and transportation
- Armrests flip up for side transfers
- Clears standard and disabled WC's for optimum versatility
- Backrest adjusts tool-free in height to suit user and is easily detached/reassembled for storage and transportation
- Supplied with a 5 litre potty with lid, which is easily removable by sliding from rear of seat
- 4 substantial (12.7cm (5")) braked castors, manoeuvre easily over different floor surfaces and can be locked in any direction for improved safety
- Self-propelling wheels are quick-release to facilitate storage and transportation and have 61cm (24") grey, solid and puncture-proof tyres requiring minimum maintenance (self-propel kit sold separately)
- Rust-proof, easy-reach aluminium brakes on self-propelling wheels allow stable positioning and are durable
- Can be cleaned easily in domestic environment or thermally decontaminated in automated washing machines

CODE	DESCRIPTION
CWCO03	Aston Shower Commode
CWCO03PACK	Aston Self Propel Kit
CWCO03HD	Aston Heavy Duty Shower Commode

CODE	WIDTH	HEIGHT	DEPTH	SEAT WIDTH	SEAT HEIGHT	SEAT DEPTH	WEIGHT CAPACITY
CWCO03	56cm (22")	97 - 102cm (38" - 40")	90cm (35")	46cm (18")	51 - 56cm (20" - 22")	43cm (17")	159kg (25st)
CWCO03PACK	69cm (27")	97 - 102cm (38" - 40")	108cm (43")	46cm (18")	51 - 56cm (20" - 22")	43cm (17")	159kg (25st)
CWCO03HD	66cm (26")	97 - 102cm (38" - 40")	90cm (35")	56cm (22")	51 - 56cm (20" - 22")	43cm (17")	260kg (40st)