

## Sentra EC Wheelchair *Bariatric*

The Sentra EC Wheelchair has been designed for the larger person to ensure their mobility experience is as comfortable as possible.



- Sturdy, carbon steel frame with durable chrome for an attractive, chip proof, maintainable finish
- Dual axle position allows transition of seat height to hemi level
- Double embossed vinyl upholstery provides comfortable support and is durable, lightweight, attractive and easy to clean
- Simple folding mechanism facilitates storage and transportation
- Footrests can be height adjusted tool free to suit user's support requirements, can be swung out of the way to facilitate transfers or removed for storage/transportation
- Padded desk style armrests provide additional comfort and allow unhindered access to work surfaces
- Armrests are easily removable for sideways transfers, transportation and storage
- Solid tyres, mounted on composite mag rear wheels provide durability and low maintenance
- 20cm (8") solid front castors also ensure longevity and minimum maintenance and are height adjustable in 3 positions
- Maintenance-free, sealed, precision bearings
- Plastic hand rims enable easy manoeuvring for the independent user
- Attendant and user accessible parking brakes on rear wheels provide additional safety
- Elevating footrests are available for users who need to raise their legs, keep them straight or who require a more comfortable foot support
- Available in 3 different seat-widths to suit varying user sizes
- Extra wide seat-widths with reinforced framework are also available – 66cm (26"), 71cm (28"), 76cm (30") with a higher weight capacity of 318kg (50st) Please enquire for details

## Generous weight capacities

CODE	DESCRIPTION
STD20ECDDA-SF	Sentra EC Wheelchair 51cm (20") Self Propel
STD22ECDDA-SF	Sentra EC Wheelchair 56cm (22") Self Propel
STD24ECDDA-SF	Sentra EC Wheelchair 61cm (24") Self Propel
STD26DFA-SF	Sentra EC Heavy Duty Wheelchair 66cm (26") Self Propel
STD28DFA-SF	Sentra EC Heavy Duty Wheelchair 71cm (28") Self Propel
STD30DFA-SF	Sentra EC Heavy Duty Wheelchair 76cm (30") Self Propel



CODE	WIDTH (OPEN)	WIDTH (CLOSED)	LENGTH	HEIGHT	SEAT WIDTH	SEAT LENGTH	SEAT HEIGHT	WEIGHT	WEIGHT CAPACITY
STD20ECDDA-SF	75cm (30")	32cm (13")	107cm (42")	91cm (36")	51cm (20")	41cm (16")	50cm (20")	22kg (49lbs)	200kg (32st)
STD22ECDDA-SF	80cm (32")	34cm (14")	107cm (42")	91cm (36")	56cm (22")	46cm (18")	50cm (20")	26kg (57lb)	200kg (32st)
STD24ECDDA-SF	85cm (34")	37cm (15")	107cm (42")	91cm (36")	61cm (24")	46cm (18")	50cm (20")	26kg (57lb)	200kg (32st)
STD26DFA-SF	89cm (35")	37cm (15")	119cm (47")	97cm (38")	66cm (26")	50cm (20")	43cm-48cm (17"-19")	40kg (88lb)	318kg (50st)
STD28DFA-SF	94cm (37")	37cm (15")	119cm (47")	97cm (38")	71cm (28")	50cm (20")	43cm-48cm (17"-19")	40kg (88lb)	318kg (50st)
STD30DFA-SF	99cm (39")	37cm (15")	119cm (47")	97cm (38")	76cm (30")	50cm (20")	43cm-48cm (17"-19")	40kg (88lb)	318kg (50st)