





















Symbols Glossary

Lexington Intl., LLC
 1040 Holland Drive,
 Boca Raton, FL 33487
 US: 1 (800) 973-4769
 Intl: +1 (561) 417-0200

SYMBOL	SYMBOL TITLE	EXPLANATORY TEXT	STANDARD TITLE	STANDARD REFERENCE
	Manufacturer	Indicates the medical device manufacturer.	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	ISO 15223-1 Reference #5.1.1 FDA Recognition # 5-117 ISO 7000 Reference #3082 FDA Recognition # 5-103
	EC Rep	Indicates the Authorized representative in the European Community.	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	ISO 15223-1 Reference #5.1.2 FDA Recognition # 5-117
	Date of Manufacture	Indicates the date when the medical device was manufactured.	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	ISO 15223-1 Reference #5.1.3 FDA Recognition # 5-117 ISO 7000 Reference #2497 FDA Recognition # 5-103
	Serial #	Indicates the manufacturer's serial # so that a specific medical device can be identified.	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	ISO 15223-1 Reference #5.1.7 FDA Recognition # 5-117 ISO 7000 Reference #2498 FDA Recognition # 5-103
	Keep dry	Indicates a medical device that needs to be protected from moisture.	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	ISO 15223-1 Reference #5.3.4 FDA Recognition # 5-117 ISO 7000 Reference #0626 FDA Recognition # 5-103
	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed.	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	ISO 15223-1 Reference #5.3.7 FDA Recognition # 5-117 ISO 7000 Reference #0632 FDA Recognition # 5-103
	Humidity limitation	Indicates the range of humidity to which the medical device can be safely exposed.	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	ISO 15223-1 Reference #5.3.8 FDA Recognition # 5-117 ISO 7000 Reference #2620 FDA Recognition # 5-103
	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	ISO 15223-1 Reference #5.4.4 FDA Recognition # 5-117 ISO 7000 Reference #0434A FDA Recognition # 5-103
	CE Mark European Conformity Conformité Européenne	Signifies European conformity (CE) mark Indicates manufacturer declaration that the product complies with applicable European regulations	Guide to the implementation of directives based on new approach and global approach	NA
	General Warning Sign	Signifies a general warning	Graphical symbols - Safety colours and safety signs - Registered safety signs	ISO 7010 Reference # W001 FDA Recognition # 5-116

LGDXXXYYEXLSR

SYMBOL	SYMBOL TITLE	EXPLANATORY TEXT	STANDARD TITLE	STANDARD REFERENCE
	Refer to instruction manual/ booklet	Signifies that the instruction manual/booklet must be read	Graphical symbols - Safety colours and safety signs - Registered safety signs	ISO 7010 Reference #M002 FDA Recognition # 5-116
	WEEE	Signifies waste from electrical and electronic equipment	Waste Electrical and Electronic Equipment Directive	WEEE Directive 2002/96/EC
	Atmospheric Pressure Limitation	Indicates the range of atmospheric pressure to which the medical device can be safely exposed	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	ISO 15223-1 Reference # 5.3.9 FDA Recognition # 5-117
	RF Transmitter	Indicates a radio frequency is transmitted	Graphical Symbols for Use on Equipment	IEC 60417 Reference # 5140 FDA Recognition # 5-102
IPX2	Protected against vertically falling water drops	Protected against water drops up to 15 degree angle	Degrees of protection provided by enclosures (IP Code)	IEC 60529
	Type BF Part		Graphical Symbols for Use on Equipment	IEC 60417 Reference # 5333 FDA Recognition # 5-102
FC	FCC Compliance Mark	Complies with limits for Class B digital device established by FCC Rules, Part 15	NA	Title 47 United States Code of Federal Regulations Part 15.19
	RoHS	Signifies restricted use of ten Hazardous Materials in the manufacture of electronic and electrical equipment	Restriction of Hazardous Substances Directive	RoHS Directive 2002/95/EC
	Laser Warning	Signifies a Laser Warning	Graphical Symbols for Use on Equipment	ISO 15223-1 Reference #5.4.4 FDA Recognition # 5-117 ISO 7000 Reference #0434A FDA Recognition # 5-103
	Low-level ESD	To indicate packages containing electrostatic sensitive devices, or to identify a device or a connector that has not been tested for immunity to electrostatic discharge.	Electrostatic sensitive devices	IEC 60417
	Country of Manufacture	To identify the country of manufacture of products. In the application of this symbol, the "xx" shall be replaced by either the two letter country code or the three letter country code defined in ISO 3166-1	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	EN 60417-6049 ISO 3166-1 Recognition Part 1:Country codes

SYMBOL	SYMBOL TITLE	EXPLANATORY TEXT	STANDARD TITLE	STANDARD REFERENCE
 <p>The symbol consists of three parts: a yellow triangle with a black border and a black laser beam icon, the text 'CAUTION LASER 3R' in a black box, and a red circle with a diagonal slash over a black laser beam icon.</p>	<p>Caution Laser Radiation Class 3R</p>	<p>Signifies a Laser Radiation Class 3R Warning for Eye Hazards</p>	<p>Graphic Symbols for Use on Equipment</p>	<p>IEC 60825</p>