

Is your vet clinic using the best and most complete hemostatic technology to stop the bleeding and keep it stopped?

ClotIt® VET is scientifically superior and financially smarter for your clinic. Here are the FACTS:

Revolutionary Technology

Protégé Biomedical's newly patented all-natural Ionic Hemostatic Mineral (IHM) technology is unlike anything you've ever seen. It stops the bleed with an ionic attraction of blood cells and platelets and an activation of the animal's coagulation cascade. This multi-action process stops bleeding in seconds and keeps it stopped. Unlike the competition, it does not create a mechanical plug, it enables the animal's own body to stop the bleeding with a natural clot.

- **Scientifically proven to stop bleeding faster** – in lab and in surgery
- In Vitro - **16 seconds** in test tube clotting test vs. 4 minutes and 23 seconds with the leading competitor
- In Vivo - liver biopsy example: leading competitor stops it in 2:18 seconds, but it starts again after 30 seconds. Stopped again, starts again after 5 and 10 minutes. Continues to bleed.
ClotIt VET: Stopped and KEPT stopped after 46 seconds¹

Total Clinic Approach

Not ONLY a surgical product. Protégé Biomedical offers custom products at accessible price points for surgery, daily clinic procedures (like nail clipping and flesh bleeding), and co-branded pet first aid kits with home-use ClotIt for emergencies and minor bleeding at home or on the go. Get the benefit of a complete line of both sterile and non-sterile products with ClotIt VET.


Guaranteed Satisfaction

100% money back guarantee. If you don't love our product for any reason, you can return it at any time for a full refund.

Liver Biopsy Clinical Study Results

ClotIt VET

Stopped bleeding after :43 seconds and **STAYED STOPPED.**




Hemablock®


Stopped bleeding after 2:18 seconds, then started again after 5 minutes, then stopped again with more treatment, then started again after 10 minutes.



ALL-NATURAL IONIC HEMOSTATIC MINERAL (IHM) TECHNOLOGY



When applied to a bleed site, IHM absorbs plasma and ionically attracts cells to begin the formation of a clot.



Factor XII initiation triggers the body's own coagulation process to build a fibrin web, stopping the bleed in seconds.

¹ Data on file.

Superior Science. Superior Clotting.

Surgical Hemostats Comparison Presented Clearly and Honestly

ClotIt VET



ClotIt



	ClotIt	Hemablock®	Vetigel™	Gelfoam®	Yunnan Baiyao	Vetspon®
Source	Mineral Based ¹	Potatoes & Cotton MPB Technology	Polymer Gel ⁴	Porcine Skin ¹¹	Chinese Herbal Blend ¹⁸	Gelatin ²⁴
Bio-Absorbable	YES time not specified ²	YES 1 week or less	YES time not specified ⁵	YES 4-6 weeks ¹²	NO not for internal use ¹⁹	YES 4-6 weeks ²⁵
Time to Hemostasis	20 seconds ³	1-2 minutes	15 seconds ⁶	Not Cited ¹³	Not Cited ²⁰	3-5 minutes ²⁶
Indication	Minor to severe bleeding ³	Profuse & diffuse specifically designed to be left in the body post-surgically	Moderate to severe bleeding ⁷	Surgical bleeding ¹⁴	Traumatic injury or internal bleeding ²¹	Surgical procedures when control of bleeding by pressure, ligature or other conventional procedures is ineffective or impractical ²⁷
Mode of Action	Factor XII activation flocculation ³	Molecular sieve	Mechanical barrier ⁸	Not fully understood ¹⁵	Release clotting granules ²²	No intrinsic hemostatic action Does not activate clotting cascade ²⁸
Usage Options	2 Delivery Options: A) long tip puffer bottle B) short tip puffer bottle	2 Delivery Options: A) sterile syringe applicator B) sterile foil pack inside resealable tube	Sterile hemostatic gel kit with a syringe ⁹	Sterile sponge ¹⁶	2 Delivery Options: A) ingestible pill B) powder ²³	2 Delivery Options: A) dental cube B) flex sheet ²⁹
Multi-Use	YES Sterile on initial use additional uses non-sterile*	YES Sterile on initial use additional uses non-sterile	NO ¹⁰	NO ¹⁷	YES Sterile on initial use additional uses non-sterile	Can be cut to desired size and shape ³⁰

* Note: ClotIt VET is available in sterile, single-dose applications for use in the OR. If reused, additional uses would be non-sterile.

ClotItVet.com

ClotIt VET

- <http://www.protegebiomedical.com/technology>
- <http://www.medicalplasticsnews.com/news/technology/suneris-have-invented-a-gel-that-stops-bleeding-in-seconds/>
- Protégé Biomedical, LLC Study Summary, PRO-008(b)-May 2015

Vetigel

- <http://www.i-runway.com/blog/magic-gel-heals-wounds-in-seconds/>
- <https://patents.google.com/patent/US20140287061A1/en>
- <https://cresilon.com/index.php/our-story/>
- <http://www.industrytap.com/plant-based-vetigel-stops-bleeding-instantly/29868>
- <https://www.youtube.com/watch?v=djLxRcU9No4>
- <http://www.i-runway.com/blog/magic-gel-heals-wounds-in-seconds/>
- <http://www.ibtimes.co.uk/vetigel-plant-based-gel-that-stops-traumatic-bleeding-wounds-15-seconds-1476464>
- <http://www.medicaldaily.com/pulse/vetigel-stops-bleeding-seconds-polymer-could-be-band-aid-future-311514>
- <https://cresilon.com/index.php/vetigel/>
- <https://www.acs.org/content/dam/acsorg/events/technology-innovation/Slides/2016-11-03-hemostatic-gel-slides.pdf> - Slide 8
- <https://cresilon.com/index.php/vetigel/>
- <https://cresilon.com/index.php/vetigel/>

Gelfoam

- http://www.pfizer.com/files/products/uspi_gelfoam_plus.pdf
- http://www.pfizer.com/files/products/uspi_gelfoam_plus.pdf
- <http://labeling.pfizer.com/ShowLabeling.aspx?id=624>
- http://www.pfizer.com/files/products/uspi_gelfoam_plus.pdf
- http://www.pfizer.com/files/products/uspi_gelfoam_plus.pdf
- https://www.pattersonvet.com/Supplies/ProductFamilyDetails/PIF_352922
- http://www.pfizer.com/files/products/uspi_gelfoam_plus.pdf

Yunnan Baiyao

- <https://www.activeherb.com/baiyao/>
- <http://www.veterinaryplace.com/dog-medicine/yunnan-baiyao-for-dogs/>
- <https://pathwithpaws.com/blog/2011/06/11/a-magic-vial-of-yunnan-baiyao/>
- <https://www.activeherb.com/baiyao/>
- <http://wendyblount.com/lufonco/32-yunnan-paiyao.pdf>
- <https://yunnan-baiyao-store.com/blogs/how-can-my-dog-benefit-from-yunnan-baiyao/45054465-how-can-yunnan-baiyao-capsules-help-with-your-sick-dog>
- <https://pathwithpaws.com/blog/2011/06/11/a-magic-vial-of-yunnan-baiyao/>

Vetspon

- <http://slideplayer.com/slide/5838/> - Slide 1
- <http://slideplayer.com/slide/5838/> - Slide 6
- <https://www.google.com/patents/US8697118>
- <http://slideplayer.com/slide/5838/> - Slide 2
- <http://slideplayer.com/slide/5838/> - Slide 1 and 2
- <http://slideplayer.com/slide/5838/> - Slide 7
- <http://slideplayer.com/slide/5838/> - Slide 2

Now that you know the FACTS, make **ClotIt® VET** your hemostatic brand of choice.

www.protegebiomedical.com
www.clotitvet.com

