

CUBRO

bock[®] v2 electric bed and Posture-Temp[™] mattress package



The electronic 4 section bock[®] v2 bed comes in 2 sizes and a choice of 2 endboard colours. Coupled with the Posture-Temp[™] mattress, this package provides comfort for clients, is easy-to-use and durable.

A range of accessories are also available.

Get in touch

For more information see
contact us on 0800 656 527
or email sales@cubro.co.nz
or visit cubro.co.nz

.bock[™] **POSTURE-TEMP[™]** **CUBRO**

bock® v2 electric bed and Posture-Temp™ mattress package

.bock



Made in Germany

bock® v2 4-section electric bed

The bock® v2 is a 4-section electric bed that combines maximum comfort and ease of use. The professional design and economy of this bed is great for daily care in the home.

Combine good looks with durability for a tough, scratch-resistant prime melamine surface on the end boards. Available in either Oyster line velvet or Dark oak wood finish.

A range of optional extras are available.



4-section bed movement



Ergo Plus hand controller is designed for ease of use with fully adjustable mounting clip to allow multi positioning. The integrated disabler function allows all electric functions to be disabled with a simple hand controller lock.



Dark oak endboards and Oyster line endboards

Posture-Temp™ mattress

The Posture-Temp™ mattress is made from unique temperature-sensitive, ultra-soft layers of open-cell, slow recovery, visco-memory foam. It warms to conform to the user's shape, maximising the body's contact area and thereby facilitating excellent comfort and support with great pressure redistribution.

POSTURE-TEMP™
CUBRO



Bock® v2 electric bed with Dark oak endboards



Posture-Temp™ mattress with visco-memory foam layer



CUBRO

sales@cubro.co.nz | 0800 656 527
cubro.co.nz

bock® v2 electric bed and Posture-Temp™ mattress package

.bock

POSTURE-TEMP™
CUBRO

Features and specifications

Package features

- Includes: 4-section electric bed and Posture-Temp™ Pressure relieving mattress
- Sprung slat base for pressure care and ventilation
- Main's isolation prevents the generation of electric fields
- All function hand controller with integrated disabler control and adjustable rail clamp
- High grade coiled cable with strain relief and anti-kink feature to protect against cable damage
- Automatic triple action knee break helps to reduce the pressure on heels
- Battery for emergency lowering
- The Super Single bed option has inbuilt compensation to reduce potential sliding



Package specifications

Code	8610	8600	8660	8650
Size	Single	Single	Super Single	Super Single
Mattress dimensions	W 900mm	900mm	1100mm	1100mm
	L 2000mm	2000mm	2000mm	2000mm
	H 150mm	150mm	150mm	150mm
Overall dimensions	W 915mm	915mm	1220mm	1220mm
	L 2160mm	2160mm	2160mm	2160mm
Endboard colour	Dark oak	Oyster line	Dark oak	Oyster line
Maximum user weight	165kg			
Warranty	Bed electronics 2 years, frame 5 years			



Optional accessories:



bock® bed extension kit
5066, 5067



Fold downside rails left/right
910T-L, 910T-R



bock® Self-help pole
5113



Clamp-on side support rail
5109-2

CUBRO

sales@cupro.co.nz | 0800 656 527
cupro.co.nz