

rollz®
takes you further



rollz

Ergonomically and stylishly designed rollators

Life is an endless discovery journey for all of us. Living life to the full, that's what Rollz is all about. Be it for a stroll in the park, a visit to an exhibition in a museum, a holiday or a shopping afternoon. Our ergonomically and stylishly designed Rollz Motion and Rollz Flex will bring you further – wherever and whenever you feel like going.

...



...

Out and about with the Rollz Motion²

The 2-in-1 rollator and wheelchair

Going out and about with the Rollz Motion without a care in the world. Walking with the rollator or sitting in the wheelchair, enjoying a lovely day in the city. Let's stroll through the old town, rest in the wheelchair while enjoying the lush park, and then walk again to that nice restaurant to end this interesting and inspiring day.



rollz motion

Product features Rollz Motion²

Stylish

Award winning Dutch design rollator and wheelchair 2-in-1 for going out in style.

Safe

Safe and controlled braking in all weather with the drum brakes.

Convenient

Resting on a comfortable cushioned seat or convert the rollator to a wheelchair.

Easy aid

Cross doorsteps and curbs easily with the threshold aid.

Stable

Walking tall, steadily and supported between the ergonomic hand grips that are adjustable in height.

Compact


The Rollz Motion is easily foldable into a compact package and fits in almost any car.

Comfortable ride

Natural shock absorption thanks to the soft tyres and foam hand grips.

Manoeuvrable

Light steering and easy motion with the manoeuvrable front wheels.



With a flick of the wrist

You can easily transform this rollator into a wheelchair and back again. Walk steadily, upright and comfortably, and rest along the way in the wheelchair if you feel like it.

...



rollz motion

Specifications

Rollz Motion²: 10.7 kg

Rollz Wheelchair package: 3.7 kg

Maximum carrying load: 125 kg / 19 st 9 lb

Minimum and maximum height user:

1.60-1.90 m / 5 ft 3" to 6 ft 3"

Unfolded H×W×L: 87-102 × 65.5 × 63 cm

Folded H×WS×L: 90-102 × 31 × 63 cm

Seat depth: 42 cm

Seat width: 45 cm

Seat height: 55 cm

Front wheel diameter: 8"

Rear wheel diameter: 12"

Warranty: 5 years on the frame and 2 years on other parts with exception of normal wear and tear.

...



Rollz Motion² Small: 10.7 kg

Rollz Wheelchair package: 3.7 kg

Maximum carrying load: 125 kg / 19 st 9 lb

Minimum and maximum height user:

1.50-1.85 m / 4 ft 11" to 6 ft 7/8"

Unfolded H×B×L: 82-95 × 65.5 × 63 cm

Folded H×B×L: 85-95 × 31 × 63 cm

Seat depth: 42 cm

Seat width: 45 cm

Seat height: 50 cm

Front wheel diameter: 8"

Rear wheel diameter: 12"

The Rollz Motion² is available in four colours: **dark purple, pebble white, island blue** and **matt black**.



“Now I dare to walk longer with the rollator, because I know that when I become tired, I can rest in the wheelchair on the go.”



“I could not believe that such a light rollator can give you so much stability and convenience.”

Out and about with the Rollz Flex

The shopper rollator

There's a clear reason why we call the Rollz Flex our shopper rollator. Whether it is the weekly market in town, getting your groceries for the weekend at the supermarket, or a nice afternoon shopping in the city, the Rollz Flex turns it into a light and comfortable experience. Walking steadily and conveniently, while the bag slowly carries more and more.



rollz flex

rollz

Product features Rollz Flex

Safe

Choice of flip brake and classic brake version for safe drum braking in any weather type.

Practical

The large water resistant bag can carry 20 kg. Including firm and removable inner bag with pocket.

Comfortable

Comfortable ride and natural suspension thanks to the soft tyres.

Assistant

Cross doorsteps and curbs easily with the threshold aid.

Stylish

Fully certified rollator with the look of a designer shopper.

Flexible

Adjustable handlebar that is not only adjustable in height, but also in three different positions. Photo shows rollator position.

Convenient

Resting on the luxurious artificial leather seat with the handlebar in the forward position as a supportive backrest.

Light and compact

The Rollz Flex weighs only 7.4 kg and is easily foldable into a compact package.

Manoeuvrable

Very light steering and easy motion with the manoeuvrable front wheels.



...

rollz flex

Multiple ways to use

The adjustable handlebar caters to your needs as regards stability and convenience: in **upright position** or **facing forward** for even more stability. The large bag is your shopping companion.

...



Details

Rollz Flex: 7.4 kg

Rollz Flex bag: 25 liter / 5.5 UK gal / 20 kg

Maximum carrying load: 125 k / 19 st 9 lb

Minimum and maximum height user:

1,55-1,95 m / 5 ft 1" - 6 ft 4 3/4"

Unfolded h×w×l: 80-100 × 62 × 60 cm

Folded h×w×l: 80-100 × 33 × 60 cm

Seat width: 44 cm

Seat height: 62 cm

Diameter wheel: 8"

Warranty:

5 years on the frame and 2 years on other parts with exception of normal wear and tear.

In rollator-position the Rollz Flex is a fully TÜV certified rollator.

...

Rollz Flex Small: 6.8 kg

Rollz Flex bag: 14 liter / 3.1 UK gal / 20 kg

Maximum carrying load: 125 kg / 19 st 9 lb

Minimum and maximum height user:

1,40-1,60 m / 4 ft 7 1/8"

Unfolded h×w×l: 70-81 × 62 × 60 cm

Folded h×w×l: 70-81 × 33 × 60 cm

Seat width: 44 cm

Seat height: 52.5 cm

Diameter wheels: 8"

The Rollz Flex bag is made out of sustainable and waterproof material and is available in two colours: **bright purple** and **denim grey**.





Use your Rollz Rollator to the fullest with our accessories.

Rollz Motion Travel cover

Due to the unique design of the Rollz travel cover, both the Rollz Motion and the wheelchair package are incorporated into one compact package.



Rollz Motion Shopper

This stylish and functional shopping bag can be clicked onto the Rollz Motion frame. You can put your shopping directly in the shopper and take it along, so you won't have to carry a shopping bag.




Rollz Motion Cane Holder

If you also use a cane, you can bring it along easily with the Rollz cane holder. Because of the clever design nearly every cane fits.



Rollz Motion Cup Holder

The Rollz cup holder is very easy to use. Because of the sliding system, the Rollz cup holder can be placed into the desired position quickly, enabling the user to take a water bottle along.



Rollz Motion Wheelchair package holder

With these holders you can also hang your wheelchair package over the seat. This gives you more depth and freedom of movement. A stick- and shopping bag holder is integrated.



rollz accessories

Rollz Motion Back support

Out on your own and left the wheelchair package at home? This back support offers you a short rest, while leaning against the soft strap.



Rollz Motion Footrests in different heights

Choose the size footrest that suits you for even more comfort. The small version is 5 cm shorter, the large version is 5 cm longer.



Rollz Flex Cane Holder

The cane holder simply slides over the handlebar. Do you need your walking stick somewhere? Then you easily slide it from its holder. The clever design fits almost any walking stick in the cane holder.



Rollz Flex Shopping hooks

If you go shopping, these hooks offer the opportunity to add another bag.



Manufacturer

Rollz International BV

info@rollz.com

www.rollz.com

rollz®