

Economical Bariatric Wheelchairs



*230 kg plus
User Weight
Limit*

*Transit
and Self
Propelling
Option*

radius
shop



AFFORDABLE BARIATRIC MOBILITY OPTIONS

These steel frame wheelchairs offer a robust design to meet rigorous demands of Community, Hire/Loan and Healthcare environments. They fold easily for storage or transportation and come with removable armrests and swing-away leg rests to assist with getting in and out of the chair - a great benefit for the user and care givers.

FEATURES

- Steel frames: Transit powder coated and Self Propelling chrome finished
- Detachable desk armrests
- Comfortable padded upholstery: Nylon covered on Transit and easy clean vinyl on Self Propelling
- Height adjustable swing-away footrests
- Folds for easy storage and transport
- 8" x 2" front castors with PU tyres
- Self Propelling rear wheel: 24" mag with PU tyres
- Transit rear wheel: 12" with PU tyres

OPTIONAL ACCESSORIES

- Anti-tippers
- Elevating leg rests


SPECIFICATIONS



Seat Width
Transit 22"
Self Propel
20", 22", 24"



Seat Depth
18"



Overall Width
29"



Overall Length
42"



Overall Height
36"



Seat to Floor
Height
19.5"




Back Height
18"



Rear Wheel
Transit 12"
Self Propel 24"



Front Castor
8" x 2" PU



Folded Width
13.5"



Total Weight*
Transit 31.6 kg
Self Propel
30 kg



User Weight
Limit
Transit 230 kg
Self Propel
250 kg

* Including Footrests.

This information is correct at the time of printing. From time to time product specifications are changed by the manufacturers without prior notice.

0800 213 313

www.radiusshop.co.nz