

AeroTrach[®] Plus^{*}

Valved Holding Chamber



Delivery and relief for your special needs patients.

Keeping the standard.


At Trudell Medical International we set the gold standard for the Valved Holding Chamber (VHC). Now we redefine that standard for special-need tracheostomy patients with the **AeroTrach Plus^{*}** VHC – a simple and cost-effective system to quickly deliver medicine where it's needed most – to your patient's lungs.

Quick relief.

The **AeroTrach Plus^{*}** VHC aerosol delivery is effective and quick – two important factors for any tracheostomy patient. Medication concentrations are delivered on inspiration only and are consistent – a benefit of the enhanced particle suspension time of the **AeroTrach Plus^{*}** VHC.

Expect innovation.

Trudell Medical International is renowned for its extensive research and development of new technology for aerosol delivery. We are proud to continue the tradition of innovation with the **AeroTrach Plus^{*}** VHC, a significant improvement in asthma relief for tracheostomy patients.




Trudell Medical International[®]

www.trudellmed.com



AeroTrach
Plus^{*}
Valved Holding Chamber



Flow Dynamic Design:
Improves respirable dose of medication. Allows exhalation away from the patient.

More efficient and comfortable medication delivery for your special needs patients.

Comfort Right Fit* Adapter

Typical 15 mm adapters are designed to connect firmly in place. The **Comfort Right Fit*** adapter means the connection will fit tightly for enhanced aerosol delivery, without excessive torque when connecting or disconnecting the device.

Inhalation Valve

With a new sensitive design for low resistance and superior responsiveness, even a low inspiratory effort will be rewarded.

Exhalation Valve

With our improved, more robust design, your patients will experience a low resistance to exhaled flow. The **Comfort Right Fit*** adapter guards against valve tampering or damage.

Chamber Design

Our new flow dynamic chamber design maximizes fine particle dosage and allows for enhanced suspension of medication in the chamber.

MDI Adapter

We've selected a new, improved polymer material for the universal MDI adapter. With this reinforced material, we've increased product durability.

Unique, clear, **Comfort Right Fit*** adapter ensures visibility

Valve and baffle system maximizes the delivery of fine particles emitted from standard pressurized MDIs



Strong polymer material for durability

Aerodynamic chamber

Universal MDI Adapter



Trudell Medical International^{*}

Trudell Medical International
725 Third Street, London,
Ontario, Canada, N5V 5G4
Phone: +1-519-455-7060
Fax: +1-519-455-7858
www.trudellmed.com

Authorised European Representative:
Trudell Medical International Europe Ltd.
BioCity Nottingham, Pennyfoot Street
Nottingham, NG1 1GF, UK
Tel: +44 115 912 4380
Fax: +44 115 912 4289

CE REGISTERED
ISO 13485
QUALITY SYSTEM[™]