

# ef IOW "rapid

Fast Nebulisation – Short inhalation time				
Drug	Fill volume	MMD	TOR	Inhalation time
	[ml]	[µm]	[g/min]	[min]
Tobramycin 300 mg	5.0	3.9	0.50	6.8
Dornase alfa 2.5 mg	2.5	3.9	0.61	2.3
Salbutamol 2.5 mg	2.5	4.3	0.68	2.1

# Worldwide mobility – small, light, silent

- Battery operation
- Rechargeable battery operation including fast recharge
- Mains adapter with 3 plugs for UK, EU, US
- Approx. weight 300 g, including rechargeable batteries
- Silent and discreet operation

# For nebuliser hygiene\*

- Can be disinfected in a vaporiser
- Can be disinfected in boiling water
- Autoclavable

\* Please apply the notes on hygiene in the actual instructions for use.

eFlow <sup>®</sup> rapid			
Aerosol Head			
Nebuliser Handset			

Order numbers.: **178G1005** 178B2603 178G8012

PARI Pharma Advancing Aerosol Therapies



### Technical data:

- Principle of nebulisation: perforated vibrating membrane
- Vibration frequency: 117 kHz
- Automatic switch off
- Battery or rechargeable battery operation: At least 90 minutes operation with fully charged rechargeable batteries
- Mains adapter: 100 - 240V~, 50/60 Hz, 2 W
- Nebuliser weight: ~ 55 g
- Total weight: ~ 300 g
- Min. fill volume: 2.0 ml
- Max. fill volume: 6.0 ml

The eFlow®*rapid* has been designed for use with all inhalation medications approved in Europe.

## More information

PARI GmbH Specialists in effective inhalation Moosstraße 3 • D- 82319 Starnberg

**PARI** Export Service Center: Tel.: 0049 (0) 8151/279-220 E-mail: info@pari.de

www.pari.de www.eFlowrapid.info

Local distributor:

06/05

© 2005 PARI GmbH · Spezialisten für effektive Inhalation · 178D0015 ·