

INSTRUCTIONS FOR USE



23111, 23113, ADJUSTABLE MOBILE SHIELD / OPTIONAL CATCHING TRAY

1) Indication for use:
Designed to
protect the
patient against
scattered
radiation, the Wolf
Mobile Shield is
ideal for chest,
skull and spine
radiography. The



screen is made of lead vinyl with a 0.5mm lead equivalency protection and is attractively covered on both sides by easy to maintain vinyl. It is 24" (60cm) wide and is adjustable from a minimum height of 33" (84cm) to a maximum height of 48"

ADJUSTABLE MOBILE SHIELD / OPTIONAL CATCHING TRAY

<u>Instructions for Use</u>

Note: It is always recommended to

- A. It's mounted on a durable chrome finished frame and rolls easily to any location on 4 conductive casters.
- B. The steel Mobile Shield Tray easily attaches (no tools necessary) to the frame of the Adjustable Mobile Shield and will help to keep your Adjustable Mobile Shield looking tidy.
- C. Turn knobs to desired height.
- 2) Routine care and cleaning instructions: Wipe with an antibacterial wipe
- 3) Precautions/Warnings: California's Proposition 65
- 4) <u>Disposal:</u> Product should be disposed of according to applicable local laws.

rev 02, 23111, 23113 8/2021