



INSTRUCTIONS FOR USE

Item # 61100: Apron Caddy

REF



- 1) Indication for use:
Wolf's Apron Caddy is designed to help transport and protect your apron. By rolling your apron and inserting it in the Caddy's protective tube, you avoid creasing the apron's internal x-ray protective lining, and potentially decreasing its usable life. Includes shoulder strap, carry handle and a handy outside pocket. Holds one adult apron. Available in Navy Blue.
- 2) Instructions for Use: Roll apron and insert it in the caddy protective tube. Hold's one adult apron.
- 3) Routine care and cleaning instructions: Spot clean
- 4) Inspection:
- 5) Precautions/Warnings:
- 6) Disposal: Product should be disposed of according to applicable local laws.

rev 01, 61100
5/2021