

INSTRUCTIONS FOR USE 80-038, 80-037, 80-124 CR/DR PLATE CARTS





1) Indication for use:

Carts are available in "Tall," "Stationary," and "Mini." They look exactly like the cassette carts except for the addition of the 3/8" padding.

a. Tall Cassette Cart

Eight 6" wide compartments with $16\frac{3}{4}$ " vertical clearance, 24" L x 18" W x $42\frac{1}{2}$ " H. Weight 120 lbs. Item # 80-038

- b. Stationary version of the tall Cassette Cart. (no casters). Item # 80-037
- c. Mini Cassette Cart

Four 6" wide compartments with $16\frac{3}{4}$ " vertical clearance, 24" L x 18" W x 36" H. Weight 95 lbs. Item # 80-024

CR/DR PLATE CARTS

2. Instructions for Use

Note: It is always recommended to

Store CR/DR plates in cart. Opening on both sides. Place CR/DR plate in compartment.

- 3) Routine care and cleaning instructions: Wipe with an antibacterial wipe
- 4) <u>Precautions/Warnings:</u> Like any device, these casters must be used with caution. The casters provide an extra measure of safety, but user vigilance is always required.
- 5) Disposal: Product should be disposed of according to applicable local laws.

rev 02, 80-038, 80-037, 80-124 8/2021