



INSTRUCTIONS FOR USE
80-046, 80-045 TOP LOADING CR/DR CARTS

REF



1) Indication for use:

- a) Like the "Tall" and "Mini" CR Carts, our Top Loading Cards feature $\frac{3}{8}$ " high impact padding on surfaces that could damage CR Plates.

Available in 2 compartments lengthwise or 4 compartments widthwise. Overall cart size 24" L x 16" W x 26" H 90 lbs.

1. Two compartments 8" W x 24" L x 13" D Item # 80-046

Four compartments 6" W x 16" L x 13" D Item # 80-045

TOP LOADING CR/DR PLATE CARTS

2. Instructions for Use

Note: It is always recommended to

Store CR/DR plates in cart by opening cover and place CR/DR plate in compartment.

- 3) Routine care and cleaning instructions: Wipe with an antibacterial wipe
- 4) Precautions/Warnings: Like any device, these casters must be used with caution. The casters provide an extra measure of safety but user vigilance is always required.
- 5) Disposal: Product should be disposed of according to applicable local laws.

rev 02, 80-046, 80-045
8/2021