

INSTRUCTIONS FOR USE 56147, 56148, 56150, 56151, 56152, 56153, DR PANEL PROTECTOR

REF



1) Indication for use:

Faceplate made of polycarbonate supporting up to 500 and X-ray lucent for high quality image contrast. Rubber floor grips prevents DR Panel Protector from slipping. Made of germ-resistant, nonporous, ease to maintain and clean polyethylene, polycarbonate and powder coated aluminum. Compatible with Cannon, Siemens, Toshiba, Phillips, Shimadzu, Quantum, Carestream, Fuji, Varian, Konica, Minolta, and AGFA digital systems.

DR PANEL

2) <u>Instructions for Use</u>

Note: It is always recommended to

- A. Simply place over the DR panel positioned on floor.
- B. Fit cable in cut out
- C. Patient steps on the DR-panel
- 3) Routine care and cleaning instructions: Wipe with an antibacterial wipe
- 4) Precautions/Warnings: California's Proposition 65
- 5) <u>Disposal:</u> Product should be disposed of according to applicable local laws.

rev 02, 56147, 56148, 56150, 56151, 56152, 56153 8/2021