

INSTRUCTIONS FOR USE

REF

50175, 50176, 50177: Digital Imaging Markers









Style	Stock Number
Digital Style R & L ¾" 2 or 3 initials ¼"	50175
Digital Style R & L ½" 2 or 3 initials ¼"	50176
Digital Style Ruler only	50177

1) Indication for use:

These markers are specifically designed for use in Teleradiology and CR imaging where the image size does not always match the actual size anatomy.

Digital Imaging Markers:

- 2) <u>Instructions for Use</u>
 - a. The image of the 3cm ruler allows you to compare the screen image to real life.
 - b. Color coded R-red, L-blue, characters are high-density lead.
- 3) Routine care and cleaning instructions: Clean with disinfectant wipes.
- 4) Precautions/Warnings: California Proposition 65
- 5) <u>Disposal:</u> Product should be disposed of according to applicable local laws.
- 6) Adverse reactions: non known

rev 02, 50175, 50176, 50177 8/2021