




INSTRUCTIONS FOR USE

REF

50160, 50161, 50162, 50163, 50164, 50165, 50166, 50167, 50167A: “Elite” Cassette Markers



1) Indication for use:

ELITE markers are color coded: R-Red / L-Blue and feature encased high density lead letters in aluminum casing for durability. Available with or without initials or numbers. Red is color for "R," blue is the color for "L." Also available in yellow and green. Markers are sealed and words are formed from individual letters.

“Elite” Lead Markers:

2) Instructions for Use

Coded red for "right," blue for "left" they are available with or without initials or numbers.

3) Routine care and cleaning instructions: Clean with disinfectant wipes.

4) Precautions/Warnings: California Proposition 65

5) Disposal: Product should be disposed of according to applicable local laws.

6) Adverse reactions: non known

rev 02, 50160, 50161, 50162, 50163, 50164, 50165, 50166, 50167, 50167A
8/2021