



INSTRUCTIONS FOR USE

REF

**21609, 21610, 216090, 2161000 WOLF XRAY UTILITY FILE
21611C, 216110, 21612C, 216120 FILE HANGING SYSTEM
25125, 25126 FILE STAK
1025118 FLOOR STANDING FILM BIN**



1) Indication for use:

- a. Utility File: is available with an opaque or translucent plexi front and in two sizes. Will hold your largest cassettes, films or jackets.
- b. File Hanging System Available in letter or legal size and your choice of translucent or opaque plexi front.
- c. File-Stack takes only 19½" x 32" of valuable floor space, stacks as high as you want. Made of sturdy 20-gauge steel.
- d. Floor Standing Film Bin: This light-proof cabinet is designed for film storage of sizes 5" X 7" to 14" X 17". Bin height is ideal for placement under a countertop.

Instructions for Use

Note: It is always recommended to

- a. Utility File:
 - i. Wall mounted via 2 keyholes.
 - ii. Place large cassettes, films or jackets
- b. File Hanging System
 - i. Wall mounted via 4 keyholes.
- c. File-Stack
 - i. Individual units are secured to each other by bolts and come with 3 adjustable dividers.
 - ii. Paneled doors slide easily at the touch of a finger for ready access. Overall finish is in attractive
 - iii. Easy to clean white baked on enamel.



- d. Floor Standing Film Bin:
 - i. It's hinged at the bottom, providing easy access without bending over.
 - ii. The non-moveable dividers provide an easy filing system.
 - iii. An optional switch is available for automatic room light cut-off when the door is opened.
 - iv. Holds 4 boxes (400 pieces) of film.
- 3) Routine care and cleaning instructions: Wipe with an antibacterial wipe
- 4) Precautions/Warnings:
 - a. Floor Standing Film Bin - Do not open in light
- 5) Disposal: Product should be disposed of according to applicable local laws.

rev 01, 21609, 21610, 21609O, 216100O, 21611C, 216110, 21612C, 216120, 25125, 25126, 25118
8/2021