




INSTRUCTIONS FOR USE

REF 56160 FOLD & STORE CD/DR TABLE TOP CASSETTE HOLDER



1) Indication for use:

Cross table device that holds all standard CR and DR plates from 8½ - 19½" width. Easily transforms from 10.5" to 3". Has a spring-loaded release button. Easy carry and storage handle. Can be stowed in a drawer or hung on the wall.

Fold & Store

2) Instructions for Use

Note: It is always recommended to

- a) Untwist lock knob to slide in cassette
- b) To open and close- utilize spring loaded release slider
- c) To carry must close unit and then use carrying handle

3) Routine care and cleaning instructions: Wipe with an antibacterial wipe

4) Precautions/Warnings:

5) Disposal: Product should be disposed of according to applicable local laws.

rev 01, 56160
8/2021