

## INSTRUCTIONS FOR USE 25515 WOLF HEAVY DUTY CASSETTE CART

REF



## 1) Indication for use:

Designed for rugged use, the heavy duty cassette cart is constructed of heavy gauge steel, with a tough washable white finish. The cart rolls about easily on 4 sturdy rubber casters that lock in position. This model has 3 levels divided into 5 compartments that hold any size cassette from 14" X 17" on down (maximum 125 cassettes).

50" H x 17" D x 20" W. 100 lbs.

## 2. <u>Instructions for Use</u>

Note: It is always recommended to

Store CR/DR plates in cart. Place CR/DR plate in the compartment.

- 3) Routine care and cleaning instructions: Wipe with an antibacterial wipe
- 4) <u>Precautions/Warnings:</u> Like any device, these casters must be used with caution. The casters provide an extra measure of safety but user vigilance is always required.
- 5) <u>Disposal:</u> Product should be disposed of according to applicable local laws.

rev 02, 25515 8/2021