


INSTRUCTIONS FOR USE



16401 UNI-HANGER 16402 SKIRT HANGER





1) <u>Indication for use:</u>

- a. <u>Uni-Hanger</u> allows you to hang any protective apron on a wall hook or in a conventional closet. The Uni-Hanger is 23½"W X 9"H.
- b. Skirt Hanger allows you to hang a skirt on any <u>Deka</u>- or <u>Penta-Rak</u>. The Skirt Hanger is 16" W x $5\frac{1}{2}$ " H.

<u>Instructions for Use</u>

Note: It is always recommended to

- a. To protect your apron on hanger after each use.
- b. To protect your protective apron use loops to hang on hanger.
- 3) Routine care and cleaning instructions: Wipe with an antibacterial wipe
- 4) Precautions/Warnings:
- 5) <u>Disposal:</u> Product should be disposed of according to applicable local laws.

rev 02, 16401, 16402 8/2021