



INSTRUCTIONS FOR USE

REF 50168, 50169, 50170, 50171, 5072, 50173, 50174: “Identifier” Lead Markers



- 1) Indication for use:
Identifier style markers feature high-density lead letters embedded in an injection molded, high-impact white polystyrene casing.

“Identifier” Lead Markers:

- 2) Instructions for Use
Coded red for "right," blue for "left" they are available with or without initials or numbers.
- 3) Routine care and cleaning instructions: Clean with disinfectant wipes.
- 4) Precautions/Warnings: California Proposition 65
- 5) Disposal: Product should be disposed of according to applicable local laws.
- 6) Adverse reactions: non known

rev 02, 50168, 50169, 50170, 50171, 5072, 50173, 50174
8/2021