

### INSTRUCTIONS FOR USE



# Lead Blockers 23308, 23309, 23310,23311, 23312, 23313, 23314, 23315



## 1) Indication for use:

**Lead Blockers** are designed for making more than one exposure on the same film. Lead Film Blockers are sized to fit half of various film sizes either length-wise or cross-wise and are made of 1/8-inch-thick lead rubber with 1 mm lead equivalency. Half-Film Blockers are available for film sizes 8 x 10, 10 x 12, 11 x 14, and 14 x 17 in.

#### **Lead Blockers:**

## 2) <u>Instructions for Use</u>

- Make multiple exposure on single film
- Cover half of film lengthwise or crosswise
- Constructed of 1/8-inch-thick lead rubber
- Provide 1mm lead protection
- Each orientation available in four sizes
- 3) Routine care and cleaning instructions: Cleans with most common disinfectants
- 4) Precautions/Warnings: California Proposition 65
- 5) <u>Disposal: Disposal:</u> Product should be disposed of according to applicable local laws.
- 6) Adverse reactions: non known

rev 01, 23308, 23309, 23310,23311, 23312, 23313, 23314, 23315

8/2021