

INSTRUCTIONS FOR USE 50150, L & R MARKER CLIP

REF



1) Indication for use:

The marker clip slips over the cassette frame and exposes "L" or "R" during patient exposure. Made of stainless steel, the Wolf Marker Clip is the industry standard. Please note that the marker clip is sold as a single clip.

L & R MARKER CLIP

2) <u>Instructions for Use</u>

Note: It is always recommended to

- A. Slips over the cassette frame and exposes "L" or "R"
- 3) Routine care and cleaning instructions: Wipe with an antibacterial wipe
- 4) Precautions/Warnings:
- 5) <u>Disposal:</u> Product should be disposed of according to applicable local laws.

rev 02, 50150 8/2021